CYSTICERCOSIS

Western blot IgG

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CYSTICERCOSIS

Western blot IgG

Immunoblot assay for
in vitro diagnostic use

#CYS-WB24G : 24 tests

#CYS-WB12G : 12 tests

#CYS-WB96G : 96 tests

INSTRUCTIONS FOR USE

Intended use

CYSTICERCOSIS IgG Western Blot (WB) is a qualitative test of serological IgG diagnosis by Immunoblot Assay of cysticercosis intended for confirmatory testing of a positive or equivocal result obtained through classic screening tests. It can be performed on sera or cerebrospinal fluid (CSF).

Principle of the test

Western Blot technique: The antigens (*Taenia solium* cysticerci extract of porcine origin), once separated by electrophoresis, are bound by electroblotting to the surface of a nitrocellulose membrane (called the transfer) cut into 24 strips numbered from 1 to 24.

Conduct of the test: Each sera (or CSF) specimen to be tested is separately incubated with a strip. The anti-cysticerci potentially present in the sample selectively bind themselves onto the antigens cysticerci. The alkaline phosphatase-anti human IgG conjugate then binds itself to the bound anti-cysticerci antibodies. Finally, the immunocomplexes react with the substrate. The antigens recognized by the anti-cysticerci antibodies of type IgG present in the samples are revealed as purple transversal bands.

Reagents supplied with the kit

italic: package of 12 tests (#CYS-WB12G) - bold: Package of 96 tests (#CYS-WB96G).

| ID | Qty | Description | Composition |
|-----|-----|---|---|
| R1 | 1 | Folder(s) of 24 (<i>12, 4x24</i>) STRIPS: precut + coloured Standards. (Each folder and each transfer is identified by a unique serial number) | Sensitized nitrocellulose. Coloured Molecular Weight (kDa): Blue: 250, Blue: 150, Blue: 100, Pink: 75, Blue: 50, Green: 37, Pink: 25, Blue: 20, Blue: 15, Yellow: 10. |
| R2 | 1 | Vial of 30 (<i>30, 125</i>) ml of SAMPLE BUFFER (Ready to use - pink solution). | Buffer + surfactant + NaN ₃ (<0.1%). |
| R3 | 1 | Vial(s) of 30 (<i>30, 2x60</i>) ml of ANTI IgG CONJUGATE (Ready to use - blue solution). | Buffer + anti-human IgG polyclonal goat sera conjugated with Alkaline Phosphatase + NaN ₃ (<0.1%) + stabilisers. |
| R5 | 1 | Vial of 30 (<i>30, 125</i>) ml of SUBSTRATE (Ready to use - opaque brown vial). | Buffer + NBT + BCIP + stabilisers. |
| R6 | 1 | Vial of 60 (<i>60, 250</i>) ml of WASH CONCENTRATE 10X BUFFER (<u>To be diluted 10 times</u> in distilled water - colourless solution). | Buffer + surfactant + NaN ₃ (<0.1%). |
| R10 | 1 | Tube of 200 (<i>200, 2x200</i>) µl of POSITIVE CONTROL SERUM (Ready to use - red cap). | Buffer + pool of human sera positive in cysticercosis serology + NaN ₃ (<0.1%) + stabilisers. |

R2, R3, R5 and R6 are common to all kits and have a unique lot number depending only on the date of their production. It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.

Additional material required but not provided

- Multi-channel polypropylene incubation trays for mini-blots (# WBPP- 08 or equivalent).
- Rocking platform for immunoblots, vacuum system for liquids (the # WBPP- 08 tubs that we supply can be emptied by simply turning them over).
- Tubes and material for drawing the samples, graduated cylinders, adapted containers. Automatic pipettes, micropipettes and disposable tips (volumes of 25 µl, 1.2 ml and 2 ml).
- Distilled or deionised water. Absorbent paper (e.g., Whatman filter paper), transparent adhesive tape.
- Latex gloves, tweezers to handle the strips, cutter or scalpel, flat transparent ruler.

Note: Our reagents can be used in an automated immunoblot processor. **Care should be taken with possible chemical contaminations of our reagents if the processor is shared with reagents from another manufacturer** (known example: contamination by the TWEEN 20), and bacterial contaminations. Reserve vials for the processor. After processing, do not place the remaining used reagents back into the original vials.

Storage and stability

Store between 2 and 8°C. The reagents from the kit are stable until the expiry date indicated on the outer box and the vial labels. Wash buffer diluted to 1/10 is stable for 2 months at +2 to +8 °C and one week at room temperature.

Precautions for use

Safety

- For *in vitro* use only. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- Positive control is a serum of human origin that has been screened and found negative for to HIV 1 and 2 antibodies and HCV antibodies, and HB antigen. However, it must be handled like a potentially infectious product.
- The substrate contains a mixture of NBT and BCIP, toxic on contact (skin and mucous membranes) and inhalation.
- The reagents contain sodium azide which can form explosive metallic salts with lead and copper. Rinse any spill with water.
- Dispose of waste (samples, tips, tubes, wash liquid, used reagent...) according to good practices used in the industry and current regulations in the country.

Precautions

- Do not use liquid reagents from different lots together.
- Use the strips in numerical order. Do not mix strips from different serial numbers; use the transfers in succession. Establish a specific distribution plan before starting the test.
- Do not touch the strips with your fingers; use tweezers.
- The reagents must be mixed well before use, particularly the concentrated wash buffer.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-channel contamination. Watch for the formation of foam or bubbles in the pipette tips (bacterial contamination of reagent vials).
- Clean incubation trays only with clear water followed by distilled water (never use detergent or bleach).
- The omission of a sample or the distribution of an inadequate volume may render the test result negative, regardless of its actual status.

Specimen collection

Aseptically collect the samples in dry tubes. A minimum of 25 µl of serum or CSF are required. In cases of CSF, using 50 µl will increase the sensitivity of the test.

Keep the samples at 2-8 °C until they are processed. If they need to be stored, freeze the samples at -20 ± 5 °C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

Preparation of reagents

Wash buffer: For 4 tests, in a clean bottle, dilute 10 ml of Wash Concentrate 10X (**R6**) in 90 ml of distilled or deionised water.

Test procedure

Nota Bene: It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.

1. Prepare a distribution plan for the samples and C+ positive control (**R10**).
Only by using this control can the test be technically validated and identification made, for a given serial number, of the specific bands developed. A C+ strip cannot be used to interpret the results of strips from a blot of a different serial number.
2. Cut the required number of strips (R1) using a scalpel and a clean and dry flat transparent ruler, keeping the blue positioning line on the strips: hold the strips firmly in place with the ruler and cut them on the side of the stain (the numbers are visible through the ruler).
3. Distribute 1.2 ml of sample buffer (R2) in each channel according to the established plan.
4. Deposit, in their numerical order, the numbered strips in the channels: Let the strips rehydrate themselves for approximately 1 minute, with the number visible at the top, by gently shaking the tray to totally immerse them in the buffer.
5. Distribute the samples and positive control(s): according to the distribution plan, at a rate of 25 µl per channel (preferably 50 µl for CSF). Gently shake the tray after each dispense.
Place the tray on a rocking platform at room temperature.
 - ✓ Serum: **Incubate for 90 min ± 5 min** at 18-25 °C.
 - ✓ CSF: **incubate overnight (12 hours +/- 2h)** at 18-25 °C. You must cover the incubation tray with a film to avoid the desiccation.
6. Wash step: Empty the contents of the channels with a Pasteur pipette or by turning the incubation tray over. Dispense 2 to 3 ml of diluted Wash Buffer in each channel. Incubate on the rocking platform for 3 min. Repeat 2 times, then empty the contents of the channels. Ensure that the strips don't turn during these steps.
7. Dispense 1.2 ml of anti IgG conjugate (R3) into each channel. Place the tray on the rocking platform.
Incubate for 60 min ± 5 min at 18-25 °C
8. Wash step: repeat step 6.
9. Distribute 1.2 ml of NBT/BCIP substrate (R5) into each of the channels. Place on the rocking platform and protect from direct light. **Incubate for 60 min ± 5 min** at 18-25 °C.

Regardless of the parameter, monitor the development of the colour. The development can be stopped if the background colour of the strip darkens to the point where reading is difficult (the quality of the wash steps has a fundamental influence on the background coloration). Note that the strips will lighten as they dry.

10. Stop the reaction by aspirating substrate with a Pasteur pipette or by turning the incubation tub over and dispensing 2 ml of distilled water in the channels. Repeat this last washing step one more time.
11. Drying the strips: With the channels still water-filled, take the strips by the numbered end using the tweezers and deposit them, with the number visible, onto a Whatman absorbent paper. Let air dry. The colour of the strips will naturally lighten while drying. Interpretation must only be performed after drying is complete.
12. Storage: Transfer the strips onto a sheet of paper, which will be used to archive them. Align the positioning lines. Keeping them in place with the flat ruler, stick the top of the strips with transparent adhesive tape.

For a good interpretation, the strips must be ordered by transfer and in their numerical order, spaced at a maximum of a few millimetres apart. It is unreliable to compare strips that are spaced far apart (e.g., no.2 with no.15). **It is dangerous** (false results) to compare strips from different kits (strips with different serial numbers).

Quality control and interpretation

The serum control (R10) provided with the kit must be systematically included in any immunoblot series. It shows the typical profile and allows for technical validation of the good conduct of the test (the bands must appear very clearly on the strip) and to calibrate precisely the position and aspect of the specific bands to allow interpretation of the results of the strips from the same transfer (same serial number).

- Description of the bands:

A positive sample can present numerous bands located between 2 and 200 kilodaltons (kDa). In practice and for specificity reasons, only the 6 to 55 kDa range is selected for the reading.

In this area, 5 bands are most often present at the following molecular weights (**kDa**): **6-8, 12, 23-26, 39, 50-55**. They are therefore called: **P6-8, P12, P23-26, P39 and P50-55**

Aspect of the bands:

*The **P6-8 and P23-26** bands can appear in the form of a large single band or a double band. The **P50-55** band traditionally presents itself in the form of a wide band with fairly fuzzy contours.*

Important points - **In Practice** (see Fig. 1):

- 1/ The 6-26 kDa and 39-55 kDa areas are the most specific and the most easily readable and interpretable.
- 2/ The intermediate area, delimited by the P23-26 and P39 bands, is not entirely specific to cysticercosis (frequent cross-reactions, particularly with other helminthiases and *P. falciparum* malaria).


Fig. 1: Examples of positive and negative results

- Interpretation:

The presence of a minimum of **2 well-defined bands** among the 5 previously described bands, P6-8, P12, P23-26, P39, and P50-55, is indicative of cysticercosis in the serum and of neuro-cysticercosis in the CSF. Above examples: “+” = neuro-cysticercosis - 8, 14, 11 = hydatidosis 7 = alveolar echinococcosis. **Note:** Strip 7 presents the non-specific “Mikado” aspect (cf. § Trouble shooting)

To validate the results, always compare the profile of the immunoblot of each sample with that of the R10 positive control. The aspect of the bands is important when interpreting the test.

Limitations of use

These serological results must be interpreted according to available information (epidemiological, clinical, imaging, biological) in order to establish a diagnosis.

Performances

- Sensitivity (Se):

The evaluation covered 79 samples (70 sera and 9 CSFs) that were positive according to clinical, epidemiological, radiological and/or serological criteria.

77 samples, including 9 CSFs, were found to be positive. Sensitivity Se = 97.5%

- Specificity (Sp):

The evaluation covered 95 samples, including 81 sera from patients suffering from the following parasitic infections: *Toxocara canis* (7), *Trichinella spiralis* (14), *Toxoplasma gondii* (7), filariasis (7), *Fasciola hepatica* (4), *Echinococcus granulosus* (14), *E. multilocularis* (14), *Schistosoma* (14) and 14 sera from patients suffering from autoimmune diseases: RF+ rheumatoid factor (7) and ANA+ anti-nuclear antibodies (7). All were found to be negative. Specificity Sp = 100%

Note: certain samples present isolated, narrow bands that should not be confused with specific bands (cf. examples on p. 5). In particular, the characteristic aspect (large and diffuse) of the **P50-55** band differentiates the narrow bands that are sometimes found at that level from echinococcosis, hydatidosis or schistosomiasis sera.

- Reproducibility:

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation with respect to specific bands is excellent.

- Interferences:

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

Trouble shooting

"The bands are pale with little contrast": Certain sera with low concentrations of antibodies may give such results.

"Shaded areas can be seen, more or less coloured, slightly diffuse": The strip was not totally submerged in one of the reagents and did not incubate correctly along its entire length. Stains can also be present where the sample was deposited if the tray was not shaken after dispensing.

"The background noise is significant, making reading very difficult": The washes were insufficient or the last incubation was too long. Ensure good test performance techniques, respect wash times and ensure water quality. Reduce the time of the last incubation.

Exceptionally, certain sera may react in a non-specific manner. The background staining can sometimes look like **streaks** (Mikado aspect, see example, Fig. 1, strip no. 7) that makes the immunoblot reading very difficult. Then, the result of the immunoblot cannot be used.

This non-specific background noise may involve only part of the strip, making the results uninterpretable for that part only.

"A precipitate appears in the solution during the last step of development": the substrate may in fact partially precipitate (black flakes) in the buffer at the end of development. This phenomenon does not alter the quality of the development which must be continued normally. The last wash with distilled water eliminates the possible solid particles present.

Bibliography

- Deckers, Nynke, et Pierre Dorny. 2010. « Immunodiagnosis of *Taenia Solium* Taeniosis/cysticercosis ». *Trends in Parasitology* 26 (3): 137–44. doi:10.1016/j.pt.2009.12.008.
- Del Brutto, Oscar H. 2012. « Diagnostic Criteria for Neurocysticercosis, Revisited ». *Pathogens and Global Health* 106 (5): 299–304. doi:10.1179/2047773212Y.0000000025.
- Dournon, Nathalie, Loic Epelboin, Marie-Charlotte Brion, Luc Paris, François Bricaire, et Eric Caumes. 2012. « Seroconversion of Neurocysticercosis Occurring After Anti-Helminthic Treatment: Neurocysticercosis With Seroconversion ». *Journal of Travel Medicine* 19 (6): 383–86. doi:10.1111/j.1708-8305.2012.00658.x.
- Garcia, Hector H, Theodore E Nash, et Oscar H Del Brutto. 2014. « Clinical Symptoms, Diagnosis, and Treatment of Neurocysticercosis ». *The Lancet Neurology* 13 (12): 1202–15. doi:10.1016/S1474-4422(14)70094-8.
- Gekeler, F, S Eichenlaub, E G Mendoza, J Sotelo, M Hoelscher, et T Löscher. 2002. « Sensitivity and specificity of ELISA and immunoblot for diagnosing neurocysticercosis ». *European journal of clinical microbiology & infectious diseases: official publication of the European Society of Clinical Microbiology* 21 (3): 227–29. doi:10.1007/s10096-002-0695-3.
- Michelet, Lorraine, Agnès Fleury, Edda Sciutto, Eric Kendjo, Gladis Fragoso, Luc Paris, et Bernard Bouteille. 2011. « Human neurocysticercosis: comparison of different diagnostic tests using cerebrospinal fluid ». *Journal of clinical microbiology* 49 (1): 195–200. doi:10.1128/JCM.01554-10.
- Raccurt, C P, P Agnamey, J Boncy, J-H Henrys, et A Totet. 2009. « Seroprevalence of human *Taenia solium* cysticercosis in Haiti ». *Journal of helminthology* 83 (2): 113–16. doi:10.1017/S0022149X09232330.
- Rodriguez, Silvia, Patricia Wilkins, et Pierre Dorny. 2012. « Immunological and Molecular Diagnosis of Cysticercosis ». *Pathogens and Global Health* 106 (5): 286–98. doi:10.1179/2047773212Y.0000000048.
- šOba, Barbara, Bojana Beović, Zala Lužnik, Miha Skvarč, et Jernej Logar. 2014. « Evidence of Human Neurocysticercosis in Slovenia ». *Parasitology* 141 (04): 547–53. doi:10.1017/S0031182013001947.
- Van Doorn, H. Rogier, Ellen Wentink-Bonnema, Rob J. Rentenaar, et Tom van Gool. 2007. « Specific Cross-Reactivity in Sera from Cystic Echinococcosis Patients in an Enzyme-Linked Immuno-electrotransfer Blot for Cysticercosis Diagnostics ». *Transactions of the Royal Society of Tropical Medicine and Hygiene* 101 (9): 948–50. doi:10.1016/j.trstmh.2007.04.021.