TOXOPLASMA

Western blot IgG IgM

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3. Garweg, J. G., de Groot-Mijnes, J. D. F. & Montoya, J. G. Diagnostic Approach to Ocular Toxoplasmosis. *Ocular Immunology and Inflammation* **19**, 255–261 (2011).
4. Garweg, J. G., Garweg, S.-D. L., Flueckiger, F., Jacquier, P. & Boehnke, M. Aqueous humor and serum immunoblotting for immunoglobulin types G, A, M, and E in cases of human ocular toxoplasmosis. *J. Clin. Microbiol.* **42**, 4593–4598 (2004).
5. Goldmann, H. & Witmer, R. [Antibodies in the aqueous humor]. *Ophthalmologica* **127**, 323–330 (1954).
6. L'ollivier, C. *et al.* Comparison of Mother and Child Antibodies That Target High-Molecular-Mass *Toxoplasma gondii* Antigens by Immunoblotting Improves Neonatal Diagnosis of Congenital Toxoplasmosis. *Clin. Vaccine Immunol.* **19**, 1326–1328 (2012).
7. Maenz, M. *et al.* Ocular toxoplasmosis past, present and new aspects of an old disease. *Prog Retin Eye Res* **39**, 77–106 (2014).
8. Magi, B. & Migliorini, L. Western blotting for the diagnosis of congenital toxoplasmosis. *New Microbiol.* **34**, 93–95 (2011).
9. Pinon, J. M. *et al.* Strategy for diagnosis of congenital toxoplasmosis: evaluation of methods comparing mothers and newborns and standard methods for postnatal detection of immunoglobulin G, M, and A antibodies. *J. Clin. Microbiol.* **39**, 2267–2271 (2001).
10. Potasman, I., Araujo, F. G. & Remington, J. S. *Toxoplasma* antigens recognized by naturally occurring human antibodies. *J. Clin. Microbiol.* **24**, 1050–1054 (1986).
11. Remington, J. S., Thulliez, P. & Montoya, J. G. Recent developments for diagnosis of toxoplasmosis. *J. Clin. Microbiol.* **42**, 941–945 (2004).
12. Rilling, V., Dietz, K., Krczal, D., Knotek, F. & Enders, G. Evaluation of a commercial IgG/IgM Western blot assay for early postnatal diagnosis of congenital toxoplasmosis. *Eur. J. Clin. Microbiol. Infect. Dis.* **22**, 174–180 (2003).
13. Robert-Gangneux, F. *et al.* Usefulness of immunoblotting and Goldmann-Witmer coefficient for biological diagnosis of toxoplasmic retinochoroiditis. *Eur. J. Clin. Microbiol. Infect. Dis.* **23**, 34–38 (2004).
14. Robert-Gangneux, F. & Darde, M.-L. Epidemiology of and Diagnostic Strategies for Toxoplasmosis. *Clinical Microbiology Reviews* **25**, 264–296 (2012).
15. Ronday, M. J., Ongkosuwito, J. V., Rothova, A. & Kijlstra, A. Intraocular anti-*Toxoplasma gondii* IgA antibody production in patients with ocular toxoplasmosis. *Am. J. Ophthalmol.* **127**, 294–300 (1999).
16. Talabani, H. *et al.* Contributions of Immunoblotting, Real-Time PCR, and the Goldmann-Witmer Coefficient to Diagnosis of Atypical Toxoplasmic Retinochoroiditis. *Journal of Clinical Microbiology* **47**, 2131–2135 (2009).
17. Tridapalli, E. *et al.* Congenital toxoplasmosis: the importance of the western blot method to avoid unnecessary therapy in potentially infected newborns. *Acta Paediatr.* **97**, 1298–1300 (2008).
18. Turunen, H. J., Leinikki, P. O. & Saari, K. M. Demonstration of intraocular synthesis of immunoglobulin G *Toxoplasma* antibodies for specific diagnosis of toxoplasmic chorioretinitis by enzyme immunoassay. *J. Clin. Microbiol.* **17**, 988–992 (1983).
19. Villard, O. *et al.* Comparison of Enzyme-Linked Immunosorbent Assay, Immunoblotting, and PCR for Diagnosis of Toxoplasmic Chorioretinitis. *Journal of Clinical Microbiology* **41**, 3537–3541 (2003).
20. Villard, O. *et al.* Serological diagnosis of *Toxoplasma gondii* infection: Recommendations from the French National Reference Center for Toxoplasmosis. *Diagn. Microbiol. Infect. Dis.* (2015).

TOXOPLASMA

Western blot IgG IgM

CE

Immunoblot assay for
in vitro diagnostic use

#TOP-WB24GM : 24 tests

#TOP-WB12GM : 12 tests

#TOP-WB96GM : 96 tests

INSTRUCTIONS FOR USE

Intended use

TOXOPLASMA WB IgG-IgM is an immunoblot assay for the Comparison of Immunological Profiles (CIP-WB) for IgG and IgM that is intended to diagnose:

- Congenital toxoplasmosis at birth (D0): CIP-WB G+M between maternal blood and cord blood.
- Congenital toxoplasmosis in post-natal monitoring (D+N): CIP-WB G+M between the cord blood at D0 and the child's blood at D+N.
- Ocular toxoplasmosis: CIP-WB IgG between the patient's serum and aqueous humour.

This assay is not intended for screening or confirming isolated serologies. For that application, use the LDBIO TOXO II IgG test (ref. TOXO II IgG WB).

Principle of the test

Western Blot technique: The antigens of *Toxoplasma gondii*, once separated by electrophoresis, are bound by electroblotting to the surface of a nitrocellulose membrane (called the transfer) cut into 24 strips numbered from 1 to 24.

Conduct of the test:

Note: The IgG or IgM immunoblot assays described below are conducted simultaneously during manipulation.

- IgG Immunoblot:

The assay consists of separately incubating, **with 2 contiguous strips from the same transfer**, the two samples (sera or aqueous humour) for which a comparison of immunological profiles is desired.

Step 1: Each sera (or aqueous humour) specimen to be tested is separately incubated with a strip. The anti-*Toxoplasma* antibodies potentially present in the sample selectively bind themselves onto the antigens of *T. gondii*.

Step 2: The alkaline phosphatase-**anti human IgG** conjugate then binds itself to the bound anti-antibodies.

Step 3: The immunocomplexes react with the substrate. The antigens recognized by the **class IgG** anti-*Toxoplasma* antibodies present in the samples are revealed as purple transversal bands.

- IgM immunoblot:

The principle of the assay is identical, but in step 2 the previous conjugate is replaced by an alkaline phosphatase-**anti-human IgM** conjugate. The colour development will therefore reveal the antigenic bands recognised by the **class IgM** anti-*Toxoplasma* antibodies present in the samples.

- Reading:

The successive comparison of pairs of IgG and then IgM (or IgA) strips allows to show the potential presence of bands that are only developed by one of the samples and not the other (cf. § Interpretation).

Reagents supplied with the kit

italic: package of 12 tests (#TOP-WB12GM) - bold: Package of 96 tests (#TOP-WB96GM).

ID	Qty	Description	Composition
R1	1	Folder(s) of 24 (<i>12</i> , 4x24) STRIPS: precut + coloured Standards. (Each folder and each transfer is identified by a unique serial number)	Sensitized nitrocellulose. Coloured Molecular Weight (kDa): Blue: 250, Blue: 150, Blue: 100, Pink: 75, Blue: 50, Green: 37, Pink: 25, Blue: 20, Blue: 15.
R2	1	Vial of 30 (<i>30</i> , 125) ml of SAMPLE BUFFER (Ready to use - pink solution).	Buffer + surfactant + NaN ₃ (<0.1%).
R3	1	Vial(s) of 30 (<i>30</i> , 60) ml of ANTI IgG CONJUGATE (Ready to use - blue solution).	Buffer + anti-human IgG polyclonal goat sera conjugated with Alkaline Phosphatase + NaN ₃ (<0.1%) + stabilisers.

R4	1	Vial(s) of 30 (30, 60) ml of ANTI-IgM CONJUGATE (ready to use - yellow solution).	Buffer + polyclonal goat anti-human IgM serum conjugated to alkaline phosphatase + NaN ₃ (< 0.1%) + stabilisers.
R5	1	Vial of 30 (30, 125) ml of SUBSTRATE (Ready to use - opaque brown vial).	Buffer + NBT + BCIP + stabilisers.
R6	1	Vial of 60 (60, 250) ml of WASH CONCENTRATE 10X BUFFER (To be diluted 10 times in distilled water - colourless solution).	Buffer + surfactant + NaN ₃ (<0.1%).

R2, R3, R5 and R6 are common to all kits and have a unique lot number depending only on the date of their production. It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.

Additional material required but not provided

- Multi-channel polypropylene incubation trays for mini-blots (# WBPP- 08 or equivalent).
- Rocking platform for immunoblots, vacuum system for liquids (the # WBPP- 08 tubs that we supply can be emptied by simply turning them over).
- Tubes and material for drawing the samples, graduated cylinders, adapted containers. Automatic pipettes, micropipettes and disposable tips (volumes of 10 µl, 25µl, 1.2 ml and 2 ml).
- Distilled or deionised water. Absorbent paper (e.g., Whatman filter paper), transparent adhesive tape.
- Latex gloves, tweezers to handle the strips, cutter or scalpel, flat transparent ruler.

Note: Our reagents can be used in an automated immunoblot processor. **Care should be taken with possible chemical contaminations of our reagents if the processor is shared with reagents from another manufacturer** (known example: contamination by the TWEEN 20), and bacterial contaminations. Reserve vials for the processor. After processing, do not place the remaining used reagents back into the original vials.

Storage and stability

Store between 2 and 8°C. The reagents from the kit are stable until the expiry date indicated on the outer box and the vial labels. Wash buffer diluted to 1/10 is stable for 2 months at +2 to +8 °C and one week at room temperature.

Precautions for use

Safety

- For *in vitro* use only. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- The substrate contains a mixture of NBT and BCIP, toxic on contact (skin and mucous

membranes) and inhalation.

- The reagents contain sodium azide which can form explosive metallic salts with lead and copper. Rinse any spill with water.
- Dispose of waste (samples, tips, tubes, wash liquid, used reagent...) according to good practices used in the industry and current regulations in the country.

Precautions

- Do not use liquid reagents from different lots together.
- Use the strips in numerical order. Do not mix strips from different serial numbers; use the transfers in succession. Establish a specific distribution plan before starting the test.
- Do not touch the strips with your fingers; use tweezers.
- The reagents must be mixed well before use, particularly the concentrated wash buffer.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-channel contamination. Watch for the formation of foam or bubbles in the pipette tips (bacterial contamination of reagent vials).
- Clean incubation trays only with clear water followed by distilled water (never use detergent or bleach).
- The omission of a sample or the distribution of an inadequate volume may render the test result negative, regardless of its actual status.

Specimen collection

Aseptically collect the samples in dry tubes. A minimum of 35 µl of serum or 10 µl of aqueous humour are required. In cases of aqueous humour, using 25 µl will increase the sensitivity of the test.

Keep the samples at 2-8 °C until they are processed. If they need to be stored, freeze the samples at -20 ± 5 °C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

Preparation of reagents

Wash buffer: For 4 tests, in a clean bottle, dilute 10 ml of Wash Concentrate 10X (R6) in 90 ml of distilled or deionised water.

Test procedure

Nota Bene: It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.

1. Carefully prepare the sample distribution plan.
It is strictly mandatory to compare a pair of samples with joined strips (contiguous numbers) from a given transfer (same serial number). . It is unreliable to compare strips that are spaced far apart (e.g., no.2 with no.15). **It is dangerous** (false results) to compare strips from different kits (strips with different serial numbers).
2. Cut the required number of strips (R1) using a scalpel and a clean and dry flat transparent ruler, keeping the blue positioning line on the strips: hold the strips firmly in place with the ruler and cut them on the side of the stain (the numbers are visible through the ruler).
3. Distribute 1.2 ml of sample buffer (R2) in each channel according to the established plan.

4. Deposit, in their numerical order, the numbered strips in the channels: Let the strips rehydrate themselves for approximately 1 minute, with the number visible at the top, by gently shaking the tray to totally immerse them in the buffer.
5. Dispense the samples according to the established distribution plan (step 1) and the following volumes:

	Serum	Aqueous humour
IgG	10µl	10 or 25µl
IgM	25µl	-

In cases of aqueous humour, using 25 µl will increase the sensitivity of the test.

Gently shake the tray after each dispense. Place the tray on a rocking platform. **Incubate for 90 min** ± 5 min at 18-25 °C.

6. Wash step: Empty the contents of the channels with a Pasteur pipette or by turning the incubation tray over. Dispense 2 to 3 ml of diluted Wash Buffer in each channel. Incubate on the rocking platform for 3 min. Repeat 2 times, then empty the contents of the channels. Ensure that the strips don't turn during these steps.
7. Dispense, according to the established distribution plan, 1.2 ml of anti-IgG conjugate (R3) or 1.2 ml of anti-IgM conjugate (R4) in each of the corresponding wells. Place the tray on the rocking platform.
Incubate for 60 min ± 5 min at 18-25 °C
8. Wash step: repeat step 6.
9. Distribute 1.2 ml of NBT/BCIP substrate (R5) into each of the channels. Place on the rocking platform and protect from direct light. **Incubate for 60 min** ± 5 min at 18-25 °C.

Regardless of the parameter, monitor the development of the colour. The development can be stopped if the background colour of the strip darkens to the point where reading is difficult (the quality of the wash steps has a fundamental influence on the background coloration). **Note that the strips will lighten as they dry.**

- It is essential to simultaneously stop the colour development of the 2 strips of a given pair for a given sub-class of antibodies, but one can independently stop that of the IgGs or IgMs (IgMs, in a lower concentration, usually develop more slowly than IgGs).
- Children's serum generally has a lower concentration of IgM. The reaction must be allowed time to correctly develop and one should not be concerned to see the maternal IgM strip darken a little more.
- Aqueous humour generally has a lower concentration of antibodies. The reaction must be allowed time to correctly develop and one should not be concerned to see the serum strips darken a little more.

10. Stop the reaction by aspirating substrate with a Pasteur pipette or by turning the incubation tub over and dispensing 2 ml of distilled water in the channels. Repeat this last washing step one more time.
11. Drying the strips: With the channels still water-filled, take the strips by the numbered end using the tweezers and deposit them, with the number visible, onto a Whatman absorbent paper. Let air dry. The colour of the strips will naturally lighten while drying. Interpretation must only be performed after drying is complete.

12. Storage: Transfer the strips onto a sheet of paper, which will be used to archive them.

- Match, side-by-side, the IgG and IgM strips of each pair of samples in order of increasing numbers, by following the established distribution plan (step 1).
- Precisely adjust the strips with the positioning line.
- Keep them in place, 2 by 2, with the flat ruler and stick the top of the strips with transparent adhesive tape.

For a good interpretation, the strips must be ordered by transfer and in their numerical order, spaced at a maximum of a few millimetres apart. It is unreliable to compare strips that are spaced far apart (e.g., no.2 with no.15). It is dangerous (false results) to compare strips from different kits (strips with different serial numbers).

Interpretation

Description of the bands:

A positive sample can present a significant number of bands located between 15 and 200 kDa. Only the bands with a molecular weight less than 120 kDa can be used to compare the profiles.

1- Interpretation CIP WB G+M (congenital toxoplasmosis)

- **At birth (mother/child pairs):**

Independently compare the IgG strips and the IgM strips.

Read the 2 contiguous strips simultaneously from top to bottom while noting any antigenic band that is **present** in the cord blood **and absent** from the maternal serum.

Any band that has a well-defined resolution, a Molecular Weight (MW) less than 120 kDa and that is *only present in the child* is evidence that the child has synthesised anti-toxoplasma antibodies, suggesting congenital toxoplasmosis.

- **During post-natal monitoring (child D0/child D+N pairs):**

Independently compare the IgG strips and the IgM strips.

Read the 2 contiguous strips simultaneously from top to bottom while noting any antigenic band that is **present** in the serum at D+N **and absent** from the cord blood.

Any band that has a well-defined resolution, an MW < 120 kDa and that is *only present at D+N* is evidence that the child has synthesised anti-toxoplasma antibodies, suggesting congenital toxoplasmosis.

NB: the indication of CIP-WB IgG/IGM in post-natal monitoring is intentionally limited to 3 months for IgG and 1 month for IgM.

Notes: The juxtaposition of the stained molecular weight standard (folder R1) allows the MW of the developed antigenic bands to be estimated (it must be cut beforehand with a ruler and scalpel, like an ordinary strip, and handled with tweezers).

A - NEGATIVE CIP
(not- infected child)

B - POSITIVE CIP
(infected child)


Fig. 1: Congenital toxoplasmosis - Examples of positive and negative results

The mother-child pair (A) corresponds to a mother infected during pregnancy but whose child is uninfected: The IgG profiles are strictly identical (transmitted IgGs); there is no other additional band present on the child's IgG and/or IgM strips: **THE CIP-WB IS NEGATIVE.**

The pair (B), congenital toxoplasmosis, corresponds to a mother infected during her pregnancy and whose child was also infected. Besides the transmitted antibodies, one perfectly notes the presence of additional bands (↔), for IgG and/or IgM, on the child's strips, corresponding to the antibodies newly synthesised by the child: **THE CIP-WB IS POSITIVE.**

2- Interpretation CIP WB IgG (ocular toxoplasmosis)

Read the 2 contiguous strips simultaneously from top to bottom while noting any antigenic band **present** in the aqueous humour **and absent** from the serum.

Any band that has a well-defined resolution, a Molecular Weight (MW) less than 120 kDa and that is *only present in the aqueous humour* evidences local synthesis of anti-toxoplasma antibodies, suggesting ocular toxoplasmosis.


Fig. 2: Ocular toxoplasmosis - Examples of positive results

3- Very important points

1. The CIP-WB IgG/IgM results must be interpreted in light of other clinical, serological, parasitological, epidemiological and medical imaging information in order to establish the diagnosis of congenital or ocular toxoplasmosis.

2. A negative CIP-WB IgG/IgM result does not rule out the diagnosis of congenital or ocular toxoplasmosis. These patients must always be monitored over time until the diagnosis of toxoplasmosis can be definitively confirmed or ruled out.
3. The aspects of the bands can vary greatly: narrow, thick, more or less coloured, intense, etc. When one takes on this technique, it is recommended that several comparisons of profiles be made with known pairs of samples in order to become familiar with their reading. At the beginning, it is also recommended that the CIP-WB reading be done independently by two individuals in the laboratory. In case of discordant interpretations, a control CIP-WB must be performed.
4. The very high molecular weight (MW) antigenic fractions are very close together in the upper part of the strip in favour of better resolution of the medium and low MW fractions. The bands of MW > 120 kDa therefore cannot be used to interpret the assay: samples only presenting such profile differences cannot be returned as positive.
5. In contrast (congenital toxoplasmosis), a “triplet” (three very easily recognisable bands) located between 75 and 100 kDa is very often found on positive CIP-WB IgMs (see “T” Fig. 1, strip No. 17 to the right).
6. At birth (congenital toxoplasmosis), one must pay particular attention to any general reinforcement of the intensity of the bands (haemoconcentration) that could suggest that additional bands exist in the cord blood. Sera only presenting such profile differences are returned as negative.
7. In contrast (congenital toxoplasmosis, ocular toxoplasmosis), significant reinforcement (often in width and intensity) of one or two isolated band(s), while all the other bands are of identical or weaker intensity, is considered to be a criterion for positivity.
8. Natural antibodies (congenital toxoplasmosis):
The immunoblot technique is extremely sensitive and the antigen used for the CIP-WB test was selected for the multiplicity of antigenic bands present on the strip. Numerous publications mention bands developed by immunoblot in individuals who have apparently never contracted toxoplasmosis. These antibodies (IgG and IgM) are only rarely detected by other techniques but are very frequently detected by immunoblot. They could be due to cross-reactions with antibodies directed against immunogens of a nature that has yet to be determined. [15, 16, 18]
This is why the indication for the TOXOPLASMA WB IgG-IgM test is reserved to comparing profiles. (To confirm IgG serologies, use the specific LDBIO TOXO II IgG test that is intended for that use)
New-borns do not present natural antibodies (other than transmitted maternal antibodies), but the likelihood of natural antibodies appearing increases with the infant’s age after 3 months; they are only rarely found between 3 and 6 months. [15]
This is why the indication for CIP-WB IgG/IgM in post-natal monitoring is intentionally limited to 3 months for IgG and 1 month for IgM: non-specific bands do in fact appear earlier for IgM.
9. “Heat Shock Protein” (congenital toxoplasmosis):
A non-specific, narrow band of weak but variable intensity can be present for IgM up to 37 kDa. It is an artefact linked to the preparation of the antigen and called “heat shock protein”. Present on both the strips of the mother-child pair, it can nevertheless sometimes appear to be more pronounced with certain sera during monitoring of the child. Do not take this band into account.
10. CIP-WB (ocular toxoplasmosis): The CIP-WB IgM is of no use in the diagnosis of ocular toxoplasmosis.

Performances

❖ Specificity - Sensitivity:

These studies were conducted by independent reference laboratories

1/ CIP-WB G+M: CONGENITAL TOXOPLASMOSIS at birth (mother/child)

		TOXOPLASMA WB IgG-IgM	
		POS	NEG
CLINICAL DATA	POS CT n = 54	41	13
	NEG CT n = 60	0	60

Table 1: Performance of CIP-WB IgG/IgM at birth (n = 114):

- Specificity = 100%
- Sensitivity = 76%
- Positive predictive value = 100%
- Negative predictive value = 83%

2/ CIP-WB G+M: CONGENITAL TOXOPLASMOSIS in post natal monitoring (child D0/D20)

Out of the 54 children previously tested at D0 (*Fig. 1*), 10 uninfected children and 12 infected children (n = 22) were monitored until D20 and retrospectively analysed with the TOXOPLASMA WB IgG-IgM test.

At D0: 4 out of 12 infected children did not present a profile different from birth (false negatives).

At D20: 1 single child remains negative.

		TOXOPLASMA WB IgG-IgM	
		POS	NEG
CLINICAL DATA	POS CT n = 12	11	1
	NEG CT n = 10	0	10

Table 2: Performance of CIP-WB IgG/IgM at D20 (n = 22):

- Specificity = 100%
- Sensitivity = 92%
- Positive predictive value = 100%
- Negative predictive value = 91%

3/ CIP-WB IgG: OCULAR TOXOPLASMOSIS (serum/aqueous humour)

The performances shown below are from the meta-analysis of four studies published by reference centres [1, 13, 16, 19].

These studies compare the performances of CIP-WB IgG to those of the Goldmann Witmer Coefficient (GWC) and those of PCR. They also show the diagnostic performances obtained by the combined association of two or three of these techniques.

These four studies all used the LDBIO test in accordance with the recommendations in the kit's instructions for use.

Sensitivity was determined on 113 patients presenting clinically proven ocular toxoplasmosis. Specificity was calculated on a control population presenting an ocular condition other than toxoplasmic infection: ocular toxocariasis (n=5), viral infection (n=10), other infections (n=4), non-infectious ocular conditions (n=126) of which cataract (n=42).

✓ Sensitivity (Se):

- The overall sensitivity of CIP-WB IgG is 62.8% (n=113), performance that is comparable to the GWC (Se=61.0%, n=113) and greater than PCR (Se=43.5%, n=92, p=0.0028).
- The combination of CIP-WB with GWC and PCR improves the sensitivity of the diagnosis:
 - CIP-WB + GWC: Se=78.1% (n=96, p=0.0082)
 - CIP-WB + GWC + PCR: 86.3% (n=95, p=0.0001)

✓ Specificity (Sp):

- The overall specificity of CIP-WB IgG is 92.8% (n=111), performance that is comparable to the GWC (Sp=94.2%, n=139) and less than PCR (Sp=100%, n=131, p=0.0009).
- The combination of the two techniques, CIP-WB IgG + GWC, slightly reduces the specificity of the diagnosis (Sp=91.1%, n=101, p=0.32). Combination with PCR does not influence the specificity.

❖ Reproducibility:

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation with respect to specific bands is excellent.

❖ Interferences:

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

Trouble shooting

"The bands are pale with little contrast": Certain sera with low concentrations of antibodies may give such results.

"Shaded areas can be seen, more or less coloured, slightly diffuse": The strip was not totally submerged in one of the reagents and did not incubate correctly along its entire length. Stains can also be present where the sample was deposited if the tray was not shaken after dispensing.

"The background noise is significant, making reading very difficult": The washes were insufficient or the last incubation was too long. Ensure good test performance techniques, respect wash times and ensure water quality. Reduce the time of the last incubation.

Exceptionally, certain sera may react in a non-specific manner. Then, the result of the immunoblot cannot be used.

This non-specific background noise may involve only part of the strip, making the results uninterpretable for that part only.

"A precipitate appears in the solution during the last step of development": the substrate may in fact partially precipitate (black flakes) in the buffer at the end of development. This phenomenon does not alter the quality of the development which must be continued normally. The last wash with distilled water eliminates the possible solid particles present.

Bibliography

1. Fekkar, A. *et al.* Comparison of immunoblotting, calculation of the Goldmann-Witmer coefficient, and real-time PCR using aqueous humor samples for diagnosis of ocular toxoplasmosis. *J. Clin. Microbiol.* **46**, 1965–1967 (2008).
2. Garweg, J. G. Determinants of immunodiagnostic success in human ocular toxoplasmosis. *Parasite Immunol.* **27**, 61–68 (2005).

3. Garweg, J. G., de Groot-Mijnes, J. D. F. & Montoya, J. G. Diagnostic Approach to Ocular Toxoplasmosis. *Ocular Immunology and Inflammation* **19**, 255–261 (2011).
4. Garweg, J. G., Garweg, S.-D. L., Flueckiger, F., Jacquier, P. & Boehnke, M. Aqueous humor and serum immunoblotting for immunoglobulin types G, A, M, and E in cases of human ocular toxoplasmosis. *J. Clin. Microbiol.* **42**, 4593–4598 (2004).
5. Goldmann, H. & Witmer, R. [Antibodies in the aqueous humor]. *Ophthalmologica* **127**, 323–330 (1954).
6. L'ollivier, C. *et al.* Comparison of Mother and Child Antibodies That Target High-Molecular-Mass *Toxoplasma gondii* Antigens by Immunoblotting Improves Neonatal Diagnosis of Congenital Toxoplasmosis. *Clin. Vaccine Immunol.* **19**, 1326–1328 (2012).
7. Maenz, M. *et al.* Ocular toxoplasmosis past, present and new aspects of an old disease. *Prog Retin Eye Res* **39**, 77–106 (2014).
8. Magi, B. & Migliorini, L. Western blotting for the diagnosis of congenital toxoplasmosis. *New Microbiol.* **34**, 93–95 (2011).
9. Pinon, J. M. *et al.* Strategy for diagnosis of congenital toxoplasmosis: evaluation of methods comparing mothers and newborns and standard methods for postnatal detection of immunoglobulin G, M, and A antibodies. *J. Clin. Microbiol.* **39**, 2267–2271 (2001).
10. Potasman, I., Araujo, F. G. & Remington, J. S. *Toxoplasma* antigens recognized by naturally occurring human antibodies. *J. Clin. Microbiol.* **24**, 1050–1054 (1986).
11. Remington, J. S., Thulliez, P. & Montoya, J. G. Recent developments for diagnosis of toxoplasmosis. *J. Clin. Microbiol.* **42**, 941–945 (2004).
12. Rilling, V., Dietz, K., Krczal, D., Knotek, F. & Enders, G. Evaluation of a commercial IgG/IgM Western blot assay for early postnatal diagnosis of congenital toxoplasmosis. *Eur. J. Clin. Microbiol. Infect. Dis.* **22**, 174–180 (2003).
13. Robert-Gangneux, F. *et al.* Usefulness of immunoblotting and Goldmann-Witmer coefficient for biological diagnosis of toxoplasmic retinochoroiditis. *Eur. J. Clin. Microbiol. Infect. Dis.* **23**, 34–38 (2004).
14. Robert-Gangneux, F. & Darde, M.-L. Epidemiology of and Diagnostic Strategies for Toxoplasmosis. *Clinical Microbiology Reviews* **25**, 264–296 (2012).
15. Ronday, M. J., Ongkosuwito, J. V., Rothova, A. & Kijlstra, A. Intraocular anti-*Toxoplasma gondii* IgA antibody production in patients with ocular toxoplasmosis. *Am. J. Ophthalmol.* **127**, 294–300 (1999).
16. Talabani, H. *et al.* Contributions of Immunoblotting, Real-Time PCR, and the Goldmann-Witmer Coefficient to Diagnosis of Atypical Toxoplasmic Retinochoroiditis. *Journal of Clinical Microbiology* **47**, 2131–2135 (2009).
17. Tridapalli, E. *et al.* Congenital toxoplasmosis: the importance of the western blot method to avoid unnecessary therapy in potentially infected newborns. *Acta Paediatr.* **97**, 1298–1300 (2008).
18. Turunen, H. J., Leinikki, P. O. & Saari, K. M. Demonstration of intraocular synthesis of immunoglobulin G *Toxoplasma* antibodies for specific diagnosis of toxoplasmic chorioretinitis by enzyme immunoassay. *J. Clin. Microbiol.* **17**, 988–992 (1983).
19. Villard, O. *et al.* Comparison of Enzyme-Linked Immunosorbent Assay, Immunoblotting, and PCR for Diagnosis of Toxoplasmic Chorioretinitis. *Journal of Clinical Microbiology* **41**, 3537–3541 (2003).
20. Villard, O. *et al.* Serological diagnosis of *Toxoplasma gondii* infection: Recommendations from the French National Reference Center for Toxoplasmosis. *Diagn. Microbiol. Infect. Dis.* (2015).