# LDBIO TOXO II IgG

## CONFIRMATION

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
# LDBIO TOXO II IgG CONFIRMATION

Immunoblot assay for  
*in vitro* diagnostic use

#TOXO II-24G : 24 tests

#TOXO II-12G : 12 tests

#TOXO II-96G : 96 tests

## INSTRUCTIONS FOR USE

### Intended use

LDBIO TOXO II IgG is a qualitative test of serological IgG diagnosis by Immunoblot Assay of toxoplasmosis intended for confirmatory testing of a positive or equivocal result obtained through classic screening tests. It can be performed on sera, cerebrospinal fluid (CSF) or aqueous humour.

### Principle of the test

Western Blot technique: The antigens of *Toxoplasma gondii*, once separated by electrophoresis, are bound by electroblotting to the surface of a nitrocellulose membrane (called the transfer) cut into 24 strips numbered from 1 to 24.

Conduct of the test: Each sera (CSF or aqueous humour) specimen to be tested is separately incubated with a strip. The anti-*Toxoplasma* antibodies potentially present in the sample selectively bind themselves onto the antigens of *T. gondii*. The alkaline phosphatase-anti human IgG conjugate then binds itself to the bound anti-*Toxoplasma* antibodies. Finally, the immunocomplexes react with the substrate. The antigens recognized by the anti-*Toxoplasma* antibodies of type IgG present in the samples are revealed as purple transversal bands.

## Reagents supplied with the kit

*italic: package of 12 tests (#TOXO II 12G) - bold: Package of 96 tests (#TOXO II 96G).*

ID	Qty	Description	Composition
R1	1	Folder(s) of 24 ( <i>12, 4x24</i> ) STRIPS: precut + coloured Standards. (Each folder and each transfer is identified by a unique serial number)	Sensitized nitrocellulose. Coloured Molecular Weight (kDa): Blue: 250, Blue: 150, Blue: 100, Pink: 75, Blue: 50, Green: 37, Pink: 25, Blue: 20, Blue: 15.
R2	1	Vial of 30 ( <i>30, 125</i> ) ml of SAMPLE BUFFER (Ready to use - pink solution).	Buffer + surfactant + NaN <sub>3</sub> (<0.1%).
R3	1	Vial(s) of 30 ( <i>30, 2x60</i> ) ml of ANTI IgG CONJUGATE (Ready to use - blue solution).	Buffer + anti-human IgG polyclonal goat sera conjugated with Alkaline Phosphatase + NaN <sub>3</sub> (<0.1%) + stabilisers.
R5	1	Vial of 30 ( <i>30, 125</i> ) ml of SUBSTRATE (Ready to use - opaque brown vial).	Buffer + NBT + BCIP + stabilisers.
R6	1	Vial of 60 ( <i>60, 250</i> ) ml of WASH CONCENTRATE 10X BUFFER (To be diluted 10 times in distilled water - colourless solution).	Buffer + surfactant + NaN <sub>3</sub> (<0.1%).
R10	1	Tube of 100 ( <i>100, 2x100</i> ) µl of POSITIVE CONTROL SERUM (Ready to use - red cap).	Buffer + pool of human sera positive in <i>Toxoplasma</i> serology + NaN <sub>3</sub> (<0.1%) + stabilisers.

**R2, R3, R5 and R6 are common to all kits and have a unique lot number depending only on the date of their production.** It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.

## Additional material required but not provided

- Multi-channel polypropylene incubation trays for mini-blots (# WBPP- 08 or equivalent).
- Rocking platform for immunoblots, vacuum system for liquids (the # WBPP- 08 tubs that we supply can be emptied by simply turning them over).
- Tubes and material for drawing the samples, graduated cylinders, adapted containers. Automatic pipettes, micropipettes and disposable tips (volumes of 10 µl, 25µl, 1.2 ml and 2 ml).
- Distilled or deionised water. Absorbent paper (e.g., Whatman filter paper), transparent adhesive tape.
- Latex gloves, tweezers to handle the strips, cutter or scalpel, flat transparent ruler.

Note: Our reagents can be used in an automated immunoblot processor. **Care should be taken with possible chemical contaminations of our reagents if the processor is shared with reagents from another manufacturer** (known example: contamination by the TWEEN 20), and bacterial contaminations. Reserve vials for the processor. After processing, do not place the remaining used reagents back into the original vials.

## Storage and stability

Store between 2 and 8°C. The reagents from the kit are stable until the expiry date indicated on the outer box and the vial labels. Wash buffer diluted to 1/10 is stable for 2 months at +2 to +8 °C and one week at room temperature.

## Precautions for use

### Safety

- For *in vitro* use only. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- Positive control is a serum of human origin that has been screened and found negative for to HIV 1 and 2 antibodies and HCV antibodies, and HB antigen. However, it must be handled like a potentially infectious product.
- The substrate contains a mixture of NBT and BCIP, toxic on contact (skin and mucous membranes) and inhalation.
- The reagents contain sodium azide which can form explosive metallic salts with lead and copper. Rinse any spill with water.
- Dispose of waste (samples, tips, tubes, wash liquid, used reagent...) according to good practices used in the industry and current regulations in the country.

### Precautions

- Do not use liquid reagents from different lots together.
- Use the strips in numerical order. Do not mix strips from different serial numbers; use the transfers in succession. Establish a specific distribution plan before starting the test.
- Do not touch the strips with your fingers; use tweezers.
- The reagents must be mixed well before use, particularly the concentrated wash buffer.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-channel contamination. Watch for the formation of foam or bubbles in the pipette tips (bacterial contamination of reagent vials).
- Clean incubation trays only with clear water followed by distilled water (never use detergent or bleach).
- The omission of a sample or the distribution of an inadequate volume may render the test result negative, regardless of its actual status.

## Specimen collection

Aseptically collect the samples in dry tubes. A minimum of 10 µl of serum, aqueous humour or CSF are required. In cases of aqueous humour or CSF, using 25 µl will increase the sensitivity of the test. Keep the samples at 2-8 °C until they are processed. If they need to be stored, freeze the samples at -20 ± 5 °C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

## Preparation of reagents

**Wash buffer:** For 4 tests, in a clean bottle, dilute 10 ml of Wash Concentrate 10X (**R6**) in 90 ml of distilled or deionised water.

## Test procedure

*Nota Bene: It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.*

1. Prepare a distribution plan for the samples and C+ positive control (**R10**).  
Only by using this control can the test be technically validated and identification made, for a given serial number, of the specific bands developed. A C+ strip cannot be used to interpret the results of strips from a blot of a different serial number.
2. Cut the required number of strips (R1) using a scalpel and a clean and dry flat transparent ruler, keeping the blue positioning line on the strips: hold the strips firmly in place with the ruler and cut them on the side of the stain (the numbers are visible through the ruler).
3. Distribute 1.2 ml of sample buffer (R2) in each channel according to the established plan.
4. Deposit, in their numerical order, the numbered strips in the channels: Let the strips rehydrate themselves for approximately 1 minute, with the number visible at the top, by gently shaking the tray to totally immerse them in the buffer.
5. Distribute the samples and positive control(s): according to the distribution plan, at a rate of 10 µl per channel (preferably 25 µl for aqueous humour or CSP). Gently shake the tray after each dispense. Place the tray on a rocking platform. **Incubate for 90 min** ± 5 min at 18-25 °C.
6. Wash step: Empty the contents of the channels with a Pasteur pipette or by turning the incubation tray over. Dispense 2 to 3 ml of diluted Wash Buffer in each channel. Incubate on the rocking platform for 3 min. Repeat 2 times, then empty the contents of the channels. Ensure that the strips don't turn during these steps.
7. Dispense 1.2 ml of anti IgG conjugate (R3) into each channel. Place the tray on the rocking platform.  
**Incubate for 60 min** ± 5 min at 18-25 °C
8. Wash step: repeat step 6.
9. Distribute 1.2 ml of NBT/BCIP substrate (R5) into each of the channels. Place on the rocking platform and protect from direct light. **Incubate for 60 min** ± 5 min at 18-25 °C.

Regardless of the parameter, monitor the development of the colour. The development can be stopped if the background colour of the strip darkens to the point where reading is difficult (the quality of the wash steps has a fundamental influence on the background coloration). Note that the strips will lighten as they dry.

10. Stop the reaction by aspirating substrate with a Pasteur pipette or by turning the incubation tub over and dispensing 2 ml of distilled water in the channels. Repeat this last washing step one more time.
11. Drying the strips: With the channels still water-filled, take the strips by the numbered end using the tweezers and deposit them, with the number visible, onto a Whatman absorbent paper. Let air dry.

The colour of the strips will naturally lighten while drying. Interpretation must only be performed after drying is complete.

12. Storage: Transfer the strips onto a sheet of paper, which will be used to archive them. Align the positioning lines. Keeping them in place with the flat ruler, stick the top of the strips with transparent adhesive tape.

For a good interpretation, the strips must be ordered by transfer and in their numerical order, spaced at a maximum of a few millimetres apart. It is unreliable to compare strips that are spaced far apart (e.g., no.2 with no.15). **It is dangerous** (false results) to compare strips from different kits (strips with different serial numbers).


## Quality control and interpretation

The serum control (R10) provided with the kit must be systematically included in any immunoblot series. It shows the typical profile and allows for technical validation of the good conduct of the test (the bands must appear very clearly on the strip) and to calibrate precisely the position and aspect of the specific bands to allow interpretation of the results of the strips from the same transfer (same serial number).

- Description of the bands:

A positive sample can present numerous bands located between 15 and 200k kilodaltons (kDa). Search for the presence of specific bands in the 30-45 kDa area for each of the tested samples with the calibration tools described above.

These bands, grouped and well isolated, are typical and can generally be found very easily.

**Fig. 1:** Examples of positive and negative results

Left and right: picture of the positive control of the kit (R10), which identifies the specific bands in the 30-45 kDa area: 30, 31, 33, 40 and 45.

- Interpretation:

The presence on the strip of a minimum of 3 bands out of specific bands 30, 31, 33, 40 and 45, and the inclusion of the band at 30 kDa, allows the assay to be interpreted as positive and to conclude that anti-*T. gondii* IgG antibodies are present in the tested sample.

A: seroconversion example. Sera 2 and 3, positive with LDBIO-TOXO II IgG, were negative with the screening technique (named *ELISA 2 IgG* below in the performance study).

B: example in neo-natal monitoring. Sera 5 and 6, positive with LDBIO-TOXO II IgG, were negative with the *ELISA 2 IgG* screening technique.

NB: other bands may be observed. They are not taken into account in the reading of the assay.

To validate the results, always compare the profile of the immunoblot of each sample with that of the R10 positive control. The aspect of the bands is important when interpreting the test.

## Limitations of use

These serological results must be interpreted according to available information (epidemiological, clinical, imaging, biological) in order to establish a diagnosis.

## Performances

The evaluation was done in a reference laboratory specialising in the diagnosis of toxoplasmosis.

**The principle of the evaluation consisted of comparing, the results of 529 sera obtained with the LDBIO-TOXO II IgG technique, the results obtained with Sabin and Feldman's DYE TEST, the results of two marketed screening techniques, "ELISA 1 IgG" and "ELISA 2 IgG", as well as the patients' clinical and biological data.**

### Threshold of the techniques used:

	NEGATIVE	EQUIVOCAL	POSITIVE
DYE TEST (IU/ml)	< 2	-	≥ 2
ELISA 1 (IU/ml)	< 4	4 - 8	≥ 8
ELISA 2 (IU/ml)	< 6	-	≥ 6
LDBIO TOXO II IgG	0	-	≥ 1

### Statistical analysis of the results:

We established the sensitivity and specificity values when possible. The correlation between the results found with the various techniques was evaluated with McNemar's CHI-2 test on matched series.

- PATIENTS:

All the analyses were done on sera stored frozen at -20°C. The samples come from 5 different groups of patients.

- Group I – DYE TEST

Study on 200 sera obtained during toxoplasmosis screening in pregnant women and tested by the DYE TEST. The "positive" sub-group corresponds to 98 sera, positive with the DYE TEST, from women immunised against *T. gondii*. This sub-group included sera with moderate titres of IgG with the DYE TEST (from 2 to 32 IU/ml) in order to test the sensitivity of LDBIO-TOXO II IgG versus the other techniques. The "negative" sub-group corresponded to 102 sera, negative with the DYE TEST, from pregnant women not immunised against toxoplasmosis. These 200 sera were tested in parallel with the LDBIO-TOXO II IgG, ELISA 1 IgG and ELISA 2 IgG techniques.

○ Group II - Seroconversions

This is a retrospective analysis of 17 serum sequences (101 samples) from patients having presented toxoplasmosis seroconversion during their pregnancy.

Each sequential series includes the last negative serum and then a series of 3 to 5 sera showing the appearance of specific IgMs and the synthesis of specific IgGs (*ELISA 2 IgG*).

○ Group III – Monitoring of uninfected children

This is a retrospective analysis of 74 samples corresponding to 20 sequences of post-natal monitoring of children born to mothers who presented toxoplasmosis seroconversion during pregnancy. Each sequence of 2 to 6 sera shows the decrease in the titre of transmitted maternal IgGs until the serology became negative with the *ELISA 2 IgG* technique (between 5 and 13 months).

○ Group IV – Monitoring of infected children

This is a retrospective analysis of 85 samples from post-natal monitoring of 30 children with a congenital infection. Serological monitoring was done by *ELISA 2 IgG*.

○ Group V – Sensitivity - Specificity (malaria and viral infections)

Study of 69 sera from patients suffering from malaria or viral infections (Table 1). These samples were tested by *ELISA 2 IgG*. (They were all negative in their IgM screening). All the negative as well as discordant results were tested with the DYE TEST.

Infectious agent (n = 69)	POSITIVE <i>ELISA 2 IgG</i> (n = 44)	NEGATIVE <i>ELISA 2 IgG</i> (n = 25)
EBV (n = 5)	0	5
VZV (n = 3)	2	1
CMV (n = 5)	2	3
HBV (n = 9)	8	1
HAV (n = 2)	0	2
HCV (n = 10)	8	2
HIV (n = 10)	6	4
MALARIA (n = 25)	18	7

**Table 1:** Different infections tested in the study

- RESULTS:

- Group I: DYE TEST

	DYE TEST	LDBIO TOXO II IgG	ELISA 1 IgG	ELISA 2 IgG
POSITIVE	98	97	61	93
NEGATIVE	102	103	114	107
EQUIVOCAL	-	-	25	-
SPECIFICITY	-	100%	100%	100%
SENSITIVITY	-	99%	85%	95%

**Table 2:** Correlation of the DYE TEST with the 3 techniques. (The ELISA 1 IgG technique presents an equivocal area).

- 4 negative ELISA 2 IgG sera are positive with the LDBIO-TOXO II IgG and Dye test
- 11 negative ELISA 1 IgG sera are positive with the LDBIO-TOXO II IgG and Dye test
- 25 ELISA 1 IgG sera are equivocal: 24 are positive with the LDBIO-TOXO II IgG and Dye Test, and 1 serum is negative with the LDBIO-TOXO II IgG and Dye test

- Group II: Seroconversions

		ELISA 2 IgG	
		POSITIVE	NEGATIVE
LDBIO TOXO II	POSITIVE	70	10
	NEGATIVE	0	21

**Table 3:** LDBIO-TOXO II IgG/ELISA 2 IgG correlation over 101 seroconversion sera.  $p = 0.0016$

For 8/17 seroconversions (47%), IgGs are screened earlier by LDBIO-TOXO II IgG.

- Groups III and IV: Monitoring of new-borns

		ELISA 2 IgG	
		POSITIVE	NEGATIVE
LDBIO TOXO II	POSITIVE	130	18
	NEGATIVE	0	11

**Table 4:** LDBIO-TOXO II IgG/TEST 2 IgG correlation over 159 sera from post-natal monitoring.  $p < 0.0001$

Uninfected children: 13 sera, corresponding to 10/20 post-natal monitoring cases (50%), negative with ELISA 2 IgG, remain positive with LDBIO-TOXO II IgG, which reveals the transmitted maternal antibodies while the ELISA 2 IgG technique no longer detects them.

Infected children: 5 sera corresponding to 3 children are discordant. One of them shows that the results of his/her serology are temporarily negative with ELISA 2 IgG. The LDBIO-TOXO II IgG test remains positive, confirming his/her infection. For the other 2 children, LDBIO TOXO II IgG test shows an earlier positive result than ELISA 2 IgG.

It is nevertheless impossible to confirm a neosynthesis of IgG as this test does not differentiate transmitted maternal antibodies from newly synthesised antibodies.

- Group V: Sensitivity and Specificity (malaria and viral infections)

		ELISA 2 IgG	
		POSITIVE	NEGATIVE
LDBIO TOXO II + DYE TEST	POSITIVE	42	2
	NEGATIVE	2	23

**Table 5:** *LDBIO-TOXO II IgG/DYE TEST/ELISA 2 IgG* correlation over 69 malaria or viral infection sera.

In this population, there is 100% agreement between *LDBIO-TOXO II IgG* and DYE TEST: these results confirm the specificity and sensitivity of the *LDBIO-TOXO II IgG* test.

The study reveals 4 discordant results with the *ELISA 2 IgG* technique, 2 false negatives (1 HIV and 1 *P. falciparum*) and 2 false positives (2 *P. falciparum*), emphasising the utility of a confirmation technique for all the results that are close to the threshold.

- CONCLUSION:

- Group I (DYE TEST):

The *LDBIO-TOXO II IgG/DYE TEST* correlation is excellent (99% sensitivity, 100% specificity). *LDBIO-TOXO II IgG* could confirm the immune status of patients presenting an equivocal result or a low antibody titre in the screening.

- Group II (seroconversions):

The sensitivity of *LDBIO-TOXO II IgG* is greater than that of *ELISA 2 IgG* ( $p = 0.0016$ ). *LDBIO-TOXO II IgG* could confirm a seroconversion earlier than *ELISA 2 IgG*.

- Group III and IV (monitoring of new-borns):

the sensitivity of *LDBIO-TOXO II IgG* is much greater than that of *ELISA 2 IgG* ( $p < 0.0001$ ). When monitoring children, *LDBIO-TOXO II IgG* could be used to confirm or set aside negative serology results.

Nevertheless, *LDBIO-TOXO II IgG* does not differentiate transmitted maternal antibodies from antibodies that are newly synthesised by the child.

- Group V (malaria and viral infect.):

The *LDBIO-TOXO II IgG/DYE TEST* correlation is excellent (100% sensitivity and 100% specificity). These results demonstrate the need to use a confirmation assay to test the samples presenting a result in the screening that is close to the threshold.

The excellent performance of the *LDBIO-TOXO II IgG* kit justifies its use in confirming the results obtained by the IgG screening techniques (results that are equivocal, slightly positive or that pose problems of interpretation).

- **Reproducibility:**

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation with respect to specific bands is excellent.

- **Interferences:**

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

## Trouble shooting

***"The bands are pale with little contrast"***: Certain sera with low concentrations of antibodies may give such results.

***"Shaded areas can be seen, more or less coloured, slightly diffuse"***: The strip was not totally submerged in one of the reagents and did not incubate correctly along its entire length. Stains can also be present where the sample was deposited if the tray was not shaken after dispensing.

***"The background noise is significant, making reading very difficult"***: The washes were insufficient or the last incubation was too long. Ensure good test performance techniques, respect wash times and ensure water quality. Reduce the time of the last incubation.

Exceptionally, certain sera may react in a non-specific manner. Then, the result of the immunoblot cannot be used.

This non-specific background noise may involve only part of the strip, making the results uninterpretable for that part only.

***"A precipitate appears in the solution during the last step of development"***: the substrate may in fact partially precipitate (black flakes) in the buffer at the end of development. This phenomenon does not alter the quality of the development which must be continued normally. The last wash with distilled water eliminates the possible solid particles present.

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