

# Interest of Toxoplasma ICT IgG-IgM as a diagnostic tool of Toxoplasma serology in cats.

Ana Patrícia Lopes<sup>1,2</sup>, Luís Cardoso<sup>1,2</sup>, Raphaël Piarroux<sup>3</sup> \*

<sup>1</sup> Department of Veterinary Sciences, School of Agrarian and Veterinary Sciences (ECAV), University of Trás-os-Montes e Alto Douro (UTAD), Vila Real, Portugal

<sup>2</sup> Animal and Veterinary Research Centre (CECAV), UTAD, Vila Real, Portugal

<sup>3</sup> LDBIO Diagnostics, Lyon, France

\* Auteur correspondant, contact : [rpiarroux@ldbiodiag.com](mailto:rpiarroux@ldbiodiag.com)

Conflict of interest: RP employee of LDBIO Diagnostics laboratories ; other authors: no conflict of interest

## Introduction

Cats, definitive hosts of *Toxoplasma gondii*, have an important role in the live cycle of the parasite. Serology can be used as a marker of environmental exposure in a feline population.

We studied the Toxoplasma ICT IgG-IgM<sup>®</sup> (LDBIO Diagnostics) as a diagnostic test of toxoplasmosis in cats.

## Material and Methods

Retrospective study on *F. catus* sera (N=139) previously tested for toxoplasmosis with 1 or 2 techniques<sup>1</sup>: Toxoscreen<sup>®</sup> (BioMérieux) + Toxo HAI<sup>®</sup> (Fumouze) (N=88) ; Toxoscreen<sup>®</sup> alone (N=51).

Use of Toxoplasma ICT IgG-IgM<sup>®</sup> according to instruction for use: 15µl serum, reading between 20 and 30 minutes (Fig 1).

Serum classification:

- If all techniques are in agreement, the result was kept.
- Discrepancies: confirmation by WB Toxo II LDBIO with A/G conjugate (Fig 2). Results of WB Toxo II was kept as serological diagnosis for all discrepancies.

Fig 1:  
Toxoplasma ICT in cats  
Positive result (Top)  
Negative result (Bottom)





Fig 2:  
WB LDBIO Toxo II in cats  
Negative results (strips 6-7)  
Positives results (strips 9-10)  
Human positive control (strip 8)



## Results

	ICT		HAI		Toxoscreen <sup>®</sup> (1/40 <sup>e</sup> )	
	Positive	Negative	Positive	Negative	Positive	Negative
<b>Positive serology</b> (N=33)	30	3	12	7	31	2
<b>Negative serology</b> (N=101)	0	101	0	65	5	96
<b>Performances</b>	Sensibility: 90,9% Spécificity: 100%		Sensibility: 63,1% Spécificity: 100%		Sensibility: 93,9% Spécificity: 95%	

Notes: 14 (10,1%) discrepant samples confirmed in WB Toxo II (10 negatives, 4 positives).

2 (1,4%) discrepant sample had insufficient volume for WB and were excluded from the study

3 (2,2%) invalid samples in ICT (no C band) were excluded from the study

## Conclusion

- Fast test, easy to use and interpret with high performances: LDBIO Toxoplasma ICT IgG-IGM<sup>®</sup> shows in this study its aptitude to diagnose veterinary toxoplasmosis, here in cats;
- Note: the use of ICT on whole blot has been validated in humans<sup>2</sup>. It remains to confirm/evaluate it in veterinary use. It allow the use of ICT as a real field test, without invasive procedures.

