TOXOPLASMA

CE0459

Western Blot IgG IgM

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Literatura

- Fekkar, A. *et al.* Comparison of immunoblotting, calculation of the Goldmann-Witmer coefficient, and real-time PCR using aqueous humor samples for diagnosis of ocular toxoplasmosis. *J. Clin. Microbiol.* **46**, 1965–1967 (2008).
- Garweg, J. G. Determinants of immunodiagnostic success in human ocular toxoplasmosis. *Parasite Immunol.* **27**, 61–68 (2005).
- Garweg, J. G., de Groot-Mijnes, J. D. F. & Montoya, J. G. Diagnostic Approach to Ocular Toxoplasmosis. *Ocular Immunology and Inflammation* **19**, 255–261 (2011).
- Garweg, J. G., Garweg, S.-D. L., Flueckiger, F., Jacquier, P. & Boehnke, M. Aqueous humor and serum immunoblotting for immunoglobulin types G, A, M, and E in cases of human ocular toxoplasmosis. *J. Clin. Microbiol.* **42**, 4593–4598 (2004).
- Goldmann, H. & Witmer, R. [Antibodies in the aqueous humor]. *Ophthalmologica* **127**, 323–330 (1954).
- L'ollivier, C. *et al.* Comparison of Mother and Child Antibodies That Target High-Molecular-Mass *Toxoplasma gondii* Antigens by Immunoblotting Improves Neonatal Diagnosis of Congenital Toxoplasmosis. *Clin. Vaccine Immunol.* **19**, 1326–1328 (2012).
- Maenz, M. *et al.* Ocular toxoplasmosis past, present and new aspects of an old disease. *Prog Retin Eye Res* **39**, 77–106 (2014).
- Magi, B. & Migliorini, L. Western blotting for the diagnosis of congenital toxoplasmosis. *New Microbiol.* **34**, 93–95 (2011).

- Pinon, J. M. *et al.* Strategy for diagnosis of congenital toxoplasmosis: evaluation of methods comparing mothers and newborns and standard methods for postnatal detection of immunoglobulin G, M, and A antibodies. *J. Clin. Microbiol.* **39**, 2267–2271 (2001).
- Potasman, I., Araujo, F. G. & Remington, J. S. Toxoplasma antigens recognized by naturally occurring human antibodies. *J. Clin. Microbiol.* **24**, 1050–1054 (1986).
- Remington, J. S., Thulliez, P. & Montoya, J. G. Recent developments for diagnosis of toxoplasmosis. *J. Clin. Microbiol.* **42**, 941–945 (2004).
- Rilling, V., Dietz, K., Krczal, D., Knotek, F. & Enders, G. Evaluation of a commercial IgG/IgM Western blot assay for early postnatal diagnosis of congenital toxoplasmosis. *Eur. J. Clin. Microbiol. Infect. Dis.* **22**, 174–180 (2003).
- Robert-Gangneux, F. *et al.* Usefulness of immunoblotting and Goldmann-Witmer coefficient for biological diagnosis of toxoplasmic retinochoroiditis. *Eur. J. Clin. Microbiol. Infect. Dis.* **23**, 34–38 (2004).
- Robert-Gangneux, F. & Darde, M.-L. Epidemiology of and Diagnostic Strategies for Toxoplasmosis. *Clinical Microbiology Reviews* **25**, 264–296 (2012).
- Runday, M. J., Ongkosuwito, J. V., Rothova, A. & Kijlstra, A. Intraocular anti-Toxoplasma gondii IgA antibody production in patients with ocular toxoplasmosis. *Am. J. Ophthalmol.* **127**, 294–300 (1999).
- Talabani, H. *et al.* Contributions of Immunoblotting, Real-Time PCR, and the Goldmann-Witmer Coefficient to Diagnosis of Atypical Toxoplasmic Retinochoroiditis. *Journal of Clinical Microbiology* **47**, 2131–2135 (2009).
- Tridapalli, E. *et al.* Congenital toxoplasmosis: the importance of the western blot method to avoid unnecessary therapy in potentially infected newborns. *Acta Paediatr.* **97**, 1298–1300 (2008).
- Turunen, H. J., Leinikki, P. O. & Saari, K. M. Demonstration of intraocular synthesis of immunoglobulin G toxoplasma antibodies for specific diagnosis of toxoplasmic chorioretinitis by enzyme immunoassay. *J. Clin. Microbiol.* **17**, 988–992 (1983).
- Villard, O. *et al.* Comparison of Enzyme-Linked Immunosorbent Assay, Immunoblotting, and PCR for Diagnosis of Toxoplasmic Chorioretinitis. *Journal of Clinical Microbiology* **41**, 3537–3541 (2003).
- Villard, O. *et al.* Serological diagnosis of Toxoplasma gondii infection: Recommendations from the French National Reference Center for Toxoplasmosis. *Diagn. Microbiol. Infect. Dis.* (2015).

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