ECHINOCOCCUS

Western Blot IgG

Test immunoblot de diagnostic *in vitro*
Technique manuelle / semi-automatisable

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ECHINOCOCCUS Western Blot IgG

INTENDED USE

ECHINOCOCCUS Western Blot (WB) IgG is a single use qualitative test of serological IgG diagnosis by Immunoblot Assay of alveolar echinococcosis and hydatidosis intended for confirmatory testing of a positive or equivocal result obtained through classic screening tests.

PRINCIPLE OF THE TEST

Western Blot technique

The antigens of *Echinococcus multilocularis* larvae, once separated by electrophoresis, are bound by electroblotting to the surface of a nitrocellulose membrane (called the transfer) cut into 24 strips numbered from 1 to 24.

Conduct of the test

Each specimen to be tested is separately incubated with a strip. The specific antibodies potentially present in the sample selectively bind themselves onto the antigens. The alkaline phosphatase-anti human IgG conjugate then binds itself to the bound antibodies. Finally, the immunocomplexes react with the substrate. The antigens recognized by the specific antibodies of type IgG present in the samples are revealed as purple transversal bands.

REAGENTS SUPPLIED

Default: package of 24 test (#ECH-WB24G)

Italic: package of 12 tests (#ECH-WB12G) – **Bold:** package of 96 tests (#ECH-WB96G)

ID	Qty	Description	Composition
R1	1	Folder(s) of 24 (<i>12, 4x24</i>) STRIPS: precut + coloured Standards. (Each folder and each transfer is identified by a unique serial number)	Sensitized nitrocellulose. Coloured Molecular Weight (kDa): Blue: 250, Blue: 150, Blue: 100, Pink: 75, Blue: 50, Green: 37, Pink: 25, Blue: 20, Blue: 15, Yellow: 10.
R2	1	Vial of 30 (<i>30, 125</i>) mL of SAMPLE BUFFER (Ready to use - pink solution).	Buffer + surfactant.
R3	1	Vial(s) of 30 (<i>30, 2x60</i>) mL of ANTI IgG CONJUGATE (Ready to use - blue solution).	Buffer + anti-human IgG polyclonal goat sera conjugated with Alkaline Phosphatase + NaN ₃ (<0.1%) + stabilisers.
R5	1	Vial of 30 (<i>30, 125</i>) mL of SUBSTRATE (Ready to use - opaque brown vial).	Buffer + NBT + BCIP + stabilisers.
R6	1	Vial of 60 (<i>60, 250</i>) mL of WASH CONCENTRATE 10X BUFFER (To be diluted 10 times in distilled water - colourless solution).	Buffer + surfactant.
R10	1	Tube of 200 (<i>200, 2x200</i>) µL of POSITIVE CONTROL SERUM (Ready to use - red cap).	Buffer + pool of human sera positive in <i>Echinococcus multilocularis</i> serology + NaN ₃ (<0.1%) + stabilisers.

R1: The letter before each strip number is specific to the parameter.

R2, R3, R5 and R6 are common to all kits and have a unique lot number depending only on the date of their production. **It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.**

R10 is calibrated in immunoblot according to a reference lot and is only dedicated to this technique.

R3, R10 (NaN₃): EUH 032 - Contact with acids liberates very toxic gas.

EUH 210 Safety data sheet available on request as well as on our website www.ldbiodiagnostics.com

ADDITIONAL MATERIAL REQUIRED BUT NOT PROVIDED

- One multi-channel polypropylene incubation trays for mini-blots (# WBPP- 08 or equivalent).
- One rocking platform for immunoblots, one vacuum system for liquids (the #WBPP-08 tubs that we supply can be emptied by simply turning them over).
- Tubes and material for drawing the samples, graduated cylinders, adapted containers. Automatic pipettes, micropipettes and disposable tips (volumes of 25 µL, 1.2 mL and 2 mL).
- Distilled or deionised water. Absorbent paper (e.g., Whatman filter paper), transparent adhesive tape.
- Gloves, tweezers to handle the strips, cutter or scalpel, flat transparent ruler.

Note: Our reagents can be used in an automated immunoblot processor. Care should be taken with possible chemical contaminations of our reagents if the processor is shared with reagents from another manufacturer (known example: contamination by the TWEEN 20), and bacterial contaminations. Reserve vials for the processor. After processing, do not place the remaining used reagents back into the original vials.

STORAGE AND STABILITY

Store between 2 and 8°C. The reagents from the kit are stable until the expiry date indicated on the outer box and the vial labels. Do not use contaminated or cloudy reagent. Wash buffer diluted to 1/10 is stable for 2 months at +2 to +8 °C and one week at room temperature.

PRECAUTIONS FOR USE

Safety

- For *in vitro* use only. For professional use only. Only for technically trained personnel. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- The positive control is a serum of human origin that has been inactivated for HIV 1 and 2, hepatitis B and hepatitis C viruses. However, it must be handled like a potentially infectious product.
- The substrate contains a mixture of NBT and BCIP, toxic on contact (skin and mucous membranes) and inhalation.
- The reagents contain sodium azide which can form explosive metallic salts with lead and copper. Rinse any spill with water.
- Dispose of waste (samples, tips, tubes, wash liquid, used reagent...) according to good practices used in the industry and current regulations in the country.
- Any serious incident must be the subject of a declaration to the manufacturer and the competent authority.

Precautions

- Read and interpret the results under direct white light.
- It is preferable to use all reagents from the same batch. If different batches are used, ensure traceability.
- Use the strips in numerical order. Do not mix strips from different serial numbers; use the transfers in succession. Establish a specific distribution plan before starting the test.
- Do not touch the strips with your fingers; use tweezers.
- The reagents must be mixed well before use, particularly the concentrated wash buffer.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-channel contamination. Watch for the formation of foam or bubbles in the pipette tips (bacterial contamination of reagent vials).
- Clean incubation trays only with distilled water (never use detergent or bleach).
- The omission of a sample or the distribution of an inadequate volume may render the test result negative or positive, regardless of its actual status.

SPECIMEN COLLECTION

Aseptically collect the samples in dry tubes. A minimum of 25 µL of serum.

Keep the samples at 2-8 °C until they are processed. If they need to be stored more than a week, freeze the samples at -20 ± 5 °C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

PREPARATION OF REAGENTS

Wash buffer: For 4 tests, in a clean bottle, dilute 10 mL of Wash Concentrate 10X (**R6**) in 90 mL of distilled or deionised water. Be careful to mix the diluted buffer well.

TEST PROCEDURE

Nota Bene: It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.

1. Prepare a distribution plan for the samples and C+ positive control (**R10**).

Only by using this control can the test be technically validated and identification made, for a given serial number, of the specific bands developed. A C+ strip cannot be used to interpret the results of strips from a blot of a different serial number.

2. Cut the required number of strips (R1) using a scalpel and a clean and dry flat transparent ruler, keeping the blue positioning line on the strips: hold the strips firmly in place with the ruler and cut them on the side of the stain (the numbers are visible through the ruler).
3. Distribute 1.2 mL of sample buffer (R2) in each channel according to the established plan.
4. Deposit, in their numerical order, the numbered strips in the channels: Let the strips rehydrate themselves at the surface of the buffer for approximately 2 minutes, with the number visible at the top, THEN gently shake the tray to totally immerse them in the buffer.
5. Distribute the samples and positive control(s): according to the distribution plan, at a rate of 25 µL per channel. Gently shake the tray after each dispense. Place the tray on a rocking platform.
Incubate for 90 minutes ± 5 minutes at 20-26 °C.
6. Wash step: Empty the contents of the channels with a Pasteur pipette or by turning the incubation tray over. Dispense 2 to 3 mL of diluted Wash Buffer in each channel. Incubate on the rocking platform for 3 min. Repeat 2 times, then empty the contents of the channels. Ensure that the strips don't turn during these steps.
7. Dispense 1.2 mL of anti IgG conjugate (R3) into each channel. Place the tray on the rocking platform.
Incubate for 60 minutes ± 5 minutes at 20-26 °C.
8. Wash step: repeat step 6.
9. Distribute 1.2 mL of NBT/BCIP substrate (R5) into each of the channels. Place on the rocking platform and protect from direct light. **Incubate for 60 minutes** ± 5 minutes at 20-26 °C.

Regardless of the parameter, monitor the development of the colour. The development can be stopped if the background colour of the strip darkens to the point where reading is difficult (the quality of the wash steps has a fundamental influence on the background coloration). Note that the strips will lighten as they dry.

10. Stop the reaction by aspirating substrate with a Pasteur pipette or by turning the incubation tub over and dispensing 2 mL of distilled or deionised water in the channels. Repeat this last washing step one more time.
11. Drying the strips: With the channels still water-filled, take the strips by the numbered end using the tweezers and deposit them, with the number visible, onto a Whatman absorbent paper. Let air dry. The colour of the strips will naturally lighten while drying. Interpretation must only be performed after drying is complete.
12. Storage: Transfer the strips onto a sheet of paper, which will be used to archive them. Align the blue positioning lines. Keeping them in place with the flat ruler, stick the top of the strips with transparent adhesive tape.

For a good interpretation, the strips must be ordered by transfer and in their numerical order, spaced at a maximum of a few millimetres apart. It is unreliable to compare strips that are spaced far apart (e.g., no.2 with no.15). **It is dangerous** (false results) to compare strips from different kits (strips with different serial numbers).

QUALITY CONTROL AND INTERPRETATION

The serum control (R10) provided with the kit must be systematically included in any immunoblot series. It shows the typical profile and allows for technical validation of the good conduct of the test (the bands must appear very clearly on the strip) and to calibrate precisely the position and aspect of the specific bands to allow interpretation of the results of the strips from the same transfer (same serial number).

Nota Bene: The positive control (R10) profile may vary according to the lot number of the reagents used. Corresponding images are available on our website www.ldbiodiagnostics.com as an example.

Description of the bands

- The reading area is located on the lower half of the strip, between 7 and 26-28 kDa. The 26-28 kDa band is so named because it can present itself in various aspects: a single narrow band (at 26 or 28 kDa), double band (26 and 28 kDa) or large band covering the entire area from 26 to 28 kDa.
- The extreme 7 and 26-28 kDa bands are used to diagnose the genus *Echinococcus* (see below: § Interpretation I).
- The intermediate bands, located between 7 and 26-28 kDa are used, when they are present, to diagnose the species *granulosus* or *multilocularis* (see below: § Interpretation II).


Fig. 1: Examples of positive and negative results

The profiles are given as examples. The strips are marked with the letter "D" specific to the parameter from batch "03023".

Interpretation

- Diagnosis of the genus:
 - Presence of the extreme 7 and/or 26-28 kDa bands
- Diagnostis of the species:
 - Profile **P1** or **P2**: *Echinococcus granulosus* (E.g)
 - Profile **P3**: *Echinococcus multilocularis* (E.m)
 - Profile **P4** or **P5**: *E. multilocularis* or *E. granulosus*

Interpretation I

Diagnosis of the genus *Echinococcus*:

Search for the presence of the 7 and/or 26-28 kDa bands for each of the samples tested with the calibration tools described above (these bands are typical and generally very easy to locate).

The presence of the extreme 7 and/or 26-28 kDa bands is required to interpret the test as being positive and to conclude that anti-*Echinococcus* IgG antibodies are present in the tested sample.

Interpretation II

Differential diagnosis of the species *E. granulosus* versus *E. multilocularis*:

This is done by searching for specific bands of one or the other of the other species in the intermediate area between 7 and 26 kDa.

- Bands common to both species: 12, 15, 20, 24 kDa
- Narrow bands only found with *E. multilocularis*: 16, 17, 18 kDa
- Band only found with *E. granulosus*: a large diffuse band at 17 kDa.

5 different profiles can be found.

- Profiles P1, P2 and P3 (found in 70% of cases) diagnose the species:

PROFILE P1: Isolated 7 kDa band only.	Echinococcus granulosus
PROFILE P2: 7 kDa band + large diffuse 17 kDa band. (NB: the 26-28 kDa band is very often also present.)	Echinococcus granulosus
PROFILE P3: 26-28 band + the narrow 16 and/or 18 kDa bands. (NB: most of the other 7, 12, 15, 17, 20 or 24 kDa bands are very often also present.)	Echinococcus multilocularis

- The last 2 profiles, P4 and P5 (found in 30% of cases), do not differentiate the 2 species *E. granulosus* and *E. multilocularis*.

PROFILE P4: Isolated 26-28 kDa band only.	NO intermediate band
PROFILE P5: Association of the 7 + 26-28 kDa bands	NO intermediate band

Note 1: The isolated presence of one or more of the intermediate bands (12, 15, 16, 17, 18, 20 or 24 kDa) cannot be considered to be specific. These bands are never found isolated in the case of echinococcosis but are always associated with the 7 kDa and/or 26-28 kDa bands.

Note 2: The bands above and, more rarely, below the 7-28 kDa area are very often present. They must not be used to interpret the assay.

Note 3: Exceptionally, the 16 kDa band appeared larger than normal in a patient infected by *E. multilocularis*. Be careful not to confuse this band with the large 17 kDa band that is specific to *E. granulosus*.

Note 4: The intermediate bands are less intense than the 7 and 26-28 kDa bands. Properly developing them often requires incubation in the substrate for 60 minutes. Do not interrupt it too soon.


Fig. 2: Additional examples of positive immunoblot samples that come from patients infected with *E. multilocularis* and *E. granulosus*.

The profiles are given as examples. The strips are marked with the letter "D" specific to the parameter from batch "03023".

These samples were specially chosen to be weak positives: all the *E.m* profiles are incomplete (with the exception of the first strip, No. 13).

It is interesting to note the opposition of the profiles that is usually found for each species:

E. multilocularis: The 26-28 kDa band often appears in the form of a double band, and it is the most intense.

E. granulosus: inversely, the most intense band is the 7 kDa band.

But this rule is not absolute (e.g., *E. m* band No. 24 - *E. g* band No. 20).

To validate the results, always compare the profile of the immunoblot of each sample with that of the R10 positive control. The aspect of the bands is important when interpreting the test.

LIMITATION OF USE

- The diagnosis of an infectious disease cannot be established on the basis of a single test result.

- Serological results must be interpreted according to available information (e.g., epidemiology, clinical, imaging, biology...) in order to establish a diagnosis. They should not be used to make a diagnosis based on their positivity alone.

PERFORMANCES (see literature references p19)

Sensitivity (Se)

A multicentre study, conducted in two independent specialised laboratories and covering 111 patient sera (50 cases of hydatidosis and 61 cases of alveolar echinococcosis identified with certainty), provided the following results:

	ECHINOCOCCUS WB IgG: obtained profiles					
	Neg	P1	P2	P3	P4	P5
Hydatidosis (n=50)	1	12	22	0	1	14
Alveolar Echinococcosis (n=61)	2	0	0	41	7	11
Total (n=111)	3	12	22	41	8	25

Table 1: Sensitivity of the test and profiles obtained

Sensitivity of the assay: **Se = 97.3% with respect to the genus *Echinococcus***
 Se = 98% with respect to the species *E. granulosus*
 Se = 96.7% with respect to the species *E. multilocularis*

Diagnosis of the species: *E. granulosus* versus *E. multilocularis*:

Table 1 above allows one to calculate an ability to discriminate between the two species of **67.6%** (P1 + P2 + P3 profiles).

Specificity (Sp) – Cross reactions

147 serum samples, corresponding to 147 patients, were tested with the **ECHINOCOCCUS WB IgG** kit by the two preceding laboratories.

Sera from patients suffering from the following were included: neuro-cysticercosis *Taenia solium* (42), *Schistosoma sp.*(42), *Fasciola hepatica* (10), *Loa loa* (6), *Trichinella spiralis* (6), *Toxocara canis* (6), *Strongyloides stercoralis* (4), *Entamoeba histolytica* (4), *Leishmania infantum* (4), *Plasmodium falciparum* (3) and the following autoimmune diseases: RF rheumatoid factor (8), ANA Anti-Nuclear Antibodies (12).

139 sera are negative, demonstrating **94.6% specificity** in this population.

The 8 cross-reactions were exclusively observed as part of:

- Cysticercosis: presence of an isolated 7 kDa band in 5/42 patients.
- Autoimmune diseases: presence of an isolated narrow band at 28 kDa in 1/8 patients (FR+) and 2/12 ANA+ patients.

Nota bene: Fasciolosis: the presence of an isolated very large band (25-30 kDa) was found in 4/10 tested patients, but it cannot be confused with the specific 26-28 band.

Conclusion

The correlation between WB Echinococcus and clinical status is excellent.

Sensitivity Se = 97.3% [CI95 91.7 - 99.3%]

Specificity Sp = 94.6% [CI95 89.2 - 97.4%]

In addition, the WB allows a differential diagnosis of positive samples with a very specific profile for *E. multilocularis* and *E. granulosus*.

E. multilocularis profile (P3 profile)

Sensitivity = 67.2% [CI95 53.9-78.4%] Specificity relative to *E. granulosus* = 100% [91.1 - 100%].

E. granulosus profile (profiles P1 and P2)

Sensitivity = 68% [CI95 53.2 - 80.1%] Specificity compared to *E. multilocularis* = 100% [92.6 - 100%]. Note: The P1 profile was however found in 5 cases (out of 42) of cysticercosis.

Confidence intervals are calculated according to the Wilson's method with continuity correction.

Reproducibility

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation with respect to specific bands is excellent.

Interferences

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

TROUBLESHOOTING

"The bands are pale with little contrast": Certain sera with low concentrations of antibodies may give such results.

"Shaded areas can be seen, more or less coloured, slightly diffuse": The strip was not totally submerged in one of the reagents and did not incubate correctly along its entire length. Stains can also be present where the sample was deposited if the tray was not shaken after dispensing.

"The background noise is significant, making reading very difficult": The washes were insufficient or the last incubation was too long. Ensure good test performance techniques, respect wash times and ensure water quality. Reduce the time of the last incubation. Exceptionally, certain sera may react in a non-specific manner. Then, the result of the immunoblot cannot be used.

This non-specific background noise may involve only part of the strip, making the results uninterpretable for that part only.

"A precipitate appears in the solution during the last step of development": the substrate may in fact partially precipitate (black flakes) in the buffer at the end of development. This phenomenon does not alter the quality of the development which must be continued normally. The last wash with distilled water eliminates the possible solid particles present.

BIBLIOGRAPHIE/BIBLIOGRAPHY

Atanasov G, Benckert C, Thelen A, Tappe D, Frosch M, Teichmann vD, Barth TFE, Wittekind C, Schubert S, et Jonas S. 2013. « Alveolar Echinococcosis-Spreading Disease Challenging Clinicians: A Case Report and Literature Review ». *World Journal of Gastroenterology: WJG* 19 (26): 4257-61. doi:10.3748/wjg.v19.i26.4257.

Auer H. 2006. « [Relevance of parasitological examinations for the clinical course, epidemiology and prevention of alveolar echinococcosis - experiences of more than two decades in Austria] ». *Wiener Klinische Wochenschrift* 118 (19-20 Suppl 3): 18-26. doi:10.1007/s00508-006-0673-3.

Bart JM, Piarroux M, Sako Y, Grenouillet F, Bresson-Hadni S, Piarroux R, et Ito A. 2007. « Comparison of several commercial serologic kits and Em18 serology for detection of human alveolar echinococcosis ». *Diagnostic microbiology and infectious disease* 59 (1): 93-95. doi:10.1016/j.diagmicrobio.2007.03.018.

Brunetti E, Kern P, Vuitton DA, et Writing Panel for the WHO-IWGE. 2010. « Expert Consensus for the Diagnosis and Treatment

- of Cystic and Alveolar Echinococcosis in Humans ». *Acta Tropica* 114 (1): 1-16. doi:10.1016/j.actatropica.2009.11.001.
- Furuya K, Kawanaka M, Yamano K, Sato N, et H Honma H. 2004. « [Laboratory evaluation of commercial immunoblot assay kit for serodiagnosis of Echinococcus infections using sera from patients with alveolar hydatidosis in Hokkaido] ». *Kansenshōgaku zasshi. The Journal of the Japanese Association for Infectious Diseases* 78 (4): 320-26.
- Liance M, Janin V, Bresson-Hadni S, Vuitton DA, Houin R, et Piarroux R. 2000. « Immunodiagnosis of Echinococcus infections: confirmatory testing and species differentiation by a new commercial Western Blot ». *Journal of clinical microbiology* 38 (10): 3718-21.
- Logar J, Soba B, et Kotar T. 2008. « Serological evidence for human cystic echinococcosis in Slovenia ». *BMC infectious diseases* 8: 63. doi:10.1186/1471-2334-8-63.
- Logar J, Soba B, Lejko-Zupanc T, et Kotar T. 2007. « Human alveolar echinococcosis in Slovenia ». *Clinical microbiology and infection: the official publication of the European Society of Clinical Microbiology and Infectious Diseases* 13 (5): 544-46. doi:10.1111/j.1469-0691.2007.01701.x.
- Makni F, Hachicha L, Mseddi F, Hammami H, Cheikhrouhou F, Sellami H, Sellami A, et al. 2007. « [Contribution of Western blotting to the diagnosis of hydatidosis] ». *Bulletin De La Société De Pathologie Exotique (1990)* 100 (3): 171-73.
- Otranto D, et Eberhard ML. 2011. « Zoonotic Helminths Affecting the Human Eye ». *Parasites & Vectors* 4: 41. doi:10.1186/1756-3305-4-41.
- Reiter-Owona I, Grüner B, Frosch M, Hoerauf A, Kern P, et Tappe D. 2009. « Serological confirmatory testing of alveolar and cystic echinococcosis in clinical practice: results of a comparative study with commercialized and in-house assays ». *Clinical laboratory* 55 (1-2): 41-48.
- Rinaldi F, Brunetti E, Neumayr A, Maestri M, Goblirsch S, et Tamarozzi F. 2014. « Cystic Echinococcosis of the Liver: A Primer for Hepatologists ». *World Journal of Hepatology* 6 (5): 293-305. doi:10.4254/wjh.v6.i5.293.
- Tamarozzi, F.; Longoni, S.S.; Vola, A.; Degani, M.; Tais, S.; Rizzi, E.; Prato, M.; Scarso, S.; Silva, R.; Brunetti, E.; et al. 2021. « Evaluation of Nine Commercial Serological Tests for the Diagnosis of Human Hepatic Cyst Echinococcosis and the Differential Diagnosis with Other Focal Liver Lesions: A Diagnostic Accuracy Study ». *Diagnostics*, 11, 167. <https://doi.org/10.3390/diagnostics11020167>
- Tappe D, Grüner B, Kern P, et Frosch M. 2008. « Evaluation of a commercial Echinococcus Western Blot assay for serological follow-up of patients with alveolar echinococcosis ». *Clinical and vaccine immunology: CVI* 15 (11): 1633-37. doi:10.1128/CVI.00272-08.
- Yamano K, Yagi K, Furuya K, Sawada Y, Honma H, et Sato N. 2005. « Active Alveolar Hydatidosis with Sero-Negativity for Antibody to the 18 kDa Antigen ». *Japanese Journal of Infectious Diseases* 58 (2): 122-24.
- Zait H, Achir I, Guerchani MK, et Hamrioui B. 2013. « [Epidemiological profile of 290 cases of human cystic echinococcosis diagnosed in the Mustapha University Hospital (Algiers) from 2006 to 2011] ». *Pathologie-Biologie* 61 (5): 193-98. doi:10.1016/j.patbio.2013.03.001.

UPDATE NOTIFICATION – Please read carefully

RELEASE DATE	VERSION	MODIFICATION SUMMARY
30/11/2022	Vs16	New address
07/12/2022	Vs17	R6 without NaN3. Strip identified by letter D. Possible use of reagents from different batches.

NF EN ISO 13485

24 Av. Joannes MASSET – 69009 LYON – FRANCE
 Tel : +33(0)4 7883 3487 – Fax : +33(0)4 7883 3430
www.ldbiodiagnostics.com – info@ldbiodiag.com