SCHISTO II CE

Western Blot IgG

Test immunoblot de diagnostic *in vitro*
Technique manuelle / semi-automatisable

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SCHISTO II Western Blot IgG

INTENDED USE

SCHISTO II Western Blot (WB) IgG is a single use qualitative test of serological IgG diagnosis by Immunoblot Assay of schistosomiasis intended for confirmatory testing of a positive or equivocal result obtained through classic screening tests.

PRINCIPLE OF THE TEST

Western Blot technique

The antigens (adult *Schistosoma mansoni* + *Schistosoma haematobium*), once separated by electrophoresis, are bound by electroblotting to the surface of a nitrocellulose membrane (called the transfer) cut into 24 strips numbered from 1 to 24.

Conduct of the test

Each specimen to be tested is separately incubated with a strip. The specific antibodies potentially present in the sample selectively bind themselves onto the antigens. The alkaline phosphatase-anti human IgG conjugate then binds itself to the bound antibodies. Finally, the immunocomplexes react with the substrate. The antigens recognized by the specific antibodies of type IgG present in the samples are revealed as purple transversal bands.

REAGENTS SUPPLIED

Default: package of 24 tests (#SCH II-WB24G)

Italic: package of 12 tests (#SCH II-WB12G) – Bold: package of 96 tests (#SCH II-WB96G)

ID	Qty	Description	Composition
R1	1	Folder(s) of 24 (<i>12</i> , 4x24) STRIPS: precut + coloured Standards. (Each folder and each transfer is identified by a unique serial number).	Sensitized nitrocellulose. Coloured Molecular Weight (kDa): Blue: 250, Blue: 150, Blue: 100, Pink: 75, Blue: 50, Green: 37, Pink: 25, Blue: 20, Blue: 15, Yellow: 10.
R2	1	Vial of 30 (<i>30</i> , 125) mL of SAMPLE BUFFER (Ready to use - pink solution).	Buffer + surfactant + NaN ₃ (<0.1%).
R3	1	Vial(s) of 30 (<i>30</i> , 2x60) mL of ANTI IgG CONJUGATE (Ready to use - blue solution).	Buffer + anti-human IgG polyclonal goat sera conjugated with Alkaline Phosphatase + NaN ₃ (<0.1%) + stabilisers.
R5	1	Vial of 30 (<i>30</i> , 125) mL of SUBSTRATE (Ready to use - opaque brown vial).	Buffer + NBT + BCIP + stabilisers.
R6	1	Vial of 60 (<i>60</i> , 250) mL of WASH CONCENTRATE 10X BUFFER (<u>To be diluted 10 times</u> in distilled water - colourless solution).	Buffer + surfactant.
R10	1	Tube of 200 (<i>200</i> , 2x200) µL of POSITIVE CONTROL SERUM (Ready to use - red cap).	Buffer + pool of human sera positive in <i>Schistosoma</i> serology + NaN ₃ (<0.1%) + stabilisers.

R1: The letter before each strip number is specific to the parameter.

R2, R3, R5 and R6 are common to all kits and have a unique lot number depending only on the date of their production. It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.

R10 is calibrated in immunoblot according to a reference lot and is only dedicated to this technique.

R3, R10 (NaN₃): EUH 032 - Contact with acids liberates very toxic gas.

EUH 210 Safety data sheet available on request as well as on our website www.ldbiodiagnostics.com

ADDITIONAL MATERIAL REQUIRED BUT NOT PROVIDED

- One multi-channel polypropylene incubation trays for mini-blot (WBPP-08 or equivalent).
- One rocking platform for immunoblots, one vacuum system for liquids (the LDBIO Diagnostics WBPP-08 tubs that we supply can be emptied by simply turning them over).
- Tubes and material for drawing the samples, graduated cylinders, adapted containers. Automatic pipettes, micropipettes and disposable tips (volumes of 25 µL, 1.2 mL and 2 mL).
- Distilled or deionised water. Absorbent paper (e.g., Whatman filter paper), transparent adhesive tape.
- Gloves, tweezers to handle the strips, cutter or scalpel, flat transparent ruler.

Note: Our reagents can be used in an automated immunoblot processor. Care should be taken with possible chemical contaminations of our reagents if the processor is shared with reagents from another manufacturer (known example: contamination by the TWEEN 20), and bacterial contaminations. Reserve vials for the processor. After processing, do not place the remaining used reagents back into the original vials.

STORAGE AND STABILITY

Store between 2 and 8°C. The reagents from the kit are stable until the expiry date indicated on the outer box and the vial labels. Do not use contaminated or cloudy reagent. Wash buffer diluted to 1/10 is stable for 2 months at +2 to +8 °C and one week at room temperature.

PRECAUTIONS FOR USE

Safety

- For *in vitro* use only. For professional use only. Only for technically trained personnel. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- The positive control is a serum of human origin that has been inactivated for HIV 1 and 2, hepatitis B and hepatitis C viruses. However, it must be handled like a potentially infectious product.
- The substrate contains a mixture of NBT and BCIP, toxic on contact (skin and mucous membranes) and inhalation.
- The reagents contain sodium azide which can form explosive metallic salts with lead and copper. Rinse any spill with water.
- Dispose of waste (samples, tips, tubes, wash liquid, used reagent...) according to good practices used in the industry and current regulations in the country.
- Any serious incident must be the subject of a declaration to the manufacturer and the competent authority.

Precautions

- Read and interpret the results under direct white light.
- It is preferable to use all reagents from the same batch. If different batches are used, ensure traceability.
- Use the strips in numerical order. Do not mix strips from different serial numbers; use the transfers in succession. Establish a specific distribution plan before starting the test.
- Do not touch the strips with your fingers; use tweezers.
- The reagents must be mixed well before use, particularly the concentrated wash buffer.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-channel contamination. Watch for the formation of foam or bubbles in the pipette tips (bacterial contamination of reagent vials).
- Clean incubation trays only with distilled water (never use detergent or bleach).
- The omission of a sample or the distribution of an inadequate volume may render the test result negative or positive, regardless of its actual status.

SPECIMEN COLLECTION

Aseptically collect the samples in dry tubes. A minimum of 25 µL of serum are required.

Keep the samples at 2-8 °C until they are processed. If they need to be stored more than a week, freeze the samples at -20 ± 5 °C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

PREPARATION OF REAGENTS

Wash buffer: For 4 tests, in a clean bottle, dilute 10 mL of Wash Concentrate 10X (R6) in 90 mL of distilled or deionised water. Be careful to mix the diluted buffer well.

TEST PROCEDURE

Nota Bene: It is recommended to perform multiparameter testing (see the LDBIO Diagnostics immunoblot range) to limit the number of vials opened and to ensure better quality control.

1. Prepare a distribution plan for the samples and C+ positive control (R10).

Only by using this control can the test be technically validated and identification made, for a given serial number, of the specific bands developed. A C+ strip cannot be used to interpret the results of strips from a blot of a different serial number.

2. Cut the required number of strips (R1) using a scalpel and a clean and dry flat transparent ruler, keeping the blue positioning line on the strips: hold the strips firmly in place with the ruler and cut them on the side of the stain (the numbers are visible through the ruler).
3. Distribute 1.2 mL of sample buffer (R2) in each channel according to the established plan.
4. Deposit, in their numerical order, the numbered strips in the channels: Let the strips rehydrate themselves at the surface of the buffer for approximately 2 minutes, with the number visible at the top, THEN gently shake the tray to totally immerse them in the buffer.
5. Distribute the samples and positive control(s): according to the distribution plan, at a rate of 25 μ L per channel. Gently shake the tray after each dispense. Place the tray on a rocking platform.
Incubate for 90 minutes \pm 5 minutes at 20-26 °C.
6. Wash step: Empty the contents of the channels with a Pasteur pipette or by turning the incubation tray over. Dispense 2 to 3 mL of diluted Wash Buffer in each channel. Incubate on the rocking platform for 3 minutes. Repeat 2 times, then empty the contents of the channels. Ensure that the strips don't turn during these steps.
7. Dispense 1.2 mL of anti IgG conjugate (R3) into each channel. Place the tray on the rocking platform.
Incubate for 60 minutes \pm 5 minutes at 20-26 °C.
8. Wash step: repeat step 6.
9. Distribute 1.2 mL of NBT/BCIP substrate (R5) into each of the channels. Place on the rocking platform and protect from direct light. **Incubate for 60 minutes \pm 5 minutes at 20-26 °C.**

Regardless of the parameter, monitor the development of the colour. The development can be stopped if the background colour of the strip darkens to the point where reading is difficult (the quality of the wash steps has a fundamental influence on the background coloration). Note that the strips will lighten as they dry.

10. Stop the reaction by aspirating substrate with a Pasteur pipette or by turning the incubation tub over and dispensing 2 mL of distilled or deionised water in the channels. Repeat this last washing step one more time.
11. Drying the strips: With the channels still water-filled, take the strips by the numbered end using the tweezers and deposit them, with the number visible, onto a Whatman absorbent paper. Let air dry. The colour of the strips will naturally lighten while drying. Interpretation must only be performed after drying is complete.
12. Storage: Transfer the strips onto a sheet of paper, which will be used to archive them. Align the blue positioning lines. Keeping them in place with the flat ruler, stick the top of the strips with transparent adhesive tape.

For a good interpretation, the strips must be ordered by transfer and in their numerical order, spaced at a maximum of a few millimetres apart. It is unreliable to compare strips that are spaced far apart (e.g., no.2 with no.15). **It is dangerous** (false results) to compare strips from different kits (strips with different serial numbers).

QUALITY CONTROL AND INTERPRETATION

The serum control (R10) provided with the kit must be systematically included in any immunoblot series. It shows the typical profile and allows for technical validation of the good conduct of the test (the bands must appear very clearly on the strip) and to calibrate precisely the position and aspect of the specific bands to allow interpretation of the results of the strips from the same transfer (same serial number).

Nota Bene: The positive control (R10) profile may vary according to the lot number of the reagents used. Corresponding images are available on our website www.ldbiodiagnostics.com as an example.

Description of the bands

A positive sample can present numerous bands located between 8 and 200k kilodaltons (kDa).

The reading area is located at the bottom of the strip, between **8 and 34 kDa**.

8 bands are most often present: P8, P9, P10, P12-13, P14-15, P18, P22-24 and P30-34 at the corresponding molecular weights (see photograph on **Fig. 1**).

The aspect of the bands can be variable. The P8, P9, P10 and P18 low-molecular-weight bands are usually narrow. The other bands can take the form of a single large band, a doublet of 2 narrower bands, or of 1 of the 2 component bands of the doublet.


Fig. 1: Examples of positive and negative results

The profiles are given as examples. The strips are marked with the letter "G" specific to the parameter from batch "06016".

Interpretation

The presence of one of the **P30-34** or **P22-24** bands is indicative of schistosomiasis.

- If it is isolated (an exceptional situation), the P30-34 band must present itself as a large band in order to be taken into account. (E.g., strip No. 4 ◀ above).
- Band P22-24 can have all the aspects: narrow, large, single or double.
- The bands that are most often found are indicated on the “C+” strip to the left of the figure 1. Numerous other bands can be present in the 8-22 kDa area.
- **The P1 and P2** profiles could be indicative of the species (see: serological schistosomiasis on p. 16).
- The **P28** band is frequent. It is **non-specific** to *Schistosoma*.

Important points:

- The 22-24 can sometimes present itself in the form of an isolated band at 22 or 24 kDa.
- The strip 13, 1 and 21 “cross-reactions” sera presented to the right of Fig. 1 correspond to cases of malaria. They were specially selected among the rare sera that, during the evaluation, presented non-specific bands in the 8-34 kDa reading area.

To validate the results, always compare the profile of the immunoblot of each sample with that of the R10 positive control. The aspect of the bands is important when interpreting the test.

LIMITATION OF USE

- The diagnosis of an infectious disease cannot be established on the basis of a single test result.
- Serological results must be interpreted according to available information (e.g., epidemiology, clinical, imaging, biology...) in order to establish a diagnosis. They should not be used to make a diagnosis based on their positivity alone.

PERFORMANCES (see literature references p. 17)

The study of **Schisto II WB IgG** performance was done on 548 different sera.

Sensitivity (Se)

184 patient sera suspected of schistosomiasis were tested according to the recommendations described in the kit insert.

Schistosomiasis was proven by a positive parasitic search (*S. haematobium* (60), *S. mansoni* (38), *S.h* + *S.m* (3) co-infection) and/or suggestive clinical data.

n = 184

Number of specific bands	1	2	3	4	5	6	7
Frequency	4%	15%	14%	15%	16%	19%	15%

Table 1: Number of specific bands present on the strip for a positive result: 95% of the immunoblots present a minimum of 2 bands.

n = 184

Nature of the specific bands (kDa)	P8	P10	P12	P15	P18	P22-24	P30-34
Frequency	37%	38%	64%	57%	52%	97%	89%

Table 2: Frequency of the presence of each of the specific bands observed on the immunoblots during our study on 184 positive samples.

n = 184	POSITIVE	NEGATIVE	Se
Reference WB	177	7	96.2%
WB SCH II	182	2	98.9%

Table 3: Sensitivity: Results compared between the new *Schisto II WB IgG* test and the previous *Schistosoma WB IgG kit* (= Reference WB).

Sensitivity = 98.9%

Differential diagnosis of the species

101 samples out of 184 corresponded to patients whose parasitological search revealed the presence of eggs in the urine, faeces and/or a rectal biopsy.

In this population, we often observed a difference in the immunological profile that appears to be related to the species responsible for the infection, *S. haematobium* or *S. mansoni*. These two types of profiles are clearly presented in the image on p. 16 (blue arrows: P1 vs P2 profiles).

n = 101	<i>S.m</i> eggs	<i>S.h</i> eggs	<i>S.m</i> + <i>S.h</i> eggs
P1 profile	9	53	0
P2 profile	27	3	2
equivocal	2	4	1

Table 4: Correlation between the parasitological search and serological diagnosis.

In this population, the immunological profile allowed a species diagnosis in 79% of cases. These data must be confirmed by more extensive studies before they can be used for a clinical diagnosis.

Note: The immunological profile cannot differentiate a *S. mansoni* infection from a *S. mansoni* + *S. haematobium* co-infection.

Specificity (Sp)

364 sera corresponding to 364 different patients were tested by following the indications presented in the kit's insert. These sera belonged to healthy patients (BD = 61), patients suffering from autoimmune pathologies, anti-nuclear antibodies (ANA = 21), rheumatoid factor (RF = 20), or various helminthiasis and other parasitic diseases: cysticercosis (53), hydatidosis (11), alveolar echinococcosis (10), fasciolosis (15), strongyloidiasis (9), toxocariasis (TXA = 41), trichinosis (TRI = 21), filariasis (FIL = 24), malaria (29), leishmaniasis (31) and amebiasis (18).

12 out of 364 samples show a characteristic "positive schistosoma" profile presenting between 2 and 7 very well-defined, specific bands. These results evidence a co-infection, confirmed by the reference WB.

6 samples show a weak cross-reaction: 4 samples present a narrow band at 24 kDa and 2 samples a pale but large band at 30-34 kDa.

Specificity calculation: if one considers 12 probable co-infections as actual positives, **Specificity = 98.3%**.

Note: Regardless of the quality of the samples, a non-specific, narrow, sometimes intense band is fairly often found to be present at 28 kDa.

Conclusion

The performances of the Schisto II kit compared to the reference WB are excellent. It allows a better identification of patients with *S. haematobium* infection than the reference.

Se = 98.9% [IC95 95.7 - 99.8%]

Sp = 98.3% [IC95 96.1 - 99.3%]

Confidence intervals are calculated according to Wilson's method with continuity correction.

Reproducibility

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation with respect to specific bands is excellent.

Interferences

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

TROUBLESHOOTING

"The bands are pale with little contrast": Certain sera with low concentrations of antibodies may give such results.

"Shaded areas can be seen, more or less coloured, slightly diffuse": The strip was not totally submerged in one of the reagents and did not incubate correctly along its entire length. Stains can also be present where the sample was deposited if the tray was not shaken after dispensing.

"The background noise is significant, making reading very difficult": The washes were insufficient or the last incubation was too long. Ensure good test performance techniques, respect wash times and ensure water quality. Reduce the time of the last incubation. Exceptionally, certain sera may react in a non-specific manner. Then, the result of the immunoblot cannot be used.

This non-specific background noise may involve only part of the strip, making the results uninterpretable for that part only.

"A precipitate appears in the solution during the last step of development": the substrate may in fact partially precipitate (black flakes) in the buffer at the end of development. This phenomenon does not alter the quality of the development which must be continued normally. The last wash with distilled water eliminates the possible solid particles present.

BIBLIOGRAPHIE/BIBLIOGRAPHY

- Guegan H, Fillaux J, Charpentier E, Robert-Gangneux F, Chauvin P, Guemas E, *et al.* 2019. « Real-time PCR for diagnosis of imported schistosomiasis ». *PLoS Negl Trop Dis* 13(9): e0007711. doi:10.1371/journal.pntd.0007711
- Bevilacqua N, Pane S, Vairo F, Nicastrì E, Paglia MG, Ame S, Sañé Schepisi M, *et al.* 2012. « Accuracy of Indirect Haemagglutination and Western Blot Assays for the Detection of Anti-Schistosoma Antibodies in Non-Severe Febrile Patients in Two Tanzanian Hospitals ». *Scandinavian Journal of Infectious Diseases* 44 (6): 453-58. doi:10.3109/00365548.2011.645505.
- Boissier J, Moné H, Mitta G, Bargues MD, Molyneux D, *et Mas-Coma S.* 2015. « Schistosomiasis Reaches Europe ». *The Lancet Infectious Diseases* 15 (7): 757-58. doi:10.1016/S1473-3099(15)00084-5.
- Brunet J, W. Pfaff A, Hansmann Y, Gregorowicz G, Pesson B, Abou-Bacar A, *et Candolfi E.* 2015. « An Unusual Case of Hematuria in a French Family Returning from Corsica ». *International Journal of Infectious Diseases: IJID: Official Publication of the International Society for Infectious Diseases* 31 (février): 59-60. doi:10.1016/j.ijid.2014.10.024.
- Cavalcanti M, Silva LF, Peralta R, Barreto M, *et Peralta JM.* 2013. « Schistosomiasis in Areas of Low Endemicity: A New Era in Diagnosis ». *Trends in Parasitology* 29 (2): 75-82. doi:10.1016/j.pt.2012.11.003.
- Colley D, Bustinduy A, Secor E, *et King CH.* 2014. « Human Schistosomiasis ». *Lancet* 383 (9936): 2253-64. doi:10.1016/S0140-6736(13)61949-2.
- De Laval F, Savini H, Biance-Valero E, *et Simon F.* 2014. « Human Schistosomiasis: An Emerging Threat for Europe ». *The Lancet* 384 (9948): 1094-95. doi:10.1016/S0140-6736(14)61669-X.
- ECDC Stockholm 2014: « Rapid risk assessment: Local transmission of *Schistosoma haematobium* in Corsica, France ».: European Centre for Disease Prevention and Control.
<http://ecdc.europa.eu/en/publications/Publications/schistosoma-haematobium-risk-assessment-France-Germany.pdf>
- Holtfreter MC, Moné H, Müller-Stöver I, Mouahid G, *et Richter J.* 2014. « *Schistosoma Haematobium* Infections Acquired in Corsica, France, August 2013 ». *Euro Surveill: Bulletin Européen Sur Les Maladies Transmissibles = European Communicable Disease Bulletin* 19 (22).
- Moné H, Holtfreter MC, Allienne JF, Mintsá-Nguéma R, Ibikounlé M, Boissier J, Berry A, Mitta G, Richter J, *et Mouahid G.* 2015. « Introgressive Hybridizations of *Schistosoma Haematobium* by *Schistosoma Bovis* at the Origin of the First Case Report of Schistosomiasis in Corsica (France, Europe) ». *Parasitology Research*, août. doi:10.1007/s00436-015-4643-4.
- Noormahomed EV, Nhacupe N, Mascaró-Lazcano C, Natane Mauaie M, Buene T, Abel Funzamo C, *et Benson C.* 2014. « A Cross-Sectional Serological Study of Cysticercosis, Schistosomiasis, Toxocariasis and Echinococcosis in HIV-1 Infected People in Beira, Mozambique ». Édité par Ana Flisser. *PLoS Neglected Tropical Diseases* 8 (9): e3121. doi:10.1371/journal.pntd.0003121.
- Sulahian A, Garin Y, Izri A, Verret C, Delaunay P, Van Gool P, *et Derouin F.* 2005. « Development and evaluation of a Western blot kit for diagnosis of schistosomiasis ». *Clinical and diagnostic laboratory immunology* 12 (4): 548-51. doi:10.1128/CDLI.12.4.548-551.2005.
- Wang W, Wang L, *et Liang YS.* 2012. « Susceptibility or Resistance of Praziquantel in Human Schistosomiasis: A Review ». *Parasitology Research* 111 (5): 1871-77. doi:10.1007/s00436-012-3151-z.

UPDATE NOTIFICATION – Please read carefully

RELEASE DATE	VERSION	MODIFICATION SUMMARY
12/08/2021	Vs 22	Removal of security warning R5 – Biblio - Contact email address – NaN3 EUH 032.
30/11/2022	Vs23	New address
21/12/2022	Vs24	R6 without NaN3. Strip identified by letter. Possible use of reagents from different batches.

NF EN ISO 13485

24 Av. Joannes MASSET – 69009 LYON – FRANCE
Tel : +33(0)4 7883 3487 – Fax : +33(0)4 7883 3430
www.ldbiodiagnostics.com – info@ldbiodiag.com