LDBIO TOXO II IgM CE0459

CONFIRMATION

Test immunoblot de diagnostic *in vitro*
Technique manuelle / semi-automatisable

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LDBIO TOXO II Western Blot IgM

INTENDED USE

LDBIO TOXO II IgM is a single use qualitative test of serological IgM diagnosis by Immunoblot Assay of toxoplasmosis intended for confirmatory testing of a positive or equivocal result obtained through classic screening tests.

PRINCIPLE OF THE TEST

Western Blot technique

The antigens of *Toxoplasma gondii*, once separated by electrophoresis, are bound by electroblotting to the surface of a nitrocellulose membrane (called the transfer) cut into 24 strips numbered from 1 to 24.

Conduct of the test

Each specimen to be tested is separately incubated with a strip. The specific antibodies potentially present in the sample selectively bind themselves onto the antigens. The alkaline phosphatase-anti human IgM conjugate then binds itself to the bound antibodies. Finally, the immunocomplexes react with the substrate. The antigens recognized by the specific antibodies of type IgM present in the samples are revealed as purple transversal bands.

REAGENTS SUPPLIED

Default: package of 24 tests (#T2M-24M)

Italic: package of 12 tests (#T2M-12M) – **Bold**: package of 96 tests (#T2M-96M)

ID	Qty	Description	Composition
R1	1	Folder(s) of 24 (<i>12</i> , 4x24) STRIPS: precut + coloured Standards. (Each folder and each transfer is identified by a unique serial number)	Sensitized nitrocellulose. Coloured Molecular Weight (kDa): Blue: 250, Blue: 150, Blue: 100, Pink: 75, Blue: 50, Green: 37, Pink: 25, Blue: 20, Blue: 15.
R2	1	Vial of 30 (<i>30</i> , 125) mL of SAMPLE BUFFER (Ready to use - pink solution).	Buffer + surfactant.
R4	1	Vial(s) of 30 (<i>30</i> , 2x60) mL of ANTI IgM CONJUGATE (Ready to use - yellow solution).	Buffer + anti-human IgM polyclonal goat sera conjugated with Alkaline Phosphatase + NaN3 (<0.1%) + stabilisers.
R5	1	Vial of 30 (<i>30</i> , 125) mL of SUBSTRATE (Ready to use - opaque brown vial).	Buffer + NBT + BCIP + stabilisers.
R6	1	Vial of 60 (<i>60</i> , 250) mL of WASH CONCENTRATE 10X BUFFER (To be diluted 10 times in distilled water - colourless solution).	Buffer + surfactant.
R10	1	Tube of 100 (<i>100</i> , 2x100) µL of POSITIVE CONTROL SERUM (Ready to use).	Buffer + pool of human sera positive in <i>Toxoplasma</i> serology + NaN3 (<0.1%) + stabilisers.

R1: The letter before each strip number is specific to the parameter.

R2, R4, R5 and R6 are common to all kits and have a unique lot number depending only on the date of their production. **It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.**

R10 is calibrated in immunoblot according to a reference lot and is only dedicated to this technique.

R3, R10 (NaN₃): EUH 032 - Contact with acids liberates very toxic gas.

EUH 210 Safety data sheet available on request as well as on our website www.ldbiodiagnostics.com

ADDITIONAL MATERIAL REQUIRED BUT NOT PROVIDED

- One multi-channel polypropylene incubation trays for mini-blots (#WBPP-08 or equivalent).
- One rocking platform for immunoblots, one vacuum system for liquids (the #WBPP-08 tubs that we supply can be emptied by simply turning them over).
- Tubes and material for drawing the samples, graduated cylinders, adapted containers. Automatic pipettes, micropipettes and disposable tips (volumes of 10 µL, 1.2 mL and 2 mL).
- Distilled or deionised water. Absorbent paper (e.g., Whatman filter paper), transparent adhesive tape.
- Gloves, tweezers to handle the strips, cutter or scalpel, flat transparent ruler.

Note: Our reagents can be used in an automated immunoblot processor. Care should be taken with possible chemical contaminations of our reagents if the processor is shared with reagents from another manufacturer (known example: contamination by the TWEEN 20), and bacterial contaminations. Reserve vials for the processor. After processing, do not place the remaining used reagents back into the original vials.

STORAGE AND STABILITY

Store between 2 and 8°C. The reagents from the kit are stable until the expiry date indicated on the outer box and the vial labels. Do not use contaminated or cloudy reagent. Wash buffer diluted to 1/10 is stable for 2 months at +2 to +8 °C and one week at room temperature.

PRECAUTIONS FOR USE

Safety

- For *in vitro* use only. For professional use only. Only for technically trained personnel. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- The positive control is a serum of human origin that has been inactivated for HIV 1 and 2, hepatitis B and hepatitis C viruses. However, it must be handled like a potentially infectious product.
- The substrate contains a mixture of NBT and BCIP, toxic on contact (skin and mucous membranes) and inhalation.
- The reagents contain sodium azide which can form explosive metallic salts with lead and copper. Rinse any spill with water.
- Dispose of waste (samples, tips, tubes, wash liquid, used reagent...) according to good practices used in the industry and current regulations in the country.
- Any serious incident must be the subject of a declaration to the manufacturer and the competent authority.

Precautions

- Read and interpret the results under direct white light.

- It is preferable to use all reagents from the same batch. If different batches are used, ensure traceability.
- Use the strips in numerical order. Do not mix strips from different serial numbers; use the transfers in succession. Establish a specific distribution plan before starting the test.
- Do not touch the strips with your fingers; use tweezers.
- The reagents must be mixed well before use, particularly the concentrated wash buffer.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-channel contamination. Watch for the formation of foam or bubbles in the pipette tips (bacterial contamination of reagent vials).
- Clean incubation trays only with distilled water (never use detergent or bleach).
- The omission of a sample or the distribution of an inadequate volume may render the test result negative or positive, regardless of its actual status.

SPECIMEN COLLECTION

Aseptically collect the samples in dry tubes. A minimum of 10 µL of serum is required.

Keep the samples at 2-8 °C until they are processed. If they need to be stored more than a week, freeze the samples at -20 ± 5 °C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

PREPARATION OF REAGENTS

Wash buffer: For 4 tests, in a clean bottle, dilute 10 mL of Wash Concentrate 10X (**R6**) in 90 mL of distilled or deionised water. Be careful to mix the diluted buffer well.

TEST PROCEDURE

Nota Bene: It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.

1. Prepare a distribution plan for the samples and C+ positive control (**R10**).

Only by using this control can the test be technically validated and identification made, for a given serial number, of the specific bands developed. A C+ strip cannot be used to interpret the results of strips from a blot of a different serial number.

2. Cut the required number of strips (R1) using a scalpel and a clean and dry flat transparent ruler, keeping the blue positioning line on the strips: hold the strips firmly in place with the ruler and cut them on the side of the strain (the numbers are visible through the ruler).
3. Distribute 1.2 mL of sample buffer (R2) in each channel according to the established plan.
4. Deposit, in their numerical order, the numbered strips in the channels: Let the strips rehydrate themselves at the surface of the buffer for approximately 2 minutes, with the number visible at the top, THEN gently shake the tray to totally immerse them in the buffer.
5. Distribute the samples and positive control(s): according to the distribution plan, at a rate of 10 µL per channel. Gently shake the tray after each dispense. Place the tray on a rocking platform.
Incubate for 90 minutes \pm 5 minutes at 20-26 °C.
6. Wash step: Empty the contents of the channels with a Pasteur pipette or by turning the incubation tray over. Dispense 2 to 3 mL of diluted Wash Buffer in each channel. Incubate on the rocking platform for 3 minutes. Repeat 2 times, then empty the contents of the channels. Ensure that the strips don't turn during these steps.

7. Dispense 1.2 mL of anti IgM conjugate (R4) into each channel. Place the tray on the rocking platform. **Incubate for 60 minutes** \pm 5 minutes at 20-26 °C.
8. Wash step: repeat step 6.
9. Distribute 1.2 mL of NBT/BCIP substrate (R5) into each of the channels. Place on the rocking platform and protect from direct light. **Incubate for 60 minutes** \pm 5 minutes at 20-26 °C.

Regardless of the parameter, monitor the development of the colour. The development can be stopped if the background colour of the strip darkens to the point where reading is difficult (the quality of the wash steps has a fundamental influence on the background coloration). Note that the strips will lighten as they dry.

10. Stop the reaction by aspirating substrate with a Pasteur pipette or by turning the incubation tub over and dispensing 2 mL of distilled or deionised water in the channels. Repeat this last washing step one more time.
11. Drying the strips: With the channels still water-filled, take the strips by the numbered end using the tweezers and deposit them, with the number visible, onto a Whatman absorbent paper. Let air dry. The colour of the strips will naturally lighten while drying. Interpretation must only be performed after drying is complete.
12. Storage: Transfer the strips onto a sheet of paper, which will be used to archive them. Align the blue positioning lines. Keeping them in place with the flat ruler, stick the top of the strips with transparent adhesive tape.

For a good interpretation, the strips must be ordered by transfer and in their numerical order, spaced at a maximum of a few millimetres apart. It is unreliable to compare strips that are spaced far apart (e.g., no.2 with no.15). **It is dangerous** (false results) to compare strips from different kits (strips with different serial numbers).

QUALITY CONTROL AND INTERPRETATION

The serum control (R10) provided with the kit must be systematically included in any immunoblot series. It shows the typical profile and allows for (1) technical validation of the good conduct of the test (the bands must appear very clearly on the strip) and (2) to calibrate precisely the position and aspect of the specific bands to allow interpretation of the results of the strips from the same transfer (same serial number).

Nota Bene: The positive control (R10) profile may vary according to the lot number of the reagents used and can in particular not present all the intermediate bands (P31, 33, 38), more rarely present in positive samples. Corresponding images are available on our website www.ldbiodiagnostics.com as an example.

Description of the bands

A positive sample can present numerous bands located between 15 and 200k kilodaltons (kDa). Search for the presence of specific bands in the 30-40 kDa area for each of the tested samples with the calibration tools described above. These bands, grouped and well isolated, are typical and can generally be found very easily.


Fig. 1: Examples of positive and negative results (Molecular weight: kDa)

The profiles are given as examples. The strips are marked with the letter "K" specific to the parameter from batch "50016".

Interpretation

The presence on the strip of a **minimum of 2 bands out of specific bands P30, P31, P33, P38 and P40, AND the inclusion of the band P30 kDa**, allows the assay to be interpreted as positive and to conclude that anti-*T. gondii* IgM antibodies are present in the tested sample.

Note Bene:

- **P30 et P40** are the most frequent bands in case of mild positive IgM serology result.
- Other bands (e.g. **P35**) may be observed. They are not taken into account when reading the test.

To validate the results, always compare the profile of the immunoblot of each sample with that of the R10 positive control. The aspect of the bands is important when interpreting the test.

LIMITATION OF USE

- The diagnosis of an infectious disease cannot be established on the basis of a single test result.
- Serological results must be interpreted according to available information (e.g., epidemiology, clinical, imaging, biology...) in order to establish a diagnosis. They should not be used to make a diagnosis based on their positivity alone.

PERFORMANCES (see literature references)

The evaluation was carried out through a multicenter study involving four reference laboratories, specialized in the diagnosis of toxoplasmosis.

The evaluation included two cohorts of pregnant women featuring a positive IgM test result: one was follow-up of seroconversions proven by the appearance of IgG (93 seroconversions / 229 samples), the other of women who had false positive IgM by at least one screening technique and whose non-specificity was then proven by iterative follow-up of the samples without IgG seroconversion (68 patients / 158 samples)

Each center used its own panel of IgM techniques allowing a comparison of the WB performance with 6 commercially available kits.

1. Study of the WB performance in confirmation of seroconversion:

Sensitivity: WB Toxo II IgM was positive for 91 of the 93 seroconversions, confirming *Toxoplasma* infection.

$$Se = 97.8\% (95CI [91.7-99.6\%])$$

Specificity: WB Toxo II IgM was negative for 61 of the 68 patients with false IgM confirming the absence of *Toxoplasma* infection.

$$Sp = 89.7\% (95CI [79.3\%-95.4\%]).$$

2. Comparative performance study of the WB Toxo II IgM, sample by sample:

All WB results (n=387) were compared with those obtained by the 6 other techniques used in the study in the 4 reference laboratories: IgM ELISAs or automates and ISAGA IgM.

Detailed results of this study have been published: "*Diagnostic Accuracy of LDBIO-Toxo II IgG and IgM Western Blot in Suspected Seroconversion in Pregnancy: A Multicentre Study*". Pathogens **2022**, 11(6), 665. doi: 10.3390/pathogens11060665 / Supplementary material File 1.

The WB shows equivalent sensitivity and higher specificity than all other techniques used.

The excellent performances of the LDBIO TOXO II IgM kit justify its use as a confirmation of results obtained by IgM screening techniques (equivocal results, weak positive results or results with difficulties of interpretation).

Reproducibility

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation with respect to specific bands is excellent.

Interferences

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

TROUBLESHOOTING

"The bands are pale with little contrast": Certain sera with low concentrations of antibodies may give such results.

"Shaded areas can be seen, more or less coloured, slightly diffuse": The strip was not totally submerged in one of the reagents and did not incubate correctly along its entire length. Stains can also be present where the sample was deposited if the tray was not shaken after dispensing.

"The background noise is significant, making reading very difficult": The washes were insufficient or the last incubation was too long. Ensure good test performance techniques, respect wash times and ensure water quality. Reduce the time of the last incubation. Exceptionally, certain sera may react in a non-specific manner. Then, the result of the immunoblot cannot be used.

This non-specific background noise may involve only part of the strip, making the results uninterpretable for that part only.

"A precipitate appears in the solution during the last step of development": the substrate may in fact partially precipitate (black flakes) in the buffer at the end of development. This phenomenon does not alter the quality of the development which must be continued normally. The last wash with distilled water eliminates the possible solid particles present.

BIBLIOGRAPHIE/BIBLIOGRAPHY

- Meroni V, Genco F, Scudeller L, Brenier-Pinchart MP, Fricker-Hidalgo H, L'Ollivier C, Paris L, Pelloux H. « Diagnostic Accuracy of LDBIO-Toxo II IgG and IgM Western Blot in Suspected Seroconversion in Pregnancy: A Multicentre Study ». *Pathogens* **2022**, 11(6), 665. doi: 10.3390/pathogens11060665
- Franck J, Garin Y, et Dumon H. « LDBio-Toxo II immunoglobulin G Western blot confirmatory test for anti-toxoplasma antibody detection ». *Journal of clinical microbiology* 46, n° 7 (juillet 2008): 2334-38. doi:10.1128/JCM.00182-08.
- Jost C, Touafek F, Fekkar A, Courtin R, Ribeiro M, Mazier D, et Paris L. « Utility of immunoblotting for early diagnosis of toxoplasmosis seroconversion in pregnant women ». *Clinical and vaccine immunology: CVI* 18, n° 11 (novembre 2011): 1908-12. doi:10.1128/CVI.05303-11.
- Robert-Gangneux F, et Darde ML. « Epidemiology of and Diagnostic Strategies for Toxoplasmosis ». *Clinical Microbiology Reviews* 25, n° 2 (1 avril 2012): 264-96. doi:10.1128/CMR.05013-11
- Villard O, Cimon B, L'Ollivier C, Fricker-Hidalgo H, Godineau N, Houze S, Paris L, Pelloux H, Villena I, et Candolfi E. « Serological Diagnosis of Toxoplasma Gondii Infection: Recommendations from the French National Reference Center for Toxoplasmosis ». *Diagnostic Microbiology and Infectious Disease*, 18 septembre 2015. doi:10.1016/j.diagmicrobio.2015.09.009
- Liesenfeld O, Press C, Montoya JG, Gill R, Isaac-Renton JL, Hedman K, Remington JS. « False-positive results in immunoglobulin M (IgM) toxoplasma antibody tests and importance of confirmatory testing: the Platelia Toxo IgM test ». *Journal of clinical microbiology*. 1997 Jan;35(1):174-8. doi: 10.1128/jcm.35.1.174-178.1997.
- Dhakal R, Gajurel K, Pomares C, Talucod J, Press CJ, Montoya JG. « Significance of a Positive Toxoplasma Immunoglobulin M Test Result in the United States ». *Journal of clinical microbiology*. 2015 Nov;53(11):3601-5. doi: 10.1128/JCM.01663-15.
- Petersen E, Borobio MV, Guy E, Liesenfeld O, Meroni V, Naessens A, Spranzi E, Thulliez P. « European multicenter study of the LIAISON automated diagnostic system for determination of Toxoplasma gondii-specific immunoglobulin G (IgG) and IgM and the IgG avidity index ». *Journal of clinical microbiology*. 2005 Apr;43(4):1570-4. doi: 10.1128/JCM.43.4.1570-1574.2005.

Roberts A, Hedman K, Luyasu V, Zufferey J, Bessières MH, Blatz RM, et al. « Multicenter evaluation of strategies for serodiagnosis of primary infection with *Toxoplasma gondii* ». European Journal of Clinical Microbiology & Infectious Diseases. 2001 Jul;20(7):467–74. doi: 10.1007/pl00011289.

Genco F, Lanzarini P, Chiaretto M, Prestla M & Meroni V. « Early diagnosis fo acute toxoplasmosis in IgG negative IgM positive pregnant women ». 25th European Congress of Clinical Microbiology & Infectious Diseases. Poster. 2015

Meroni V, Genco F, Corcione A ,Scudeller L, Brenier-Pinchart MP, Fricker-Hidalgo H, et al. « Diagnostic accuracy of toxoplasma western blot test in suspected seroconversion in pregnancy : a multicentric study ». International Congress On Congenital Toxoplasmosis. Poster. 2019.

UPDATE NOTIFICATION – Please read carefully

RELEASE DATE	VERSION	MODIFICATION SUMMARY
29/06/2022	Vs 02	Performance: reference to publication instead of table – update kits reference
30/11/2022	Vs03	New address
02/01/2023	Vs04	R6 without NaN3. Strip identified by letter. Possible use of reagents from different batches.
06/02/2023	Vs05	Indication of the band P35 (not specific)
20/09/2023	Vs06	Name and references in blue color

NF EN ISO 13485

24 Av. Joannes MASSET – 69009 LYON – FRANCE
 Tel : +33(0)4 7883 3487 – Fax : +33(0)4 7883 3430
www.ldbiodiagnostics.com – info@ldbiodiag.com