LEISHMANIA

Western Blot IgG

Test immunoblot de diagnostic *in vitro*
Technique manuelle / semi-automatisable

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LEISHMANIA Western Blot IgG

INTENDED USE

LEISHMANIA Western Blot (WB) IgG is a single use qualitative test of serological IgG diagnosis by Immunoblot Assay of leishmaniasis intended for confirmatory testing of a positive or equivocal result obtained through classic screening tests.

PRINCIPLE OF THE TEST

Western Blot technique

The antigens of *Leishmania infantum*, once separated by electrophoresis, are bound by electroblotting to the surface of a nitrocellulose membrane (called the transfer) cut into 24 strips numbered from 1 to 24.

Conduct of the test

Each specimen to be tested is separately incubated with a strip. The specific antibodies potentially present in the sample selectively bind themselves onto the antigens. The alkaline phosphatase-anti human IgG conjugate then binds itself to the bound antibodies. Finally, the immunocomplexes react with the substrate. The antigens recognized by the specific antibodies of type IgG present in the samples are revealed as purple transversal bands.

REAGENTS SUPPLIED

Default: package of 24 tests (#LES-WB24G)

Italic: package of 12 tests (#LES-WB12G) – **Bold:** package of 96 tests (#LES-WB96G)

ID	Qty	Description	Composition
R1	1	Folder(s) of 24 (12, 4x24) STRIPS: precut + coloured Standards. (Each folder and each transfer is identified by a unique serial number)	Sensitized nitrocellulose. Coloured Molecular Weight (kDa): Blue: 250, Blue: 150, Blue: 100, Pink: 75, Blue: 50, Green: 37, Pink: 25, Blue: 20, Blue: 15, Yellow: 10.
R2	1	Vial of 30 (30, 125) mL of SAMPLE BUFFER (Ready to use - pink solution).	Buffer + surfactant.
R3	1	Vial(s) of 30 (30, 2x60) mL of ANTI IgG CONJUGATE (Ready to use - blue solution).	Buffer + anti-human IgG polyclonal goat sera conjugated with Alkaline Phosphatase + NaN3 (<0.1%) + stabilisers.
R5	1	Vial of 30 (30, 125) mL of SUBSTRATE (Ready to use - opaque brown vial).	Buffer + NBT + BCIP + stabilisers.
R6	1	Vial of 60 (60, 250) mL of WASH CONCENTRATE 10X BUFFER (To be diluted 10 times in distilled water - colourless solution).	Buffer + surfactant + NaN3 (<0.1%).
R10	1	Tube of 200 (200, 2x200) µL of POSITIVE CONTROL SERUM (Ready to use - red cap).	Buffer + pool of human sera positive in <i>Leishmania</i> serology + NaN3 (<0.1%) + stabilisers.

R1: The letter before each strip number is specific to the parameter.

R2, R3, R5 and R6 are common to all kits and have a unique lot number depending only on the date of their production. **It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.**

R10 is calibrated in immunoblot according to a reference lot and is only dedicated to this technique.

R3, R10 (NaN₃): EUH 032 - Contact with acids liberates very toxic gas.

EUH 210 Safety data sheet available on request as well as on our website www.ldbiodiagnostics.com

ADDITIONAL MATERIAL REQUIRED BUT NOT PROVIDED

- One multi-channel polypropylene incubation trays for mini-blots (#WBPP-08 or equivalent).
- One rocking platform for immunoblots, one vacuum system for liquids (the #WBPP-08 tubs that we supply can be emptied by simply turning them over).
- Tubes and material for drawing the samples, graduated cylinders, adapted containers. Automatic pipettes, micropipettes and disposable tips (volumes of 25 µL, 1.2 mL and 2 mL).
- Distilled or deionised water. Absorbent paper (e.g., Whatman filter paper), transparent adhesive tape.
- Gloves, tweezers to handle the strips, cutter or scalpel, flat transparent ruler.

Note: Our reagents can be used in an automated immunoblot processor. Care should be taken with possible chemical contaminations of our reagents if the processor is shared with reagents from another manufacturer (known example: contamination by the TWEEN 20), and bacterial contaminations. Reserve vials for the processor. After processing, do not place the remaining used reagents back into the original vials.

STORAGE AND STABILITY

Store between 2 and 8°C. The reagents from the kit are stable until the expiry date indicated on the outer box and the vial labels. Do not use contaminated or cloudy reagent. Wash buffer diluted to 1/10 is stable for 2 months at +2 to +8 °C and one week at room temperature.

PRECAUTIONS FOR USE

Safety

- For *in vitro* use only. For professional use only. Only for technically trained personnel. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- The positive control is a serum of human origin that has been inactivated for HIV 1 and 2, hepatitis B and hepatitis C viruses. However, it must be handled like a potentially infectious product.
- The substrate contains a mixture of NBT and BCIP, toxic on contact (skin and mucous membranes) and inhalation.
- The reagents contain sodium azide which can form explosive metallic salts with lead and copper. Rinse any spill with water.
- Dispose of waste (samples, tips, tubes, wash liquid, used reagent...) according to good practices used in the industry and current regulations in the country.
- Any serious incident must be the subject of a declaration to the manufacturer and the competent authority.

Precautions

- Read and interpret the results under direct white light.
- It is preferable to use all reagents from the same batch. If different batches are used, ensure traceability.
- Use the strips in numerical order. Do not mix strips from different serial numbers; use the transfers in succession. Establish a specific distribution plan before starting the test.
- Do not touch the strips with your fingers; use tweezers.
- The reagents must be mixed well before use, particularly the concentrated wash buffer.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-channel contamination. Watch for the formation of foam or bubbles in the pipette tips (bacterial contamination of reagent vials).
- Clean incubation trays only with distilled water (never use detergent or bleach).
- The omission of a sample or the distribution of an inadequate volume may render the test result negative or positive, regardless of its actual status.

SPECIMEN COLLECTION

Aseptically collect the samples in dry tubes. A minimum of 25 µL of serum are required.

Keep the samples at 2-8 °C until they are processed. If they need to be stored more than a week, freeze the samples at -20 ± 5 °C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

PREPARATION OF REAGENTS

Wash buffer: For 4 tests, in a clean bottle, dilute 10 mL of Wash Concentrate 10X (**R6**) in 90 mL of distilled or deionised water. Be careful to mix the diluted buffer well.

TEST PROCEDURE

Nota Bene: It is recommended to perform multiparameter testing (see the LDBIO immunoblot range) to limit the number of vials opened and to ensure better quality control.

1. Prepare a distribution plan for the samples and C+ positive control (**R10**).

Only by using this control can the test be technically validated and identification made, for a given serial number, of the specific bands developed. A C+ strip cannot be used to interpret the results of strips from a blot of a different serial number.

2. Cut the required number of strips (R1) using a scalpel and a clean and dry flat transparent ruler, keeping the blue positioning line on the strips: hold the strips firmly in place with the ruler and cut them on the side of the strain (the numbers are visible through the ruler).
3. Distribute 1.2 mL of sample buffer (R2) in each channel according to the established plan.
4. Deposit, in their numerical order, the numbered strips in the channels: Let the strips rehydrate themselves at the surface of the buffer for approximately 2 minutes, with the number visible at the top, THEN gently shake the tray to totally immerse them in the buffer.
5. Distribute the samples and positive control(s): according to the distribution plan, at a rate of 25 µL per channel. Gently shake the tray after each dispense. Place the tray on a rocking platform.
Incubate for 90 minutes ± 5 minutes at 20-26 °C.
6. Wash step: Empty the contents of the channels with a Pasteur pipette or by turning the incubation tray over. Dispense 2 to 3 mL of diluted Wash Buffer in each channel. Incubate on the rocking platform for 3 minutes. Repeat 2 times, then empty the contents of the channels. Ensure that the strips don't turn during these steps.
7. Dispense 1.2 mL of anti IgG conjugate (R3) into each channel. Place the tray on the rocking platform.
Incubate for 60 minutes ± 5 minutes at 20-26 °C.
8. Wash step: repeat step 6.
9. Distribute 1.2 mL of NBT/BCIP substrate (R5) into each of the channels. Place on the rocking platform and protect from direct light. **Incubate for 60 minutes** ± 5 minutes at 20-26 °C.

Regardless of the parameter, monitor the development of the colour. The development can be stopped if the background colour of the strip darkens to the point where reading is difficult (the quality of the wash steps has a fundamental influence on the background coloration). Note that the strips will lighten as they dry.

10. Stop the reaction by aspirating substrate with a Pasteur pipette or by turning the incubation tub over and dispensing 2 mL of distilled or deionised water in the channels. Repeat this last washing step one more time.
11. Drying the strips: With the channels still water-filled, take the strips by the numbered end using the tweezers and deposit them, with the number visible, onto a Whatman absorbent paper. Let air dry. The colour of the strips will naturally lighten while drying. Interpretation must only be performed after drying is complete.
12. Storage: Transfer the strips onto a sheet of paper, which will be used to archive them. Align the blue positioning lines. Keeping them in place with the flat ruler, stick the top of the strips with transparent adhesive tape.

For a good interpretation, the strips must be ordered by transfer and in their numerical order, spaced at a maximum of a few millimetres apart. It is unreliable to compare strips that are spaced far apart (e.g., no.2 with no.15). **It is dangerous** (false results) to compare strips from different kits (strips with different serial numbers).

QUALITY CONTROL AND INTERPRETATION

The serum control (R10) provided with the kit must be systematically included in any immunoblot series. It shows the typical profile and allows for technical validation of the good conduct of the test (the bands must appear very clearly on the strip) and to calibrate precisely the position and aspect of the specific bands to allow interpretation of the results of the strips from the same transfer (same serial number).

Nota Bene: The positive control (R10) profile may vary according to the lot number of the reagents used. Corresponding images are available on our website www.ldbiodiagnostics.com as an example.

Description of the bands

A positive sample can present numerous bands located between 8 and 200 kilodaltons (kDa). Certain are specific to leishmaniasis. However, the difficulty in locating them in the middle of the other bands without a definite specificity is a disadvantage of using them.

Search for the presence of 14 and 16 kDa bands for each of the tested samples with the calibration tools described above. These bands, located at the bottom of the strip and generally well isolated, are generally very easy to interpret.


Fig. 1: Examples of positive and negative results

The profiles are given as examples. The strips are marked with the letter "C" specific to the parameter from batch "02007".

Interpretation

The presence on the strip of the 14 kDa and/or 16 kDa antigenic band allows the test to be interpreted as positive and to conclude that anti-*Leishmania* IgG antibodies are present in the tested sample.

To validate the results, always compare the profile of the immunoblot of each sample with that of the R10 positive control. The aspect of the bands is important when interpreting the test.

LIMITATION OF USE

- The diagnosis of an infectious disease cannot be established on the basis of a single test result.
- Serological results must be interpreted according to available information (e.g., epidemiology, clinical, imaging, biology...) in order to establish a diagnosis. They should not be used to make a diagnosis based on their positivity alone.

A negative serological result does not rule out the diagnosis of visceral leishmaniasis, particularly in immuno-suppressed patients. Any suspicion of leishmaniasis must automatically result in performing a parasitological search for protozoans.

PERFORMANCES (see literature references p15)

The **LEISHMANIA WB IgG** test was the subject of a comparative study with the IFA and ELISA techniques in an independent laboratory.

Sensitivity (Se)

	IFA	ELISA	WB
POSITIVE	41	40	51
NEGATIVE	10	11	0

Table 1: 51 sera from subjects suffering from progressive visceral leishmaniasis were tested by the 3 techniques. Unlike WB, ELISA and IFA presented false-negative results, in particular in immuno-compromised patients (HIV)

	IFA	ELISA	WB
POSITIVE	0	0	15
NEGATIVE	20	20	5

Table 2: 20 sera from healthy living patients in an endemic area and presenting a positive skin sensitivity test were tested in parallel with the three techniques: the sensitivity of IFA and ELISA is insufficient to detect very low levels of antibodies.

Specificity (Sp)

	IFA	ELISA	WB
POSITIVE	0	0	0
NEGATIVE	30	30	30

Table 3: 30 sera from healthy living adult patients in a non-endemic area were tested in parallel with the immunoblot, ELISA and IFA in an independent laboratory: the specificity of the three techniques was 100%.

Nota bene: False-positive reactions are frequently found, regardless of the technique, in patients suffering from trypanosomiasis (*Trypanosoma cruzi*).

Conclusion

The WB has an excellent sensitivity that allows it to effectively detect patients with visceral leishmaniasis even in an immunodepression context.

It allows the detection of asymptomatic carriers while having an excellent specificity, as demonstrated by the absence of positive sera in non-endemic patients.

Se = 100% [IC95 91.3 - 100%]

Sp = 100% [IC95 79.9 - 100%]

Confidence intervals are calculated according to Wilson's method with continuity correction.

Reproducibility

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation with respect to specific bands is excellent.

Interferences

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

TROUBLESHOOTING

"The bands are pale with little contrast": Certain sera with low concentrations of antibodies may give such results.

"Shaded areas can be seen, more or less coloured, slightly diffuse": The strip was not totally submerged in one of the reagents and did not incubate correctly along its entire length. Stains can also be present where the sample was deposited if the tray was not shaken after dispensing.

"The background noise is significant, making reading very difficult": The washes were insufficient or the last incubation was too long. Ensure good test performance techniques, respect wash times and ensure water quality. Reduce the time of the last incubation. Exceptionally, certain sera may react in a non-specific manner. Then, the result of the immunoblot cannot be used.

This non-specific background noise may involve only part of the strip, making the results uninterpretable for that part only.

"A precipitate appears in the solution during the last step of development": the substrate may in fact partially precipitate (black flakes) in the buffer at the end of development. This phenomenon does not alter the quality of the development which must be continued normally. The last wash with distilled water eliminates the possible solid particles present.

BIBLIOGRAPHIE/BIBLIOGRAPHY

- Aoun O, Mary C, Roqueplo C, Marié JL, Terrier O, Levieuge A, et Davoust B. 2009. « Canine Leishmaniasis in South-East of France: Screening of Leishmania Infantum Antibodies (western Blotting, ELISA) and Parasitaemia Levels by PCR Quantification ». *Veterinary Parasitology* 166 (1-2): 27-31. doi:10.1016/j.vetpar.2009.08.006.
- Biglino A, Bolla C, Concialdi E, Trisciuglio A, Romano A, et Ferroglio E. 2010. « Asymptomatic Leishmania Infantum Infection in an Area of Northwestern Italy (Piedmont Region) Where Such Infections Are Traditionally Nonendemic ». *Journal of Clinical Microbiology* 48 (1): 131-36. doi:10.1128/JCM.00416-09.
- Cota GF, De Sousa MR, Nogueira Demarqui F, et Rabello A. 2012. « The Diagnostic Accuracy of Serologic and Molecular Methods for Detecting Visceral Leishmaniasis in HIV Infected Patients: Meta-Analysis ». *PLoS Neglected Tropical Diseases* 6 (5): e1665. doi:10.1371/journal.pntd.0001665.
- Deniau M, Cañavate C, Faraut-Gambarelli F, et Marty P. 2003. « The Biological Diagnosis of Leishmaniasis in HIV-Infected Patients ». *Annals of Tropical Medicine & Parasitology* 97 (Supplement-1): 115-33. doi:10.1179/000349803225002598.
- Ferroglio E, Centaro E, Mignone W, et Trisciuglio A. 2007. « Evaluation of an ELISA Rapid Device for the Serological Diagnosis of Leishmania Infantum Infection in Dog as Compared with Immunofluorescence Assay and Western Blot ». *Veterinary Parasitology* 144 (1-2): 162-66. doi:10.1016/j.vetpar.2006.09.017.
- Kallel K, Ammari L, Kaouech E, Belhadj S, Anane S, Kilani B, et Chaker E. 2007. « [Asymptomatic bearing of Leishmania infantum among Tunisian HIV infected patients] ». *Pathologie-biologie* 55 (10): 521-24. doi:10.1016/j.patbio.2007.07.017.

- Lachaud L, Dedet JP, Marty P, Faraut F, Buffet P, Gangneux JP, Ravel J, Bastien P, et Working Group for the Notification of Human Leishmanioses in France. 2013. « Surveillance of Leishmaniasis in France, 1999 to 2012 ». *Euro Surveillance: Bulletin Européen Sur Les Maladies Transmissibles = European Communicable Disease Bulletin* 18 (29): 20534.
- Marty P, Lelievre A, Quaranta JF, Rahal A, Gari-Toussaint M, et Le Fichoux Y. 1994. « Use of the Leishmanin Skin Test and Western Blot Analysis for Epidemiological Studies in Visceral Leishmaniasis Areas: Experience in a Highly Endemic Focus in Alpes-Maritimes (France) ». *Transactions of the Royal Society of Tropical Medicine and Hygiene* 88 (6): 658-59.
- Marty P, Lelièvre A, Quaranta JF, Suffia I, Eulalio M, Gari-Toussaint M, Le Fichoux Y, et Kubar J. 1995. « Detection by Western Blot of Four Antigens Characterizing Acute Clinical Leishmaniasis due to *Leishmania Infantum* ». *Transactions of the Royal Society of Tropical Medicine and Hygiene* 89 (6): 690-91.
- Mary C, Lamouroux D, Dunan S, et Quilici M. 1992. « Western Blot Analysis of Antibodies to *Leishmania Infantum* Antigens: Potential of the 14-kD and 16-kD Antigens for Diagnosis and Epidemiologic Purposes ». *The American Journal of Tropical Medicine and Hygiene* 47 (6): 764-71.
- Pomares C, Despierres L, Del Giudice P, Delaunay P, Michel G, Ferrua B, et Marty P. 2012. « Western Blot Analysis as an Aid for the Diagnosis of Cutaneous Leishmaniasis due to *Leishmania Major* ». *Transactions of the Royal Society of Tropical Medicine and Hygiene* 106 (7): 452-54. doi:10.1016/j.trstmh.2012.03.001.
- Ready P. 2014. « Epidemiology of Visceral Leishmaniasis ». *Clinical Epidemiology*, mai, 147. doi:10.2147/CLEP.S44267.
- Saghrouni F, Khammari I, Kaabia N, Bouguila J, Ben Abdeljelil J, Fathallah A, Amri F, et Ben Saïd M. 2011. « Asymptomatic carriage of *Leishmania* in family members of patients with visceral leishmaniasis in Central Tunisia ». *Pathologie-biologie*, décembre. doi:10.1016/j.patbio.2011.11.001.
- Solano-Gallego L, Miró G, Koutinas A, Cardoso L, Pennisi MG, Ferrer L, Bourdeau P, Oliva G, Baneth G, et null The LeishVet Group. 2011. « LeishVet Guidelines for the Practical Management of Canine Leishmaniosis ». *Parasites & Vectors* 4: 86. doi:10.1186/1756-3305-4-86.
- Van Griensven J, Carrillo E, López-Vélez R, Lynen L, et Moreno J. 2014. « Leishmaniasis in Immunosuppressed Individuals ». *Clinical Microbiology and Infection* 20 (4): 286-99. doi:10.1111/1469-0691.12556.

UPDATE NOTIFICATION – Please read carefully

RELEASE DATE	VERSION	MODIFICATION SUMMARY
02/08/2021	Vs 15	Removal of security warning R5 - Contact email address – NaN3 EUH 032.
24/10/2022	Vs16	R6 without NaN3. Strip identified by letter C. Possible use of reagents from different batches.
30/11/2022	Vs17	New address

NF EN ISO 13485

24 Av. Joannes MASSET – 69009 LYON – FRANCE
 Tel : +33(0)4 7883 3487 – Fax : +33(0)4 7883 3430
www.ldbiodiagnostics.com – info@ldbiodiag.com