TOXOPLASMA CE0459

Western Blot IgA

Test immunoblot de diagnostic *in vitro*
Technique manuelle / semi-automatisable

In vitro diagnostic Immunoblot assay
Semi-automated / manual technique

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These Instructions for Use are **complementary** to the TOXOPLASMA WB IgG-IgM Kit (**#TOP-WB24G, #TOP-WB12G, #TOP-WB96G**) and anti-IgA conjugate (R8) (**#WB-IA30, #WB-IA60**) Instructions for Use.

The CIP-IgA must not be performed without first reading both products Instructions for Use.

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Find more information and translated IFU on our website at www.ldbiodiagnostics.com

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TOXOPLASMA Western Blot IgA

INTENDED USE: extension of application for the TOXOPLASMA WB IgG-IgM kit

The Comparison of Immunological Profiles (CIP-WB) IgA (CIP-IgA) is intended for the diagnosis of ocular toxoplasmosis by comparison of the immunological profile of the patient's serum and aqueous humor. It must be performed **in addition** to the CIP-IgG described in the **TOXOPLASMA WB IgG-IgM kit**.

PRINCIPLE OF THE TEST

Western Blot technique

The antigens of *Toxoplasma gondii*, once separated by electrophoresis, are bound by electroblotting to the surface of a nitrocellulose membrane (called the transfer) cut into 24 strips numbered from 1 to 24.

Conduct of the test

The CIP-IgA should be carried out in a complementary way to the CIP-IgG, as part of a search for ocular toxoplasmosis.

The test consists in incubating separately, **with 2 contiguous strips from the same transfer**, the two samples (sera, aqueous humor) whose immunological profiles are to be compared and then to reveal them by the successive action of the Alkaline Phosphatase – **anti-human IgA conjugate (R8)** and the substrate. The antigens recognized by the **specific IgA** class antibodies present in the samples are thus revealed in the form of purple transverse bands. The reading is made by successive comparison of the pairs of IgG and IgA strips, allowing to note the possible presence of strips revealed only by one of the samples and not by the other (Cf. § Interpretation).

REAGENTS SUPPLIED

Refer to Toxoplasma WB IgG-IgM and anti-IgA conjugate (R8) instructions for use.

ADDITIONAL MATERIAL REQUIRED BUT NOT PROVIDED

- CIP-IgA is performed using the **TOXOPLASMA WB IgG-IgM kit (reagents R1, R2, R5 and R6)** and the **LDBIO Diagnostics IgA conjugate (R8)** provided. Do not use any other test strips or reagents.
- Refer to the **TOXOPLASMA WB IgG-IgM Kit Instructions for Use**.

Note: Although our reagents can be used on an automated immunoblot processor, incubation times of CIP-IgA are different from those of other LDBIO assays and it is therefore recommended to perform CIP-IgA by manual technique.

STORAGE AND STABILITY

Refer to the **TOXOPLASMA WB IgG-IgM Kit** and **anti-IgA conjugate (R8)** Instructions for Use.

PRECAUTIONS FOR USE

Please refer to precautions for use of **TOXOPLASMA WB IgG-IgM**. The following instructions also applies.

- Do not mix 2 vials of IgA conjugate of different batch number.
- It is preferable to use all reagents from the same batch. If different batches are used, ensure traceability.

SPECIMEN COLLECTION

A minimum of 25µl of serum and aqueous humor is required. In the particular case of aqueous humor, the use of 50µl will increase the sensitivity of the test (Cf. § Procedure).

Keep the samples at 2-8 °C until they are processed. If they need to be stored more than a week, freeze the samples at -20 ± 5 °C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

TEST PROCEDURE

Note Bene: It is recommended to use CIP-IgA in addition to CIP IgG. If the volumes of aqueous humor do not allow it, it is recommended to prefer CIP-IgG to CIP-IgA.

Follow steps 1 to 4 of the procedure described in the Instructions for Use **TOXOPLASMA WB IgG-IgM**.

5. Distribute the samples according to the established distribution plan (step 1). **Use 25µl serum and 25-50µl (if possible) aqueous humor.** Gently shake the cell after each deposition. Place on an oscillating stirrer. Cover the incubation cell with parafilm or aluminum foil.

Incubate overnight (15h +/- 3h) at 20-26 °C.

Follow steps 6 to 12 described in the **TOXOPLASMA WB IgG-IgM** Instructions for Use.

- It is essential to simultaneously stop the color development of the 2 strips of a given pair for a given sub-class of antibodies, but one can independently stop that of the IgGs or IgAs (IgAs, in a lower concentration, usually develop more slowly than IgGs).
- Aqueous humor generally has a lower concentration of antibodies. The reaction must be allowed time to correctly develop and one should not be concerned to see the serum strips darken a little more.

Match, side-by-side, the IgG and IgA strips of each pair of samples in order of increasing numbers, by following the established distribution plan (step 1).

It is strictly mandatory to perform the comparison of a pair of samples using joint strips (contiguous numbers) from the same transfer (same serial number). It is unreliable to compare strips that are spaced far apart (e.g., no.2 with no.15). **It is dangerous** (false results) to compare strips from different kits (strips with different serial numbers).

QUALITY CONTROL AND INTERPRETATION

Description of the bands

A positive sample can present a significant number of bands located between 15 and 200 kDa. Only the bands with a molecular weight less than 120 kDa can be used to compare the profiles.

Interpretation

CIP-WB IgA (ocular toxoplasmosis)

Read the 2 contiguous strips simultaneously from top to bottom noting any antigenic strips present in the aqueous humor and absent from the serum.

Any well-defined resolution band, with a Molecular Weight (MW) less than 120 kDa and present only in the aqueous humor, indicates local synthesis of anti-toxoplasmosis antibodies, in favor of ocular toxoplasmosis.

Very important points

1. Results of CIP-WB IgA must be interpreted in conjunction with other clinical status, serological, parasitological, epidemiological, and medical imaging information in order to make the diagnosis of ocular toxoplasmosis.
2. A negative CIP-WB IgA result does not rule out the diagnosis of ocular toxoplasmosis. These patients should necessarily be followed up over time until the diagnosis of toxoplasmosis can be definitively confirmed or ruled out.
3. The bands can vary widely in appearance: thin, thick, more or less colored, intense... It is recommended when using this technique to carry out a few profile comparisons with known pairs of samples in order to become familiar with their reading. It is also recommended that the CIP-WB reading be initially performed independently by two people in the laboratory. In the event of a discrepancy in interpretation, a control IPC-WB should be performed.
4. The very high molecular weight (MW) antigenic fractions are very close together in the upper part of the strip in favor of better resolution of the medium and low MW fractions. The bands of MW > 120 kDa therefore cannot be used to interpret the assay: samples only presenting such profile differences cannot be returned as positive.

LIMITATION OF USE

- The diagnosis of an infectious disease cannot be established on the basis of a single test result.
- Serological results must be interpreted according to available information (e.g., epidemiology, clinical, imaging, biology...) in order to establish a diagnosis. They should not be used to make a diagnosis based on their positivity alone.

PERFORMANCES (see literature references)

Specificity - Sensitivity

Protocol

Performances of the CIP-IgA are based on two cohorts of patients from two reference French University Hospitals. One of the cohorts has been published [Mathis et al, 2018]. One study involved 87 patients: 42 ocular toxoplasmosis and 45 controls with ocular diseases. The other involved 59 patients: 24 cases of ocular toxoplasmosis and 35 controls.

The performance of the WB is given in terms of sensitivity and specificity. The 95% confidence intervals are calculated using the Wilson's method with continuity correction.

Results

	Cases	Controls	Total
WB+	21	0	21
WB-	45	80	125
Total	66	80	146

Table 1: Performances of CIP-WB IgA

WB IgA Sensibility = 31.8% [21.2-44.6%] WB IgA Specificity = 100% [94.3-100%]

Note that 7 sera (10.6% [4.7-21.2%]) were positive for IgA alone and would not have been detected by a CIP-IgG.

Conclusion

The comparative WB IgA profile shows excellent specificity in both case series (100% [94.3-100%]).

It has a low sensitivity (31.8% [21.2-44.6%]) but can be positive while WB IgG is negative in about 10% (10.6% [4.7-21.2%]) of cases of ocular toxoplasmosis. The complementary search for IgA, in addition to IgG, in the context of suspected ocular toxoplasmosis thus improves the sensitivity of the WB profile comparison to 75.7% (50/66 [63.4%-85.1%]).

Reproducibility

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation with respect to specific bands is excellent.

Interferences

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

TROUBLESHOOTING

Refer to **TOXOPLASMA WB IgG-IgM** Instructions for Use.

BIBLIOGRAPHIE/BIBLIOGRAPHY

- Fekkar, A. *et al.* Comparison of immunoblotting, calculation of the Goldmann-Witmer coefficient, and real-time PCR using aqueous humor samples for diagnosis of ocular toxoplasmosis. *J. Clin. Microbiol.* **46**, 1965–1967 (2008).
- Garweg, J. G. Determinants of immunodiagnostic success in human ocular toxoplasmosis. *Parasite Immunol.* **27**, 61–68 (2005).
- Garweg, J. G., de Groot-Mijnes, J. D. F. & Montoya, J. G. Diagnostic Approach to Ocular Toxoplasmosis. *Ocular Immunology and Inflammation* **19**, 255–261 (2011).
- Mathis, T. *et al.* Comparison of Immunoblotting (IgA and IgG) and the Goldmann-Witmer Coefficient for Diagnosis of Ocular Toxoplasmosis in Immunocompetent Patients. *Br J Ophthalmol.* 102 (10): 1454-58.
- Robert-Gangneux, F. *et al.* Usefulness of immunoblotting and Goldmann-Witmer coefficient for biological diagnosis of toxoplasmic retinochoroiditis. *Eur. J. Clin. Microbiol. Infect. Dis.* **23**, 34–38 (2004).
- Talabani, H. *et al.* Contributions of Immunoblotting, Real-Time PCR, and the Goldmann-Witmer Coefficient to Diagnosis of Atypical Toxoplasmic Retinochoroiditis. *Journal of Clinical Microbiology* **47**, 2131–2135 (2009).
- Villard, O. *et al.* Serological diagnosis of *Toxoplasma gondii* infection: Recommendations from the French National Reference Center for Toxoplasmosis. *Diagn. Microbiol. Infect. Dis.* (2015).

UPDATE NOTIFICATION – Please read carefully

RELEASE DATE	VERSION	MODIFICATION SUMMARY
09/08/2021	Vs 04	Night time incubation - Contact email address
28/09/2022	Vs 05	R8 indication – LDBIO address – Test procedure presentation – Use of reagents from different batches.

NF EN ISO 13485

24 Av. Joannes Masset – 69009 LYON – FRANCE
 Tel : +33(0)4 7883 3487 – Fax : +33(0)4 7883 3430
www.ldbiodiagnostics.com – info@ldbiodiag.com