PEO 

ICT IgG-IgM

Test immuno-chromatographique de diagnostic *in vitro*  
Technique manuelle

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# PEO ICT IgG-IgM

## Bird breeder's lung disease ICT IgG-IgM

### INTENDED USE

**PEO ICT IgG-IgM** is a rapid immunochromatographic serological diagnostic test proposed as a screening test for bird-breeder's lung disease (BBD), extrinsic allergic alveolitis (EAA) with avian antigens.

### PRINCIPLE OF THE TEST

**PEO ICT IgG-IgM** is a single use unitary test for a qualitative diagnostic. It is based on the principle of the homogeneous sandwich (immunological reaction of two same antigen epitopes with the two binding sites of a bivalent antibody).

Inside the cassette, the device is composed of:

- A nitrocellulose strip on which are spread two reactive bands: the antigens (squab crop milk) of the "test" band (T band) and the rabbit gamma globulins of the "control" band (C band),
- A fiberglass support (conjugate pad) which is impregnated of **black latex** particles coupled with squab crop milk antigens ("test" latex = T latex) and **blue latex** particles coupled with goat anti-rabbit IgG ("control" latex = C latex).

The test is run by successively dispensing the serum sample and an eluting solution (called the eluent) in the "sample well" of the cassette. Adding the eluent starts the concomitant migration (chromatography) of the serum and the latex particles. This migration is completed in 20-30 minutes.

If specific antibodies (IgG and/or IgM) are present in the sample, a complex is formed between the T latex and the patient's antibodies which is then captured by the T band. It results in the appearance of a **black line**: the test is positive.

The direct capture of the C latex by the C band results in the appearance of a blue line, meaning that the chromatography performed well. The appearance of this **blue line** is systematic and independent of the serological status of the patient.

Both letters "T" and "C" are printed on the cassette, materializing the position of the corresponding reading area.

### KIT COMPONENTS

ID	Description	Packaging	
		20 tests	100 tests
R1	Bag (sealed + zip) of 10 ready-to-use cassettes + a desiccant	2	10
R2**	Dropper bottle of 3 mL elution buffer	1	5
	Instructions for use	1	1

#### \*\* Eluent R2

- **Hazard pictograms**
- **Signal word:** Warning!

Code	Hazard statements
H317	May cause an allergic skin reaction
Code	Prevention
P261	Avoid breathing spray
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves and protective clothing
Intervention	
P305+P351+P338 P337+P313	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention
Disposal	
P501	Dispose of contents and container in accordance with relevant regulations

- **Hazard statements**
- **Precautionary statements**

EUH 210 Safety data sheet available on request as well as on our website [www.ldbiodiagnostics.com](http://www.ldbiodiagnostics.com).

## STORAGE AND STABILITY

- **Store the original sealed bag between 2 and 8°C.** Cassettes can be used until the expiry date written on the bag label. Do not freeze. Do not use after the expiration date.
- **The first opening of a bag of 10 tests must occur at least 15 minutes after** putting the bag at room temperature in order to avoid condensation in the bag.
- **Let the eluent at least 15 minutes at room temperature before use.**
- **After the first opening** of a bag, keep it **at room temperature (18-30°C)**, carefully closed (zip closure), the desiccant packet inside. After opening, the cassettes can be used **for up to 2 months**.
- The eluent is stable up to 2 months at room temperature (18-30°C) and until expiration date (as written on the kit) if kept between 2 and 8°C.

## PRECAUTIONS FOR USE

### Safety

- For *in vitro* use only. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- For professional use only. Only for technically trained personnel.
- All serum samples must be considered as potentially infectious and handled with care.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- Dispose of waste (samples, tips, tubes, cassettes, used reagent...) according to good practices used in the industry and current regulations in the country.
- Any serious incident must be declared to the manufacturer and the competent authority.

### Precautions

- Read and interpret the results under direct white light.
- Do not use eluent from another lot number.
- Do not use cassettes from two different lot numbers in the same run.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-cassette contamination.
- Do not use reagents after their expiration date.

- The omission of a sample or the distribution of an inadequate volume may render the test result positive or negative, regardless of its actual serological status.
- Only use cassettes carefully stored in their closed bag, the desiccant packet inside.

## SAMPLE COLLECTION AND PREPARATION

- The test can be done with serum.
- Sample collection must be sterile.
- On tubes with gel, do not collect gel: it could lead to false positives.
- Avoid hemolysis as much as possible.
- Keep the samples at 2-8 °C until they are processed. If they need to be stored, freeze the samples below -15°C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

## TEST PROCEDURE

Additional material required: micropipette and disposable tips for dispensing volumes of 15µL, timer.

**If the bags of 10 tests are kept at 2-8°C, let them at least 15 minutes at room temperature** before opening: the temperature of the bag must reach the room temperature to avoid condensation in the bag.

1. Take out the desired number of cassettes, then close carefully the bag with the zip closure (the desiccant packet inside) while pressing out as much air as possible. Close and store the bag at room temperature for up to 2 months.
2. **Identify each cassette** with the reference of each sample to be tested. Do not work with runs of more than 10 cassettes. Two successive runs of 10 cassettes must be separated by a few minutes in order to be able to make the reading at the indicated times (the use of use 2 timers is recommended).
3. Use a micropipette with a disposable tip to dispense 15µL of serum in the sample well. Do so for all cassettes before moving to the next step.
4. Dispense **4 drops** of the eluent of the kit. **Do not use the eluent of another lot number.** Keep the dropper vertically while dispensing. Close the dropper after use.
5. Start the timer when the eluent is dispensed in all the cassette of the run.

## READING AND INTERPRETATION

The reading must be done near a window or under direct light (example: a desk lamp). Avoid shadows on the reading area.

The reading must be done between 20 and 30 minutes after starting the timer.

***Do not take into account the results from readings after 30 minutes.***

- **Positive test:** 2 lines, a **black "T"** and a blue **"C"** appear in the corresponding areas. Every **"T"** line must be considered positive, even of very weak intensity. For very weak lines, make the reading with the eye vertically above the reading area.
- **Negative test:** No black line appears. Only the blue **"C"** line is visible.
- **Equivocal test:** In very rare cases, a faint, diffuse, grey line can appear on the **"T"** band. This result should be considered negative but controlled on another sample or technique.
- **Invalid test:** The **"C"** line does not appear. Read once again the instruction and repeat the test. If the problem persists, contact the manufacturer or your distributor.

Notes: This is a qualitative test. Intensity of the black line does not reflect the quantity of specific antibodies in the sample.

## QUALITY CONTROL

The blue “C” line allows the validation of the good running of the test. However, it is recommended to incorporate from times to times a known weak positive sample in a run.

## TEST LIMITATIONS

- Do not use too old serum sample with this technique. It is recommended to use samples frozen for less than 2 years.
- The positivity can be caused by the presence of IgG and/or IgM directed against the infectious agent, the test does not distinguish the type of antibodies present.
- The use of other body fluids (urine, CSF, saliva...) has not been validated.
- Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care. Very hemolytic samples can hide a weak positive test, due to important red background from hemolysis.
- Do not dispense 3 or 5 drops of eluent.
- **Do not use eluent outside the one given with the cassettes (same lot number).**
- The diagnosis of bird breeder’s disease cannot be established based on the results of a single test.
- Serological results must be interpreted according to available information (e.g., epidemiology, clinical, imaging, biology, etc.) in order to establish a diagnosis.

## PERFORMANCES

### Sensitivity, Specificity

The evaluation included two complementary cohorts of patients:

- A prospective cohort of 185 sera characterized by an automated screening technique in a specialized biology laboratory and whose positives were sent for confirmation to a reference laboratory for Bird-Breeder’s lung disease diagnosis (16 BBD, 5 sensitized, 164 negative).
- A complementary cohort of 63 characterized sera from the same reference laboratory (38 BBD, 10 sensitized, 15 negative).

The purpose of the evaluation was to determine the performance of PEO ICT IgG-IgM, in particular to detect BBD patients and sensitized patients, and the specificity of the test. Confidence intervals are calculated according to Wilson's method with correction for continuity.

<i>Patients</i>	PEO ICT IgG-IgM N=248	
	Positive	Negative
<i>BBD</i> <i>N=54</i>	48	6
<i>Sensitized</i> <i>N=15</i>	6	9
<i>Negative</i> <i>N=179</i>	26	153

LDBIO PEO ICT IgG-IgM was positive in 54/69 samples, 48/54 of which were BBD, while being negative in 153/179 negative samples. Its sensitivity is therefore 78.3% (IC95 [66.4-86.9%]) for all positives and 88.9% (IC95 [76.7-95.4%]) for BBD patients. Its specificity is 85.5% (CI95 [79.2-90.1%]).

These performances allow LDBIO ICT IgG IgM to be used as a screening test for Bird-Breeder's disease EAA. However, it cannot differentiate between sensitized and BBD patients and must therefore be confirmed by a confirmatory test, such as the PEO WB IgG test.

## Reproducibility

Inter-series and inter-lot reproducibility were tested. Correlation is excellent.

## Interferences

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

## BIBLIOGRAPHIE/BIBLIOGRAPHY

- Goudswaard J, Noordzij A, Stam JW. Pigeon IgA: a major antigen in pigeon breeder's disease. *Immunol Commun.* 1978;7(6):661–8.
- McSharry C, Anderson K, Boyd G. A review of antigen diversity causing lung disease among pigeon breeders. *Clin Exp Allergy J Br Soc Allergy Clin Immunol.* 2000 Sep;30(9):1221–9.
- Chan AL, Juarez MM, Leslie KO, Ismail HA, Albertson TE. Bird fancier's lung: a state-of-the-art review. *Clin Rev Allergy Immunol.* 2012 Aug;43(1-2):69–83.
- Nademi Z, Todryk S, Baldwin C. Characteristics of antibody responses in Pigeon Fanciers' Lung. *Mol Immunol.* 2013 Jun;54(2):227–32.
- Ohtani Y, Hisauchi K, Sumi Y, Miyashita Y, Sawada M, Miyake S, et al. Sequential changes in bronchoalveolar lavage cells and cytokines in a patient progressing from acute to chronic bird Fancier's lung disease. *Intern Med Tokyo Jpn.* 1999 Nov;38(11):896–9.
- Toubas D, Aubert D, Villena I, Foudrinier F, Chemla C, Pinon JM. Use of co-immunoelectrodiffusion to detect presumed disease-associated precipitating antibodies, and time-course value of specific isotypes in bird-breeder's disease. *J Immunol Methods.* 2003 Jan 15;272(1-2):135–45.
- Lacasse Y, Girard M, Cormier Y. Recent advances in hypersensitivity pneumonitis. *Chest.* 2012 Jul;142(1):208–17.

### UPDATE NOTIFICATION – Please read carefully

RELEASE DATE	VERSION	MODIFICATION SUMMARY
13/05/2022	Vs 01	Conception
30/11/2022	Vs 02	New address
23/07/2024	Vs03	Erratum : BBD versus AAE – serum use only

NF EN ISO 13485

24 Av. Joannes MASSET – 69009 LYON – FRANCE  
Tel : +33(0)4 7883 3487 – Fax : +33(0)4 7883 3430  
[www.ldbiodiagnostics.com](http://www.ldbiodiagnostics.com) – [info@ldbiodiag.com](mailto:info@ldbiodiag.com)