# ASPERGILLUS

CE

## ICT IgG-IgM

Test immuno-chromatographique de diagnostic *in vitro*  
Technique manuelle

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# ASPERGILLUS ICT IgG-IgM

## INTENDED USE

**ASPERGILLUS ICT IgG-IgM** is a rapid test based on the immuno-chromatography technology (lateral flow), allowing the *simultaneous detection* of both IgG and IgM class anti-*Aspergillus* antibodies in human sera.

## PRINCIPLE OF THE TEST

**ASPERGILLUS ICT IgG-IgM** is a single use unitary test for a qualitative diagnostic. It is based on the principle of the homogeneous sandwich (immunological reaction of two same antigen epitopes with the two binding sites of a bivalent antibody).

Inside the cassette, the device is composed of:

- a nitrocellulose strip on which are spread two reactive bands: the antigens (purified from an *Aspergillus fumigatus* culture) of the “test” band (T band) and the rabbit gamma globulins of the “control” band (C band),
- a fiberglass support (conjugate pad) which is impregnated of **black latex** particles coupled with *A. fumigatus* antigens (“test” latex = T latex) and **blue latex** particles coupled with goat anti-rabbit IgG (“control” latex = C latex).

The test is run by successively dispensing the serum sample and an eluting solution (called the eluent) in the “sample well” of the cassette. Adding the eluent starts the concomitant migration (chromatography) of the serum and the latex particles. This migration is completed in 20-30 minutes.

If specific antibodies (IgG and/or IgM) are present in the sample, a complex is formed between the T latex and the patient’s antibodies which is then captured by the T band. It results in the appearance of a **black line**: the test is positive.

The direct capture of the C latex by the C band results in the appearance of a blue line, meaning that the chromatography performed well. The appearance of this **blue line** is systematic and independent of the serological status of the patient.

Both letters “T” and “C” are printed on the cassette, materializing the position of the corresponding reading area.

## KIT COMPONENTS

ID	Description	Packaging	
		20 tests	100 tests
R1	Bag (sealed + zip) of 10 ready-to-use cassettes + a desiccant	2	10
R2**	Dropper bottle of 3 mL elution buffer	1	5
	Instructions for use	1	1

### \*\* Eluent R2

- **Hazard pictograms**
- **Signal word:** Warning!

○ **Hazard statements**

Code	Hazard statements
H302	Harmful if swallowed
H317	May cause an allergic skin reaction

○ **Precautionary statements**

Code	Prevention
P261	Avoid breathing spray
P264	wash hands thoroughly after handling
P270	Do not eat, drink or smoke when using the product
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves and protective clothing
	<b>Intervention</b>
P301+ P312+P330	IF SWALLOWED: rinse mouth. Do NOT induce vomiting. Call a POISON CENTER or doctor/physician if you feel unwell
P302+P352 P333+P313 P363	IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention Wash contaminated clothing before reuse
	<b>Disposal</b>
P501	Dispose of contents and container in accordance with relevant regulations

○ **Hazards not otherwise classified (HNOC)**

EUH 032 Contact with acids liberates very toxic gas.

EUH 210 Safety data sheet available on request as well as on our website [www.ldbiodiagnostics.com](http://www.ldbiodiagnostics.com).

## STORAGE AND STABILITY

- **Store the original sealed bag between 2 and 8°C.** Cassettes can be used until the expiry date written on the bag label. Do not freeze. Do not use after the expiration date.
- **The first opening of a bag of 10 tests must occur at least 15 minutes after** putting the bag at room temperature in order to avoid condensation in the bag.
- **Let the eluent remaining at least 15 minutes at room temperature before use.**
- **After the first opening** of a bag, keep it **at room temperature (18-30°C)**, carefully closed (zip closure), the desiccant packet inside. After opening, the cassettes can be used **for up to 2 months**.
- The eluent is stable up to 2 months at room temperature (18-30°C) and until expiration date (as written on the kit) if kept between 2 and 8°C.

## PRECAUTIONS FOR USE

### Safety

- For *in vitro* use only. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- For professional use only. Only for a technically trained technician.
- All serum samples must be considered as potentially infectious and handled with care.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- Dispose of waste (samples, tips, tubes, cassettes, used reagent...) according to good practices used in the industry and current regulations in the country.
- Any serious incident must be declared to the manufacturer and the competent authority.

## Precautions

- Read and interpret the results under direct white light.
- Do not use eluent from another lot number.
- Do not use cassettes from two different lot numbers in the same run.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-cassette contamination.
- Do not use reagents after their expiration date.
- The omission of a sample or the distribution of an inadequate volume may render the test result positive or negative, regardless of its actual serological status.
- Only use cassettes carefully stored in their closed bag, the desiccant packet inside.

## SERUM COLLECTION AND PREPARATION

- The test can be done with either serum or plasma.
- Sample collection must be sterile and can be done either on dry tube or with heparin, citrate or EDTA.
- On tubes with gel, do not collect gel: it could lead to false positives.
- Avoid hemolysis as much as possible.
- Keep the samples at 2-8 °C until they are processed. If they need to be stored, freeze the samples below -15°C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

## TEST PROCEDURE

Additional material required: micropipette and disposable tips for dispensing volumes of 15µL, timer.

**If the bags of 10 tests are kept at 2-8°C, let them at least 15 minutes at room temperature** before opening: the temperature of the bag must reach the room temperature to avoid condensation in the bag.

1. Take out the desired number of cassettes, then close carefully the bag with the zip closure (the desiccant packet inside) while pressing out as much air as possible. **Close and store the bag at room temperature for up to 2 months.**
2. **Identify each cassette** with the reference of each sample to be tested. Do not work with runs of more than 10 cassettes. Two successive runs of 10 cassettes must be separated by a few minutes in order to be able to make the reading at recommended times (use 2 timers).
3. Use a micropipette with a disposable tip to dispense 15µL of serum or plasma in the sample well. Do so for all cassettes before moving to the next step.
4. Dispense **4 drops** of the eluent of the kit. **Do not use the eluent of another lot number.** Keep the dropper vertically while dispensing. Close the dropper after use.
5. Start the timer when the eluent is dispensed in all the cassette of the run.

## READING AND INTERPRETATION

The reading must be done near a window or under direct light (example: a desk lamp). Avoid shadows on the reading area.

The reading must be done between 20 and 30 minutes after starting the timer.

**Do not take into account the results from readings after 30 minutes.**

- **Positive test:** 2 lines, a black "T" and a blue "C" appear in the corresponding areas. Every "T" line must be considered positive, even of very weak intensity. For very weak lines, make the reading with the eye vertically above the reading area.
- **Negative test:** No black line appears. Only the blue "C" line is visible.
- **Equivocal test:** In very rare cases, a faint, diffuse, grey line can appear on the "T" band. This result should be considered negative but controlled on another sample or technique.

- **Invalid test:** The “C” line does not appear. Read once again the instruction and repeat the test. If the problem persists, contact the manufacturer or your distributor.

Notes: This is a qualitative test. Intensity of the black line does not reflect the quantity of anti-*Aspergillus* antibodies in the sample.

Positivity of the test is proof of the contact of the patient with the infectious agent but doesn't prejudge the contact date or the clinical status of the patient.

## QUALITY CONTROL

- The blue “C” line allows the validation of the good running of the test.
- However, it is recommended to incorporate from times to times a known weak positive sample in a run.

Note: With some sera a second band black and attached to the blue band can be seen (=Duplication of the “C” band). The second band does not impair the result of the test which must only be read on the “T” line.

## TEST LIMITATIONS

- Do not use too old serum sample with this technique. It is recommended to use samples frozen for less than 2 years.
- The positivity can be caused by the presence of IgG and/or IgM directed against the infectious agent, the test does not distinguish the type of antibodies present.
- The use of other body fluids (urine, CSF, saliva, whole blood...) has not been validated.
- The use of hemolytic, icteric or lipidic samples is not recommended. However, no interference in the reaction was recorded. Still, very hemolytic samples can hide a weak positive test, due to important red background from hemolysis.
- Do not dispense 3 or 5 drops of eluent.
- **Do not use eluent outside the one given with the cassettes (same lot number).**
- The diagnosis of an infectious disease cannot be established based on the results of a single test.
- Serological results must be interpreted according to available information (e.g., epidemiology, clinical, imaging, biology, etc.) in order to establish a diagnosis.

## PERFORMANCES

### Sensitivity, Specificity

#### Material and methods

The evaluation of the kit was done prospectively in a reference hospital laboratory specialized in the diagnosis of *Aspergillus* related diseases. ASPERGILLUS ICT IgG-IgM was compared to two commercially available techniques used routinely in the laboratory: an ELISA technique (Orgentec®, France) and an ImmunoElectroPhoresis (IEP; Sebia®, France).

Each patient's file was studied and classified according to diagnostic consensus. The following disease categories were used: *Aspergillus* colonization, Allergic BronchoPulmonary Aspergillosis (ABPA), Chronic Pulmonary Aspergillosis (CPA) and all others proven *Aspergillus* diseases pooled together (acute or sub-acute aspergillosis, non-pulmonary aspergillosis). Patients not responding to any classification were classified as controls.

Sensitivity and specificity of each test was calculated, with their respective 95% Confidence Intervals (95CI) (Wilson's method).

#### Results

263 samples were gathered during the study: 44 cases and 219 controls.

The 44 cases were classified as follow: 6 ABPA, 23 colonization, 11 CPA, 4 others (acute and sub-acute aspergillosis).

Table 1 : Sensitivity

N = 44	ELISA	IEP	LDBIO ICT
Positive	20	35	40
Negative	15	9	4
Equivocal	9	0	0
Se [IC95]	57.1% [41-72%]	79.5% [66-89%]	90.9% [79-96%]

Table 2 : Specificity

N = 219	ELISA	IEP	LDBIO ICT
Positive	20	25	8
Negative	10	194	211
Equivocal	189	0	0
Sp [IC95]	92.6% [86-94%]	88.6% [84-92%]	96.3% [86-94%]

**Tables 1 and 2:** Results of the evaluation for all techniques on all serums (n=263). For Elisa, equivocal results were excluded from sensitivity/specificity calculations.

## Conclusion

Correlation between clinical status and LDBIO ICT was very good:

**Sensitivity Se = 90.9% [95CI: 79 - 96%]**

**Specificity Sp = 96.3% [95CI: 93 - 98%]**

## Reproducibility

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation is excellent.

## Interferences

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

## BIBLIOGRAPHIE/BIBLIOGRAPHY

Hunter E, Wilopo B, Richardson M, Kosmidis C, Denning D. 2021. « Effect of patient immunodeficiencies on the diagnostic performance of serological assays to detect *Aspergillus*-specific antibodies in chronic pulmonary aspergillosis ». *Respiratory Medicine* 178 (2021) 106290. doi:10.1016/j.rmed.2020.106290.

Kwizera R, Katende A, Teu A, Apolot D, Worodria W, Kirenga B, Bongomin F. 2019. « Algorithm-aided diagnosis of chronic pulmonary aspergillosis in low- and middle-income countries by use of a lateral flow device ». *European Journal of Clinical Microbiology & Infectious Diseases* (2020) 39:1–3. doi:10.1007/s10096-019-03782-x.

Thornton CR. 2019. « Detection of the ‘Big Five’ mold killers of humans: *Aspergillus*, *Fusarium*, *Lomentospora*, *Scedosporium* and *Mucormycetes* ». *Advances in Applied Microbiology*. doi: 10.1016/bs.aambs.2019.10.003.

Wilopo B, Richardson M, Denning D. 2019. « Diagnostic Aspects of Chronic Pulmonary Aspergillosis: Present and New Directions ». *Current Fungal Infection Reports* (2019) 13:292–300. doi:10.1007/s12281-019-00361-7.

Denning D, Riniotis K, Dobrashian R, et Sambatakou H. 2003. « Chronic Cavitary and Fibrosing Pulmonary and Pleural Aspergillosis: Case Series, Proposed Nomenclature Change, and Review ». *Clinical Infectious Diseases: An Official Publication of the Infectious Diseases Society of America* 37 Suppl 3 (octobre): S265-80. doi:10.1086/376526.

De Pauw B, Walsh TJ, Donnelly JP, Stevens DA, Edwards J, Calandra T, Pappas P, et al. 2008. « Revised Definitions of Invasive Fungal Disease from the European Organization for Research and Treatment of Cancer/Invasive Fungal Infections Cooperative Group and the National Institute of Allergy and Infectious Diseases Mycoses Study Group (EORTC/MSG) Consensus Group ». *Clinical Infectious Diseases: An Official Publication of the Infectious Diseases Society of America* 46 (12): 1813-21. doi:10.1086/588660.

Hohl T, et Feldmesser M. 2007. « *Aspergillus Fumigatus*: Principles of Pathogenesis and Host Defense ». *Eukaryotic Cell* 6 (11): 1953-63. doi:10.1128/EC.00274-07.

Oliva A, Flori P, Hennequin C, Dubus JC, Reynaud-Gaubert M, Charpin D, Vergnon JM, et al. 2015. « Evaluation of the *Aspergillus* Western Blot IgG Kit for Diagnosis of Chronic Aspergillosis ». *Journal of Clinical Microbiology* 53 (1): 248-54. doi:10.1128/JCM.02690-14.

Persat F, 2012. « [Aspergillus serology, from yesterday to today for tomorrow] ». *Journal De Mycologie Médicale* 22 (1): 72-82. doi:10.1016/j.mycmed.2012.01.004. <https://doi.org/10.3389/fcimb.2019.00012>

Piarroux R, Romain T, Martin A, Vainqueur D, Vitte J Lachaud L et al. 2019. « Multicenter Evaluation of a Novel Immunochromatographic Test for Anti-aspergillus IgG Detection ». *Frontiers in cellular and infection microbiology*. 2019;9:12. <https://doi.org/10.3389/fcimb.2019.00012>

Stucky Hunter E, Page I, Richardson M, Denning D. « Evaluation of the LDBio *Aspergillus* ICT lateral flow assay for serodiagnosis of allergic bronchopulmonary aspergillosis ». *Plos One*. 15(9):e0238855 <https://doi.org/10.1371/journal.pone.0238855>

Stucky Hunter E, Richardson M, Denning D. « Evaluation of LDBio *Aspergillus* ICT Lateral Flow Assay for IgG and IgM Antibody Detection in Chronic Pulmonary Aspergillosis ». *Journal of Clinical Microbiology* 2019;57(9):e00538–19. <https://doi.org/10.1128/JCM.00538-19>.

Zmeili O S, et Soubani AO. 2007. « Pulmonary Aspergillosis: A Clinical Update ». *QJM: Monthly Journal of the Association of Physicians* 100 (6): 317-34. doi:10.1093/qjmed/hcm035.

#### UPDATE NOTIFICATION – Please read carefully

RELEASE DATE	VERSION	MODIFICATION SUMMARY
15.10/2021	Vs 06	Biblio - addition EUH 032
20/11/2022	Vs07	New address
30/05/2023	Vs08	Toxicity update

NF EN ISO 13485

24 Av. Joannes MASSET – 69009 LYON – FRANCE  
 Tel : +33(0)4 7883 3487 – Fax : +33(0)4 7883 3430  
[www.ldbiodiagnostics.com](http://www.ldbiodiagnostics.com) – [info@ldbiodiag.com](mailto:info@ldbiodiag.com)