SCHISTOSOMA

CE

ICT IgG-IgM

Test immuno-chromatographique de diagnostic *in vitro*
Technique manuelle

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SCHISTOSOMA ICT IgG-IgM

INTENDED USE

SCHISTOSOMA ICT IgG-IgM is a rapid test based on the immuno-chromatography technology (lateral flow), allowing the *simultaneous detection* of both IgG and IgM class anti-*Schistosoma* antibodies in human sera.

PRINCIPLE OF THE TEST

SCHISTOSOMA ICT IgG-IgM is a single use unitary test for a qualitative diagnostic. It is based on the principle of the homogeneous sandwich (immunological reaction of two same antigen epitopes with the two binding sites of a bivalent antibody).

Inside the cassette, the device is composed of:

- a nitrocellulose strip on which are spread two reactive bands: the antigens (adult *Schistosoma mansoni*) of the “test” band (T band) and the rabbit gamma globulins of the “control” band (C band),
- a fiberglass support (conjugate pad) which is impregnated of **red latex** particles coupled with *S. mansoni* antigens (“test” latex = T latex) and **blue latex** particles coupled with goat anti-rabbit IgG (“control” latex = C latex).

The test is run by successively dispensing the serum sample and an eluting solution (called the eluent) in the “sample well” of the cassette. Adding the eluent starts the concomitant migration (chromatography) of the serum and the latex particles. This migration is completed in 20-30 minutes.

If specific antibodies (IgG and/or IgM) are present in the sample, a complex is formed between the T latex and the patient’s antibodies which is then captured by the T band. It results in the appearance of a **red line**: the test is positive.

The direct capture of the C latex by the C band results in the appearance of a blue line, meaning that the chromatography performed well. The appearance of this **blue line** is systematic and independent of the serological status of the patient.

Both letters “T” and “C” are printed on the cassette, materializing the position of the corresponding reading area.

KIT COMPONENTS

ID	Description	Packaging	
		20 tests	100 tests
R1	Bag (sealed + zip) of 10 ready-to-use cassettes + a desiccant	2	10
R2**	Dropper bottle of 2 mL elution buffer	1	5
	Instructions for use	1	1

** Eluent R2

- **Hazard pictograms**
- **Signal word:** Warning!

Code	Hazard statements
H317	May cause an allergic skin reaction
Code	Prevention
P261	Avoid breathing spray
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves and protective clothing
Intervention	
P305+P351+P338 P337+P313	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention
Disposal	
P501	Dispose of contents and container in accordance with relevant regulations

EUH 210 Safety data sheet available on request as well as on our website www.ldbiodiagnostics.com.

STORAGE AND STABILITY

- **Store the original sealed bag between 2 and 8°C.** Cassettes can be used until the expiry date written on the bag label. Do not freeze. Do not use after the expiration date.
- **The first opening of a bag of 10 tests must occur at least 15 minutes after** putting the bag at room temperature in order to avoid condensation in the bag.
- **Let the eluent remaining at least 15 minutes at room temperature before use.**
- **After the first opening** of a bag, keep it **at room temperature (18-30°C)**, carefully closed (zip closure), the desiccant packet inside. After opening, the cassettes can be used **for up to 2 months**.
- The eluent is stable up to 2 months at room temperature (18-30°C) and until expiration date (as written on the kit) if kept between 2 and 8°C.

PRECAUTIONS FOR USE

Safety

- For *in vitro* use only. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- For professional use only. Only for a technically trained technician.
- All serum samples must be considered as potentially infectious and handled with care.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- Dispose of waste (samples, tips, tubes, cassettes, used reagent...) according to good practices used in the industry and current regulations in the country.
- Any serious incident must be declared to the manufacturer and the competent authority.

Precautions

- Read and interpret the results under direct white light.
- Do not use eluent from another lot number.
- Do not use cassettes from two different lot numbers in the same run.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-cassette contamination.
- Do not use reagents after their expiration date.
- The omission of a sample or the distribution of an inadequate volume may render the test result positive or negative, regardless of its actual serological status.
- Only use cassettes carefully stored in their closed bag, the desiccant packet inside.

SERUM COLLECTION AND PREPARATION

- The test can be done with either serum or plasma.
- Sample collection must be sterile and can be done either on dry tube or with heparin (**do not use plasma from citrate or EDTA samples**).
- On tubes with gel, do not collect gel: it could lead to false positives.
- Avoid hemolysis as much as possible.
- Keep the samples at 2-8 °C until they are processed. If they need to be stored, freeze the samples below -15°C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.

TEST PROCEDURE

Additional material required: micropipette and disposable tips for dispensing volumes of 30µL, timer.

If the bags of 10 tests are kept at 2-8°C, let them at least 15 minutes at room temperature before opening: the temperature of the bag must reach the room temperature to avoid condensation in the bag.

1. Take out the desired number of cassettes, then close carefully the bag with the zip closure (the desiccant packet inside) while pressing out as much air as possible. **Close and store the bag at room temperature for up to 2 months.**
2. **Identify each cassette** with the reference of each sample to be tested. Do not work with runs of more than 10 cassettes. Two successive runs of 10 cassettes must be separated by a few minutes in order to be able to make the reading at recommended times (use 2 timers).
3. Use a micropipette with a disposable tip to dispense 30µL of serum or plasma in the sample well. Do so for all cassettes before moving to the next step.
4. Dispense **3 drops** of the eluent of the kit. **Do not use the eluent of another lot number.** Keep the dropper vertically while dispensing. Close the dropper after use.
5. Start the timer when the eluent is dispensed in all the cassette of the run.

READING AND INTERPRETATION

The reading must be done near a window or under direct light (example: a desk lamp). Avoid shadows on the reading area.

The reading must be done between 20 and 30 minutes after starting the timer.

Do not take into account the results from readings after 30 minutes.

- **Positive test:** 2 lines, a red "T" and a blue "C" appear in the corresponding areas. Every "T" line must be considered positive, even of very weak intensity. For very weak lines, make the reading with the eye vertically above the reading area.
- **Negative test:** No red line appears. Only the blue "C" line is visible.
- **Equivocal test:** In very rare cases, a faint, diffuse, grey line can appear on the "T" band. This result should be considered negative but controlled on another sample or technique.
- **Invalid test:** The "C" line does not appear. Read once again the instruction and repeat the test. If the problem persists, contact the manufacturer or your distributor.

Notes: This is a qualitative test. Intensity of the red line does not reflect the quantity of specific antibodies in the sample.

Positivity of the test is proof of the contact of the patient with the infectious agent but doesn't prejudice the contact date or the clinical status of the patient.

QUALITY CONTROL

The blue “C” line allows the validation of the good running of the test. However, it is recommended to incorporate from times to times a known weak positive sample in a run.

TEST LIMITATIONS

- Do not use too old serum sample with this technique. It is recommended to use samples frozen for less than 2 years.
- The positivity can be caused by the presence of IgG and/or IgM directed against the infectious agent, the test does not distinguish the type of antibodies present.
- The use of other body fluids (urine, CSF, saliva, whole blood...) has not been validated.
- The use of hemolytic, icteric or lipidic samples is not recommended (§ interferences). Very hemolytic samples can hide a weak positive test, due to important red background from hemolysis.
- Do not dispense 2 or 4 drops of eluent.
- **Do not use eluent outside the one given with the cassettes (same lot number).**
- The diagnosis of an infectious disease cannot be established based on the results of a single test.
- Serological results must be interpreted according to available information (e.g., epidemiology, clinical, imaging, biology, etc.) in order to establish a diagnosis.

PERFORMANCES

Sensitivity (Se)

The retrospective evaluation covered 354 samples.

166 samples were clinically documented (eggs, biopsy, clinical features...) from *S. haematobium* and *S. mansoni* infected patients (approximate 1:1 ratio).

Sample's nature	Positive ICT	Negative ICT	Se %
Clinical schistosomiasis (n=166)	164	2	98.8

Table 1: Sensitivity of *SCHISTOSOMA* ICT IgG-IgM on samples from schistosomiasis clinical cases. The 2 negative ICT results were also negative in the reference immunoblot *SCHISTO* II WB IgG.

43 samples were obtained from the recent epidemic in Corsica (hybridization of *S. haematobium* by *S. bovis*). 145 samples came from routine diagnostic samples featuring positive or equivocal results in one or several screening tests (HIA, IFA, ELISA).

Sample's nature	Positive ICT	Negative ICT	Se %
Corsican epidemic : positive samples (n=43)	39	4	90.7
Serological schistosomiasis (n=145)	136	9	93.8
TOTAL (n=188)	175	13	93.1

Fig. 2: Sensitivity of *SCHISTOSOMA* ICT IgG-IgM on serological schistosomiasis samples. The 188 samples were confirmed positive in reference immunoblot *SCHISTO* II WB IgG.

Sensitivity of **SCHISTOSOMA ICT IgG-IgM** (n=354): **Se = 95.8% [93.0 - 97.5%]** (confidence interval according to Wilson's method with continuity correction)

Specificity (Sp)

The evaluation covered 275 samples, including 53 blood donors, 181 sera from patients suffering from the following parasitic infections: *Echinococcus granulosus* (14), *E. multilocularis* (8), malaria (26), filariasis (7), *Stongiloïdes stercoralis* (8), *Fasciola hepatica* (8), *Trichinella spiralis* (9), *Toxocara canis* (26), cysticercosis (45), *Leishmania infantum* (30) and 41 sera from patients suffering from autoimmune diseases: rheumatoid factor RF+ (20), anti-nuclear antibodies ANA+ (21).

Sample's nature	Negative ICT	Positive ICT	Sp%
Blood donor (n=53)	53	0	100.0
Cystic echinococcosis (n=14)	11	3	72.7
Alveolar echinococcosis (n=8)	8	0	100.0
Malaria (n=26)	24	2	91.7
Filariasis (n=7)	7	0	100.0
Anguillulosis (n=8)	7	1	85.7
Fasciolosis (n=8)	7	1	85.7
Trichinellosis (n=9)	9	0	100.0
Toxocariasis (n=26)	26	0	100.0
Cysticercosis (n=45)	35	10	71.4
Leishmaniosis (n=30)	27	3	88.9
Rhumatoid factor (n=20)	20	0	100.0
Anti-nuclear antibodies (n=21)	20	1	95.0
TOTAL (n=275)	254	21	92.4

Table 3: Comparative specificity of **SCHISTOSOMA ICT IgG-IgM** on samples from blood donors and patients suffering from parasitic infections or auto-immune diseases. The 275 samples were found to be negative in reference Immunoblot **SCHISTO II WB IgG**.

Two helminth infections, cystic echinococcosis and cysticercosis, are giving cross reactions with **SCHISTOSOMA ICT IgG-IgM**: This highlights the interest of the immunoblot in making the confirmation of the serological diagnostic.

Specificity of **SCHISTOSOMA ICT IgG-IgM** (n=275): **Sp = 92.4% [88.4 - 95.1%]** (confidence interval according to Wilson's method with continuity correction).

Conclusion

Correlation between clinical and serological results and ICT are very good.

Sensitivity Se = 95.8% [93.0 - 97.5%]

Sp specificity = 92.4% [88.4 - 95.1%]

However, as with any serological result, a positive result does not differentiate an active infection from a treated previous infection.

Reproducibility

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum correlation is excellent.

Interferences

Even though no particular cross-reaction has been observed with haemolysed or lipidic sera, it is recommended to interpret the results from the use of such samples with care. **Icteric sera:** Supplementation trials have shown a possible cross-reaction (false positive) for bilirubin concentration above 100 µmol/L.

BIBLIOGRAPHIE/BIBLIOGRAPHY

- Beltrame A, Guerriero M, Angheben A, Gobbi F, Requena-Mendez A, Zammarchi L, *et al.* 2017. « Accuracy of parasitological and immunological tests for the screening of human schistosomiasis in immigrants and refugees from African countries: An approach with Latent Class Analysis ». *PLoS Negl Trop Dis* 11(6): e0005593. doi:10.1371/journal.pntd.0005593.
- Bevilacqua N, Pane S, Vairo F, Nicastrì E, Paglia M, Ame S, Sañé Schepisi M, *et al.* 2012. « Accuracy of Indirect Haemagglutination and Western Blot Assays for the Detection of Anti-Schistosoma Antibodies in Non-Severe Febrile Patients in Two Tanzanian Hospitals ». *Scandinavian Journal of Infectious Diseases* 44 (6): 453–58. doi:10.3109/00365548.2011.645505.
- Brunet J, Pfaff A, Hansmann Y, Gregorowicz G, Pesson B, Abou-Bacar A, *et al.* Candolfi E. 2015. « An Unusual Case of Hematuria in a French Family Returning from Corsica ». *International Journal of Infectious Diseases: IJID: Official Publication of the International Society for Infectious Diseases* 31 (février): 59–60. doi:10.1016/j.ijid.2014.10.024.
- Cavalcanti M, Silva L, Peralta R, Barreto M, *et al.* Peralta J. 2013. « Schistosomiasis in Areas of Low Endemicity: A New Era in Diagnosis ». *Trends in Parasitology* 29 (2): 75–82. doi:10.1016/j.pt.2012.11.003.
- Colley D, Bustinduy A, Secor E, *et al.* King CH. 2014. « Human Schistosomiasis ». *Lancet* 383 (9936): 2253–64. doi:10.1016/S0140-6736(13)61949-2.
- De Laval F, Savini H, Biance-Valero E, *et al.* Simon F. 2014. « Human Schistosomiasis: An Emerging Threat for Europe ». *The Lancet* 384 (9948): 1094–95. doi:10.1016/S0140-6736(14)61669-X.
- ECDC Stockholm 2014: « Rapid risk assessment: Local transmission of Schistosoma haematobium in Corsica, France ». European Centre for Disease Prevention and Control.
<http://ecdc.europa.eu/en/publications/Publications/schistosoma-haematobium-risk-assessment-France-Germany.pdf>
- Holtfreter MC, Moné H, Müller-Stöver I, Mouahid G, *et al.* Richter J. 2014. « Schistosoma Haematobium Infections Acquired in Corsica, France, August 2013 ». *Euro Surveillance: Bulletin Européen Sur Les Maladies Transmissibles = European Communicable Disease Bulletin* 19 (22).
- Moné H, Holtfreter MC, Allienne JF, Mintsá-Nguéma R, Ibikounlé M, Boissier J, Berry A, Mitta G, Richter J, *et al.* Mouahid G. 2015. « Introgressive Hybridizations of Schistosoma Haematobium by Schistosoma Bovis at the Origin of the First Case Report of Schistosomiasis in Corsica (France, Europe) ». *Parasitology Research*, août. doi:10.1007/s00436-015-4643-4.
- Noormahomed EV, Nhacupe N, Mascaró-Lazcano C, Natane Mauaie M, Buene T, Abel Fonzamo C, *et al.* Benson C. 2014. « A Cross-Sectional Serological Study of Cysticercosis, Schistosomiasis, Toxocariasis and Echinococcosis in HIV-1 Infected People in Beira, Mozambique ». Édité par Ana Flisser. *PLoS Neglected Tropical Diseases* 8 (9): e3121. doi:10.1371/journal.pntd.0003121.
- Sulahian A, Garin Y, Izri A, Verret C, Delaunay P, Van Gool P, *et al.* Derouin F. 2005. « Development and evaluation of a Western blot kit for diagnosis of schistosomiasis ». *Clinical and diagnostic laboratory immunology* 12 (4): 548–51. doi:10.1128/CDLI.12.4.548-551.2005.

UPDATE NOTIFICATION – Please read carefully

RELEASE DATE	VERSION	MODIFICATION SUMMARY
16/04/2021	Vs 08	Typo correction : color of the line (fr/en) – sample volume (en) – “Toxoplasma” (en)
14/06/2021	Vs 09	Add reference for 100 tests kit – contact email address – - Biblio – Bilirubin interference - updated eluent toxicity
30/11/2022	Vs10	New address
02/11/2024	Vs11	Colour of product name and references

NF EN ISO 13485

24 Av. Joannes MASSET – 69009 LYON – FRANCE
 Tel : +33(0)4 7883 3487 – Fax : +33(0)4 7883 3430
www.ldbiodiagnostics.com – info@ldbiodiag.com