TOXOPLASMA

CE0459

ICT IgG-IgM

Test immuno-chromatographique de diagnostic *in vitro*
Technique manuelle

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TOXOPLASMA ICT IgG-IgM

INTENDED USE

TOXOPLASMA ICT IgG-IgM is a single-use, qualitative immunochromatographic test (lateral flow). This rapid test is intended for professional use as a serological screening test for toxoplasmosis by simultaneous detection of anti-Toxoplasma gondii IgG and IgM, in any person with suspected toxoplasmosis, except in the context of organ transplantation. It can be performed on human serum, plasma or whole blood.

PRINCIPLE OF THE TEST

TOXOPLASMA ICT IgG-IgM is a single use unitary test for a qualitative diagnostic. It is based on the principle of the homogeneous sandwich (immunological reaction of two same antigen epitopes with the two binding sites of a bivalent antibody).

Inside the cassette, the device is composed of:

- A nitrocellulose strip on which are spread two reactive bands: the antigens (*Toxoplasma gondii*) of the “test” band (T band) and the rabbit gamma globulins of the “control” band (C band),
- A fiberglass support (conjugate pad) which is impregnated of **black latex** particles coupled with *T. gondii* antigens (“test” latex = T latex) and **blue latex** particles coupled with goat anti-rabbit IgG (“control” latex = C latex).

The test is run by successively dispensing the sample (serum, plasma or whole blood) and an eluting solution (called the eluent) in the “sample well” of the cassette. Adding the eluent starts the concomitant migration (chromatography) of the serum and the latex particles. This migration is completed in 20-30 minutes.

If specific antibodies (IgG and/or IgM) are present in the sample, a complex is formed between the T latex and the patient’s antibodies which is then captured by the T band. It results in the appearance of a **black line**: the test is positive.

The direct capture of the C latex by the C band results in the appearance of a blue line, meaning that the chromatography performed well. The appearance of this **blue line** is systematic and independent of the serological status of the patient.

Both letters “T” and “C” are printed on the cassette, materializing the position of the corresponding reading area.

KIT COMPONENTS

ID	Description	Packaging	
		20 tests	100 tests
R1	Bag (sealed + zip) of 10 ready-to-use cassettes + a desiccant	2	10
R2**	Dropper bottle of 3 mL elution buffer	1	5
	Instructions for use	1	1

** Eluent R2

- **Hazard pictograms**
- **Signal word:** Warning!

Code	Hazard statements
H317	May cause an allergic skin reaction
Code	Prevention
P261	Avoid breathing spray
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves and protective clothing
Intervention	
P305+P351+P338 P337+P313	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention
Disposal	
P501	Dispose of contents and container in accordance with relevant regulations

- **Hazard statements**
- **Precautionary statements**

EUH 210 Safety data sheet available on request as well as on our website www.ldbiodiagnostics.com.

STORAGE AND STABILITY

- **Store the original sealed bag between 2 and 8°C.** Cassettes can be used until the expiry date written on the bag label. Do not freeze. Do not use after the expiration date.
- **The first opening of a bag of 10 tests must occur at least 15 minutes after** putting the bag at room temperature in order to avoid condensation in the bag.
- **Let the eluent at least 15 minutes at room temperature before use.**
- **After the first opening** of a bag, keep it **at room temperature (18-30°C)**, carefully closed (zip closure), the desiccant packet inside. After opening, the cassettes can be used **for up to 2 months**.
- The eluent is stable up to 2 months at room temperature (18-30°C) and until expiration date (as written on the kit) if kept between 2 and 8°C.

PRECAUTIONS FOR USE

Safety

- Safety data sheet available on request as well as on our website www.ldbiodiagnostics.com.
- For *in vitro* use only. Handle according to Good Laboratory Practices and consider any reagent and any sample as potentially toxic and/or infectious.
- For professional use only. Only for technically trained personnel.
- All human samples must be considered as potentially infectious and handled with care.
- Wear a lab coat, gloves and glasses; do not drink, eat or smoke in the laboratory. Do not mouth the pipettes.
- Dispose of waste (samples, tips, tubes, cassettes, used reagent...) according to good practices used in the industry and current regulations in the country.
- Any serious incident must be declared to the manufacturer and the competent authority.

Precautions

- Read and interpret the results under direct white light.
- Do not use eluent from another lot number.
- Do not use cassettes from two different lot numbers in the same run.
- Close the vials after use; do not use if a substance was accidentally introduced in the reagents. Do not use reagent from a vial that presents signs of leakage. Do not use cloudy or precipitated solution.
- Use only disposable pipette tips. Avoid any inter-cassette contamination.
- Do not use reagents after their expiration date.

- The omission of a sample or the distribution of an inadequate volume may render the test result positive or negative, regardless of its actual serological status.
- Only use cassettes carefully stored in their closed bag, the desiccant packet inside.

SAMPLE COLLECTION AND PREPARATION

- The test can be done with either serum or plasma.
 - Sample collection must be sterile and can be done either on dry tube or with heparin, citrate or EDTA.
 - On tubes with gel, do not collect gel: it could lead to false positives.
 - Avoid hemolysis as much as possible.
 - Keep the samples at 2-8 °C until they are processed. If they need to be stored, freeze the samples below -15°C. Do not use a contaminated sample. Avoid freezing and thawing the samples repeatedly.
- Test can also be performed on whole blood
 - In that case, sampling must be performed on a previously cleaned skin according to applicable good sampling practices.

TEST PROCEDURE

Serum or plasma

Additional material required: micropipette and disposable tips for dispensing volumes of 15µL, timer.

If the bags of 10 tests are kept at 2-8°C, let them at least 15 minutes at room temperature before opening: the temperature of the bag must reach the room temperature to avoid condensation in the bag.

1. Take out the desired number of cassettes, then close carefully the bag with the zip closure (the desiccant packet inside) while pressing out as much air as possible. Close and store the bag at room temperature for up to 2 months.
2. **Identify each cassette** with the reference of each sample to be tested. Do not work with runs of more than 10 cassettes. Two successive runs of 10 cassettes must be separated by a few minutes in order to be able to make the reading at the indicated times (the use of 2 timers is recommended).
3. Use a micropipette with a disposable tip to dispense 15µL of serum or plasma in the sample well. Do so for all cassettes before moving to the next step.
4. Dispense **4 drops** of the eluent of the kit. **Do not use the eluent of another lot number.** Keep the dropper vertically while dispensing. Close the dropper after use.
5. Start the timer when the eluent is dispensed in all the cassette of the run.

Whole blood

Additional material required: sterile lancets, capillaries or micropipette and disposable tips for dispensing volumes of 30µL, timer.

If the bags of 10 tests are kept at 2-8°C, let them at least 15 minutes at room temperature before opening: the temperature of the bag must reach the room temperature to avoid condensation in the bag.

1. Do not work in series but one patient at a time
2. After removing the cassette, close carefully the bag with the zip closure (the desiccant packet inside) while pressing out as much air as possible. Close and store the bag at room temperature for up to 2 months.
3. **Identify the cassette.** Do not work in series. Two successive samples must be separated by a few minutes in order to be able to make the reading at indicated times (the use of 2 timers is recommended).
4. With an adequate sterile lancet, proceed to fingerprick (middle or ring fingers recommended): sample 30µl of blood (1/2 a 60µl capillary tube for instance) and dispense immediately in the sample well by capillarity.

5. Dispense **4 drops** of the eluent of the kit. **Do not use the eluent of another lot number.** Keep the dropper vertically while dispensing. Close the dropper after use.
6. Start the timer after the eluent has been dispensed.

READING AND INTERPRETATION

The reading must be done near a window or under direct light (example: a desk lamp). Avoid shadows on the reading area.

The reading must be done between 20 and 30 minutes after starting the timer.

Do not take into account the results from readings after 30 minutes.

- **Positive test:** 2 lines, a **black "T"** and a blue **"C"** appear in the corresponding areas. Every **"T"** line must be considered positive, even of very weak intensity. For very weak lines, make the reading with the eye vertically above the reading area.
- **Negative test:** No black line appears. Only the blue **"C"** line is visible.
- **Equivocal test:** In very rare cases, a faint, diffuse, grey line can appear on the **"T"** band. This result should be considered negative but controlled on another sample or technique.
- **Invalid test:** The **"C"** line does not appear. Read once again the instruction and repeat the test. If the problem persists, contact the manufacturer or your distributor.

Notes: This is a qualitative test. Intensity of the black line does not reflect the quantity of anti-*Toxoplasma* antibodies in the sample.

Positivity of the test is proof of the contact of the patient with the infectious agent but doesn't prejudge the contact date or the clinical status of the patient.

QUALITY CONTROL

The blue **"C"** line allows the validation of the good running of the test. However, it is recommended to incorporate from times to times a known weak positive sample in a run.

TEST LIMITATIONS

- Do not use too old serum sample with this technique. It is recommended to use samples frozen for less than 2 years.
- The positivity can be caused by the presence of IgG and/or IgM directed against the infectious agent, the test does not distinguish the type of antibodies present.
- The use of other body fluids (urine, CSF, saliva...) has not been validated.
- The use of hemolytic, icteric or lipidic samples is not recommended. However, no interference in the reaction was recorded. Still, very hemolytic samples can hide a weak positive test, due to important red background from hemolysis.
- Do not dispense 3 or 5 drops of eluent.
- **Do not use eluent outside the one given with the cassettes (same lot number).**
- The diagnosis of an infectious disease cannot be established based on the results of a single test.
- Serological results must be interpreted according to available information (e.g., epidemiology, clinical, imaging, biology, etc.) in order to establish a diagnosis.

PERFORMANCES (see literature references)

Sensitivity, Specificity (serum/plasma)

Material and methods

The evaluation was done in a Reference Laboratory specialising in the diagnosis of toxoplasmosis and member of the French National Reference Centres network (CNR) for toxoplasmosis.

The serological status (positive or negative) was the final result of the routine diagnostic in the laboratory from the patient's biological and epidemiological data. The discordant biological results were solved by the use of confirmatory techniques LDBIO-TOXO II IgG, ISAGA IgM and Sabin and Feldman's DYE TEST).

The principle of the evaluation consisted of comparing the results of 486 samples (including 67 cord blood samples) obtained with the LDBIO ICT and the results obtained with 3 commercial serological tests: ELISA IgG, ELISA IgM (Abbott Architect) and Indirect Hemagglutination Assay – IHA (Fumouze).

The results of ELISA IgG and ELISA IgM were compiled (ELISA G+M) in order to have a better comparison as both LDBIO ICT and IHA detect simultaneously IgG and IgM class antibodies.

Two studies were conducted in parallel:

- A prospective study on 356 samples from the routine activity of the laboratory.
- A retrospective study on 130 samples selected for their specific profile and spread into 5 groups:
 - Group 1 – seroconversion: retrospective analysis of 9 serum sequence (30 samples) from patients having presented a toxoplasmosis seroconversion during their pregnancy.
 - Group 2 – monitoring of uninfected children: This is a retrospective analysis of 15 samples corresponding to 5 sequences of post-natal monitoring of children born to mothers who presented a toxoplasmosis seroconversion during pregnancy.
 - Group 3 – false positive with ELISA: 4 false positive with ELISA IgG and 11 false positive ELISA IgM including 1 false positive and 6 equivocal results by ISAGA IgM.
 - Group 4 – particular populations:
 - residual IgM (n=5),
 - weak IgG level (n=20), selected in ELISA IgG with antibodies level between 1.1 and 3.3 UI/mL (note: of these 20 samples, 2 were therefore considered negative and 16 equivocal according to the criteria of the ELISA used in the study)
 - strong IgG level (n=10),
 - serological reactivation (n=10),
 - sera from multi-organ transplant donor (n=10, 9 positives, 1 negative).
 - Group 5 – potential cross reactions: sera negative for toxoplasmosis but with anti-CMV and/or anti-EBV (n=10) and sera with rheumatoid factor (RF, n=5).

Results

- Whole population (prospective + retrospective studies): n=486

	ELISA G+M	IHA	LDBIO ICT		ELISA G+M	IHA	LDBIO ICT
Positive	151	175	189	Positive	16	2	7
Negative	9	14	0	Negative	276	295	290
Equivocal	29	0	0	Equivocal	5	0	0
Se	94.6%	92.9%	100%	Sp	94.6%	99.3%	97.7%

Table 1: results of the different techniques compared to the serological status (n=486) for the entire study. Equivocal results in ELISA were excluded from sensitivity and specificity calculations.

- Prospective study: n=356

	ELISA G+M	IHA	LDBIO ICT		ELISA G+M	IHA	LDBIO ICT
Positive	88	96	105	Positive	1	0	6
Negative	4	9	0	Negative	249	251	245
Equivocal	13	0	0	Equivocal	1	0	0
Se	95.7%	91.4%	100%	Sp	99.6%	100%	97.6%

Table 2: results of the different techniques compared to the serological status (n=486) for the prospective study. Equivocal results in ELISA were excluded from sensitivity and specificity calculations.

- Retrospective study n=130
 - Group 1 – seroconversion:
Of the 9 sequences, 1 is detected earlier by the LDBIO ICT. All other results were concordant for all three techniques.
 - Group 2 – monitoring of uninfected children:
LDBIO ICT allowed the detection of a serological residual positivity on a serum already negative with the two other techniques. All other results were concordant for all three techniques.
 - Group 3 – false positive ELISA (NEGATIVE serological status) n=15:

	IHA		ISAGA IgM			LDBIO ICT	
	POS	NEG	POS	EQUIVOCAL	NEG	POS	NEG
False positive ELISA IgG	0	4	NA	NA	NA	0	4
False positive ELISA IgM	2	9	1	6	4	0	11

Table 3: Compared performances of the different techniques for the false positive serum in ELISA.

- Group 4 – particular populations:
With the exception of the weak IgG level samples, the results were concordant between all the techniques. No prozone effect was observed in any of the three techniques, even with 2000UI/mL IgG level.
For the weak IgG level (POSITIVE serological status), the results are as follow: n=20

	ELISA G+M			IHA		LDBIO ICT	
	POS	EQUIVOCAL	NEG	POS	NEG	POS	NEG
IgG between 1.1 and 3.3 UI/mL	2	16	2	18	2	20	0

Table 4: Performances of the different techniques for the weak IgG level group

- Group 5 – potential cross reactions:
Of the 10 sera positive for EBV and/or CMV, none were found positive with ICT and IHA, but one was found equivocal with ELISA. Of the 5 sera positive for rheumatoid factor, one was false positive in ICT and equivocal in ELISA but correctly identified by IHA. All others samples were correctly identified.
There may be a potential cross reaction with this population, but a study on a greater number of serums is needed in order to conclude.

Note: Weak false positive reactions have been observed on ~10% samples positive for *Leishmania* serology. This could explain some of the nonspecific reactions that were recorded.

Conclusions

Correlation between the serological status and LDBIO ICT is excellent:

Sensitivity Se=100%, [IC95: 95.6 – 100 %]

Specificity Sp=97.7%, [IC95: 94.9 – 99.0 %]

Confidence intervals were done using Aggresti-Coull method.

LDBIO ICT has a low detection limit (<2UI/mL) without equivocal results.

The analysis of the 67 cord blood samples did not give any discordant results between ICT and the conclusion, showing that LDBIO ICT could be use on such samples.

A positive result in LDBIO ICT cannot discriminate between IgM, IgG or the simultaneous presence of IgG and IgM.

Sensitivity, Specificity (whole blood)

Material and methods

This study was carried out by the National Institute of Hygiene in Morocco. 349 patients (female sex, >15 years old) were enrolled at 3 sites. Patients had a capillary sample taken through a lancet and a 60µl collection capillary, half filled. The sample was used for the point-of-care **Toxoplasma ICT IgG-IgM** whole blood test. They also had a blood sample taken from venipuncture in order to perform a **Toxoplasma ICT IgG-IgM** test on serum in the laboratory after centrifugation. A **WB LDBIO Toxo II IgG** confirmation test was also performed.

Results

		Toxoplasma ICT IgG-IgM serum		Total
		Pos	Neg	Total
Toxoplasma ICT IgG-IgM whole blood	Pos	108	1	109
	Neg	4	236	240
	Total	112	237	349

Table 6: Results of whole blood evaluation

The correlation between serum results and WB is 100%.

Conclusion

The performance of the kit in point-of-care use with capillary and whole blood samples is excellent with a sensitivity of 96.4% [95CI 91.1-99.0%] and a specificity of 99.6% [95CI 97.7-99.9%].

Confidence intervals are calculated according to the Aggresti-Coull method.

Reproducibility

Inter-series and inter-lot reproducibility were tested. In both cases, the serum to serum, and capillarity sampling to capillarity sampling, correlation is excellent.

Interferences

Even though no particular cross-reaction has been observed with haemolysed, icteric or lipidic sera, it is recommended to interpret the results from the use of such samples with care.

BIBLIOGRAPHIE/BIBLIOGRAPHY

- Begeman I, Lykins J, Zhou Y, Shiun Lai B, Levigne P, El Bissati K, Boyer K, *et al.* 2017. « Point-of-Care Testing for Toxoplasma Gondii IgG/IgM Using Toxoplasma ICT IgG-IgM Test with Sera from the United States and Implications for Developing Countries ». *PLoS Neglected Tropical Diseases* 11 (6):e0005670. <https://doi.org/10.1371/journal.pntd.0005670>.
- Chapey E, Wallon M, et Peyron F. 2017. « Evaluation of the LDBIO Point of Care Test for the Combined Detection of Toxoplasmic IgG and IgM ». *Clinica Chimica Acta; International Journal of Clinical Chemistry* 464 (janvier):200-201. <https://doi.org/10.1016/j.cca.2016.10.023>.
- Franck J, Garin Y, et Dumon H. 2008. « LDBio-Toxo II Immunoglobulin G Western Blot Confirmatory Test for Anti-Toxoplasma Antibody Detection ». *Journal of Clinical Microbiology* 46 (7):2334-38. <https://doi.org/10.1128/JCM.00182-08>.
- Jost C, Touafek F, Fekkar A, Courtin R, Ribeiro M, Mazier D, et Paris L. 2011. « Utility of Immunoblotting for Early Diagnosis of Toxoplasmosis Seroconversion in Pregnant Women ». *Clinical and Vaccine Immunology* 18 (11):1908-12. <https://doi.org/10.1128/CVI.05303-11>.
- Khammari I, Saghrouni F, Lakhal S, Bouratbine A, Ben Said M, et Boukadida J. 2014. « A New IgG Immunoblot Kit for Diagnosis of Toxoplasmosis in Pregnant Women ». *The Korean Journal of Parasitology* 52 (5):493-99. <https://doi.org/10.3347/kjp.2014.52.5.493>.
- Khammari I, Saghrouni F, Yaacoub A, Gaied Meksi S, Ach H, Garma L, Fathallah A, et Ben Saïd M. 2013. « IgG Western Blot for Confirmatory Diagnosis of Equivocal Cases of Toxoplasmosis by EIA-IgG and Fluorescent Antibody Test ». *The Korean Journal of Parasitology* 51 (4):485-88. <https://doi.org/10.3347/kjp.2013.51.4.485>.
- Mahinc C, Flori P, Delaunay E, Guillerme C, Charaoui S, Raberin H, Hafid J, et L'Ollivier C. 2017. « Evaluation of a New ICT Test (LDBIO Diagnostics) to Detect Toxoplasma IgG and IgM: Comparison with the Routine Architect Technique. » *Journal of Clinical Microbiology*, septembre, JCM.01106-17. <https://doi.org/10.1128/JCM.01106-17>.
- Villard O, Cimon B, L'Ollivier C, Fricker-Hidalgo H, Godineau N, Houze S, Paris L, Pelloux H, Villena I, et Candolfi E. 2016. « Help in the Choice of Automated or Semiautomated Immunoassays for Serological Diagnosis of Toxoplasmosis: Evaluation of Nine Immunoassays by the French National Reference Center for Toxoplasmosis ». *Journal of Clinical Microbiology* 54 (12):3034-42. <https://doi.org/10.1128/JCM.01193-16>.

UPDATE NOTIFICATION – Please read carefully

RELEASE DATE	VERSION	MODIFICATION SUMMARY
01/10/2020	Vs 10	Layout - Reagents supplied – Precaution for use – Stability – addition of whole blood use
10/06/2021	Vs 11	Add reference for 100 tests kit – contact email address - updated eluent toxicity
30/11/2022	Vs12	New address
13/08/2024	Vs13	Updated test indication - SDS § safety - color designation and references

NF EN ISO 13485

24 Av. Joannes MASSET – 69009 LYON – FRANCE
 Tel : +33(0)4 7883 3487 – Fax : +33(0)4 7883 3430
www.ldbiodiagnostics.com – info@ldbiodiag.com